

Cleaning and Care

Remove soft components and wipe shell clean with a damp cloth. Liner and straps can be hand wash with mild detergent. Triple rinse, squeeze and air dry. Avoid twisting the material to remove water.

Clinical Fitting Instructions

Select size according to leg circumference. Measure leg circumference 6" above mid-patella, at knee center, and at 6" below mid-patella.

Adult Sizing Guide

Size	Thigh	Knee Joint	Calf
Small	15"-18"	11"-12.75"	12.5"-13"
Medium	16"-20"	12.75"-15.25"	13.5"-16"
Large	18"-24"	15.25"-16.625"	14.5"-18"
XLarge	18"-26"	16.625"-18.125"	14.5"-19"

Liability/Warranty

This device is warrantied for a period of one year. The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use. Details of the warranty are also available at www.ottobockus.com

Manufacturer's Limited Warranty

Otto Bock HealthCare, LP (OttoBock) warrants all of its devices, to the original purchaser, to be free from defects in materials and workmanship. This warranty applies, subject to normal wear and tear, when the devices are used as intended, without unapproved modifications, following all OttoBock instructions and requirements; and when they are fitted by or under the direct supervision of certified/licensed practitioners. This Limited Warranty does not cover device damage caused by accidents, neglect, misuse or operation beyond capacity, parts damaged by improper installation, substitution of parts not approved by OttoBock, or any alteration or repair by others that, in OttoBock's judgment, materially or adversely affect the device.

The duration of this Limited Warranty varies by product types and is effective from the date of delivery to the end-user. Please refer to www.ottobockus.com or call 800 328 4058 for questions. OttoBock's sole obligation under this Limited Warranty shall be to repair, replace, refabricate the device at no charge, or refund the cost of the device to the original purchaser, at OttoBock's sole discretion.

THE EXPRESS WARRANTIES SET FORTH ABOVE ARE IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, AND ALL SUCH OTHER WARRANTIES ARE HEREBY DISCLAIMED AND EXCLUDED BY OTTOBOCK. IN NO EVENT SHALL OTTOBOCK'S LIABILITY OF ANY KIND INCLUDE ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, EVEN IF OTTOBOCK SHALL HAVE BEEN ADVISED OF THE POSSIBILITY OF SUCH POTENTIAL LOSS OR DAMAGE.

Some states do not allow the exclusion of incidental or consequential damages, thereby rendering the aforementioned limitation in applicable to certain original purchasers.

Otto Bock HealthCare LP
 Two Carlson Parkway North, Suite 110
 Minneapolis MN, 55447
 P 800 328 4058 · F 800 962 2549
 Made in India

Ottobock LTE Sport Brace
 Instructions for Use



50K301=* Ottobock LTE Sport Brace

Thank you for using this Ottobock device. The following information will help you fit and care for this device. Please review it before using this device and ask for clarification if needed.

Indications for Use:

Indicated for support of the ACL, PCL, MCL, and LCL following injury or surgical reconstruction.

Warning

This device does not prevent injury and is not intended to reduce or eliminate the risk of injury.

Warning

Always consult with your physician or healthcare provider before making changes to the brace.

Warning

Carefully read all instructions for use before using or fitting this product.

Warning

Proper rehabilitation and activity modification are also a part of a safe treatment program.

Warning

If you experience pain, swelling, or any unusual reaction, contact your physician or healthcare provider immediately.

CAUTION

This device is offered under order of a physician or other qualified healthcare provider.

CAUTION

This device is intended for single patient use and is not intended to be reused on a second patient.

Application Instructions

1. Unfasten straps by disengaging the hook and loop fasteners on the strap end and pulling the strap through the d-ring.
2. With the leg bent between 45° and 60° (Fig. 1) position the brace on the leg so the hinges are aligned with the top of the knee cap.

For PCL patients: With the leg positioned in full extension, position the brace on the leg so the hinges are aligned with the top of the knee cap.

3. Secure the upper calf strap first (Fig. 2). Pull the strap through the d-ring and secure using the hook and loop fastener. The strap should be snug. Do not over tighten the strap.
4. Secure the lower calf strap (Fig. 3). Pull the strap through the d-ring and secure using the hook and loop fastener. The strap should be snug. Do not over tighten the strap.
5. Secure the lower thigh strap (Fig. 4). Pull the strap through the d-ring and secure using the hook and loop fastener. The strap should be snug. Do not over tighten the strap.
6. Secure the upper thigh strap (Fig. 5). Pull the strap through the d-ring and secure using the hook and loop fastener. The strap should be snug. Do not over tighten the strap.
7. Secure the front shin strap being careful not to pull the brace forward. Pull the strap through the d-ring and secure using the hook and loop fastener (Fig. 6). The strap should be snug. Do not over tighten the strap.
8. Adjust strap tensions as needed and ensure the hinges remain aligned with the top of the knee cap.



Fig. 1



Fig. 2



Fig. 3



Fig. 4



Fig. 5



Fig. 6