





50R300 Dyneva

DE Gebrauchsanweisung	6
EN Instructions for use	11
FR Instructions d'utilisation	15
IT Istruzioni per l'uso	20
ES Instrucciones de uso	25
PT Manual de utilização	30
NL Gebruiksaanwijzing	35
SV Bruksanvisning	40
DA Brugsanvisning	44
NO Bruksanvisning	49
FI Käyttöohje	53
PL Instrukcja użytkowania	58
HU Használati utasítás	62
CS Návod k použití	67
SK Návod na používanie	72
TR Kullanma talimatı	76
ZH 使用说明书	80
KO 사용 설명서	84

	Größe / Size	Abstand / Distance	
		cm	inch
	XS	a ~ 27; b ~ 30	a ~ 10.6; b ~ 11.8
	S	a ~ 30; b ~ 33	a ~ 11.8; b ~ 13.0
	M	a ~ 33; b ~ 36	a ~ 13.0; b ~ 14.2
	L	a ~ 36; b ~ 39	a ~ 14.2; b ~ 15.4
	XL	a ~ 39; b ~ 42	a ~ 15.4; b ~ 16.5
	XXL	a ~ 42; b ~ 45	a ~ 16.5; b ~ 17.7

Material	PA, TPU, PU, TPE-E, PE-Schaum / PE foam, Federstahl / Spring steel, Elastan / Spandex, PE, Edelstahl / stainless steel
-----------------	--

INFORMATION

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INFORMATION

Last update: 2016-04-21

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use contain important information for fitting and applying the 50R300 Dyneva lumbar spine flexion brace.

2 Product description

Item	Component (see fig. 1)	Item	Component (see fig. 1)
1	Connecting rods	2	Hook-and-loop closure for belt
3	Tension belt	4	Abdominal strap closure
5	Abdominal strap segments	6	Lower tube frame
7	Sacral pad	8	Spring adapter with abdominal strap element
9	Cover caps	10	Upper tube frame
11	Shoulder pads		

The product is delivered completely assembled (see fig. 2).

3 Application

3.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the back and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

3.2 Indications

With pain syndromes of the lumbar spine caused by:

- Lumbar spinal stenosis
- Lumbar facet joint syndrome
- Lumbar radicular syndrome (anterior thigh pain, lumboschialgia)
- Lumbar instability (spondylolisthesis, spondylolysis, rotational listhesis, degenerative spondylolisthesis)
- Lumbar spondylarthrosis
- Intervertebral disc prolapse and protrusion in the lumbar spine area
- General degenerative changes with reduced intervertebral disc height

The indication must be determined by the physician.

3.3 Contraindications

3.3.1 Absolute Contraindications

None known.

3.3.2 Relative Contraindications

Report to and discuss with the physician any of the following: contractures of the spine, restricted chest breathing ability, skin diseases/injuries, inflammation, raised scars with swelling, reddening, and warming in the area where the brace is fitted, tumors, lymphatic congestion – including unclear swelling of soft tissue distal to the medical aid, and dysesthesia in the legs.

3.4 Mechanism of Action

The three-point principle and the spring mechanism result in an effective and statistically significant reduction of the muscle-force-induced compression of the movement segments and the small vertebral joints in the lumbar spine. Spine deformation occurs in both the sagittal and the frontal planes. This applies for moderate as well as maximum loads. The brace effectively, significantly and systematically reduces the load on the lumbar spine.

4 Safety

4.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

4.2 General safety instructions

CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

CAUTION

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

NOTICE

Contact with ointments, lotions, or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to ointments, lotions, or other products that contain oils or acids.

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below $-10\text{ }^{\circ}\text{C}$ ($14\text{ }^{\circ}\text{F}$) or above $+50\text{ }^{\circ}\text{C}$ ($122\text{ }^{\circ}\text{F}$) (e.g. sauna, excessive sunlight, drying on a radiator).

NOTICE

Contact with saltwater or water containing chlorine/soap, fresh water or abrasive substances (e.g. sand)

Damage and premature product wear

- ▶ Following any contact with the substances identified above, promptly clean the product in accordance with the section "Cleaning".

5 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

5.1 Selecting the Size

> Required tools: 743S1=60 diameter or 743E10 measuring tool

- 1) Measure the width of the torso 4 cm below the shoulder blades (see fig. 3).
- 2) Measure the width of the torso at the level of the iliac crest (see fig. 4).
- 3) Determine the size of the tube frames using the measured dimensions and the size chart on page 5.
- 4) **INFORMATION: Determine the size of the brace (XS/S/M/L/XL/XXL) on the basis of the width of the torso below the shoulder blades. If necessary, order the lower tube frame one size smaller or bigger and then replace it (reference number 29R394=XS/S/M/L/XL/XXL).**

5.2 Adaptation

NOTICE

Improper changes

Damage to the product due to overloading of the material and improper fit of the product due to breakage of load-bearing components

- ▶ Do not make any improper changes to the product.

- 1) **NOTICE! Make sure not to damage the pins when attaching or removing the pads.**

Position the brace frame on the back (see fig. 5):

Sacral pad on the sacrum

Shoulder pads approx. 4 cm below the shoulder blades

Optional: To increase wearer comfort, attach and adapt the side pads, which can be ordered separately, to the ends of the tube frames. One set of side pads (29R390) is required for the upper and lower tube frames (see fig. 6).

INFORMATION: To position the shoulder pads vertically, push the upper tube frame inside the spring adapters. The maximum pull-out height is indicated on the tube frame (see fig. 7). For horizontal positioning of the shoulder pads, move them horizontally. Do this by sliding a flat screwdriver under the cover caps, loosening the cover caps and pulling off the shoulder pads. Insert the shoulder pad pin in one of the three holes (see fig. 8).

- 2) After vertical and horizontal positioning, put on the cover caps.

- 3) Check the fit of the tube frames on the patient and make any necessary adjustments by bending the tube frame or replacing the lower tube frame. The replacement size may only be one size smaller or bigger (reference number 29R394=S/M/L/XL/XXL).
- 4) Determine the number of abdominal strap segments required (see fig. 9). Pay attention to the minimum and maximum overlapping of the strap end elements (see fig. 10).
- 5) Remove unnecessary abdominal strap segments by pulling out the connecting rods (see fig. 11).
- 6) Reconnect the abdominal strap segments that are being used with a connecting rod.
- 7) Apply the frame of the brace to the back.
- 8) If necessary, adjust the angle of the abdominal strap at the spring elements after loosening the screws and removing the covers (see fig. 12).
- 9) Pull evenly on the two abdominal strap closures towards the front (see fig. 13).
- 10) Close the abdominal strap (see fig. 14).
- 11) Attach the hook-and-loop closures for the tension belts to the waist belt and determine the optimal length of the tension belts (see fig. 15).
NOTICE! Make sure that the angle of the springs does not exceed 20° when the tension belts are taut.
- 12) Shorten the length of the tension belts with scissors (see fig. 16).

5.3 Application

⚠ CAUTION

Direct skin contact with the product

Skin irritation due to chafing or sweat formation

- ▶ Do not wear this product directly on the skin.

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the brace is applied properly and fits correctly.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- 1) Place the brace frame in such a way that the sacral pad comes to be positioned on the sacrum, and the two shoulder pads below the shoulder blades (see fig. 5).
- 2) Evenly pull the two abdominal strap closures to the front and close them (see fig. 17).
- 3) Evenly tighten the tension belts and attach them to the positions identified during the fitting by means of hook-and-loop closures (see fig. 18).

5.4 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the brace regularly:

- Remove the textile component from the brace.
- Close the hook-and-loop closures.
- Wash the textile component in warm water at 30 °C/104 °F with a standard mild detergent. Rinse thoroughly.
- Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).

6 Disposal

Dispose of the product in accordance with national regulations.

7 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

7.1 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorized modification of the product.

7.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The manufacturer therefore drew up the declaration of conformity on its own responsibility in accordance with Annex VII of the directive.

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Otto Bock HealthCare GmbH