

17B82*

DE	Gebrauchsanweisung	4
EN	Instructions for use	9
FR	Instructions d'utilisation	15
IT	Istruzioni per l'uso	21
ES	Instrucciones de uso	26
PT	Manual de utilização	33
NL	Gebruiksaanwijzing	38
SV	Bruksanvisning	44
DA	Brugsanvisning	49
NO	Bruksanvisning	55
PL	Instrukcja użytkowania	60
HU	Használati utasítás	66
CS	Návod k použití	72
TR	Kullanma talimatı	77
EL	Οδηγίες χρήσης	82
RU	Руководство по применению	89

1

2

INFORMATION

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INFORMATION

Date of last update: 2021-05-10

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.

- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide you with important information on the processing of the 17B82* orthotic hip joint.

2 Product description

2.1 Available sizes

Article number	Orthotic hip joint	System width	Material
17B82=L/R16	Top contoured, bottom straight	16 mm	Stainless steel
17B82=L/R20		20 mm	

2.2 Components/design

Scope of delivery (see fig. 1)			
Item	Quantity [pcs.]	Designation	Reference number
1	1	Upper joint section	
2	1	Screws	501T7=7.5x9xM5
3	1	Joint screw	501A6=5x5xM5
4	1	Ball bearing	509K1=5x16x5
5	1	Lower joint section	

Not included in scope of delivery			
Item	Quantity [Piece (s)]	Designation	Article number
(see fig. 1), item 6	1	Side bar	17B38*, 17B39*, 17B40*, 17B41*, 17B84*
not illustrated	1	Lamination bar	17Y128*, 17Y129*

3 Intended use

3.1 Indications for use

The orthotic hip joint is intended **exclusively** for orthotic fittings of the lower limbs on one patient.

3.2 Indications

- Partial or total paralysis of the leg muscles
- Orthopaedic diseases of the lower limbs

Indications must be determined by the physician.

3.3 Lifetime

The product is designed for a lifetime of **3 years** when used as intended and assembled professionally.

3.4 Qualification

Patients may be fitted with the product only by trained qualified personnel. The qualified personnel must be familiar with the handling of the various techniques, materials, machines and tools.

4 Safety

4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

4.2 Safety instructions

 **CAUTION**

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

 **CAUTION**

Incorrect alignment or assembly

Injuries due to changes in or loss of functionality

- ▶ The product may only be installed by trained, qualified personnel.
- ▶ Observe the alignment and assembly instructions.

⚠ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

⚠ CAUTION

Clothing getting caught in the orthotic joint

Risk of injury due to loss of orthotic joint functionality

- ▶ Make sure your clothing does not get caught in the orthotic joint.

⚠ CAUTION

Excessive strain due to use on more than one patient

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

NOTICE

Thermal overloading of the product

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at temperatures above **300 °C (570 °F)**.
- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +60 °C (140 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

5 Preparing the product for use

CAUTION

Repeated bending at the same position

Risk of injury due to breakage of the bar, change in or loss of functionality

- ▶ Avoid repeated bending at the same position.

Shaping the bars

- > **Required tools:** 711S5 bending irons, 711S3 horizontal bending irons, 711S2 bending pusher set.
- 1) Use the 711S5 bending irons for bending.
- 2) Use the 711S3 horizontal bending irons or the 711S2 bending pusher set for horizontal bending.
- 3) Use a bending radius over 10 mm (see fig. 2).
- 4) Remove grooves and burrs by grinding.
- 5) For final assembly, secure the joint screw with Loctite 241.

Gluing in bars

Required materials: 636W28* special adhesive, degreasing cleaner

The bars must be glued in the insertion zones to ensure stability.

- 1) Clean the connecting surfaces with a degreasing cleaner.
- 2) Apply the special adhesive to the connecting surfaces.
- 3) Insert the bars.
- 4) Secure with the enclosed screws.
- 5) Allow to cure for at least **4 hours**.

INFORMATION: Final bonding strength will be reached after 16 hours.

Improving corrosion protection

Polish or sinter the product for improved corrosion protection. Ottobock recommends the 618T40* sintering powder for surface coating.

- ▶ Do not heat the material longer than **5 minutes** at max. **150 °C** for surface coating.

Using protective devices

- ▶ Use joint protectors to protect the product against soiling.

6 Cleaning

Promptly clean the product after contact with water containing salt, chlorine or soap, or if it gets dirty.

- 1) Rinse the product with clear fresh water.

- 2) Dry the product with a cloth or allow it to air dry. Avoid exposure to direct heat (e. g. from an oven or radiator).

7 Maintenance

INFORMATION

The product may be exposed to increased loads by the patient.

- ▶ Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

Spare parts are listed in the section "Components/design".

Only use 633F7 special lubricant.

Replacing the ball bearing

NOTICE

Canting the ball bearing

Mechanical damage

- ▶ Do not cant or damage the ball bearing while inserting it into the orthotic joint.

- 1) Loosen the screws and disassemble the orthotic joint.
- 2) Remove the ball bearing from the orthotic joint and replace it with a new one.
- 3) Assemble the orthotic joint and tighten the screws.
- 4) For the final assembly, secure the screw connection with Loctite 241.

8 Disposal

Dispose of the product in accordance with national regulations.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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