
4E50-2

DE	Gebrauchsanweisung (Fachpersonal).....	3
EN	Instructions for use (Qualified personnel).....	9
FR	Instructions d'utilisation (Personnel spécialisé).....	14
IT	Istruzioni per l'uso (Personale tecnico specializzato).....	20
ES	Instrucciones de uso (Personal técnico especializado).....	26
PT	Manual de utilização (Pessoal técnico).....	32
NL	Gebruiksaanwijzing (Vakmensen).....	38
SV	Bruksanvisning (Fackpersonal).....	44
DA	Brugsanvisning (Faguddannet personale).....	50
PL	Instrukcja użytkowania (Personel fachowy).....	56
HU	Használati utasítás (szakszemélyzet).....	62
CS	Návod k použití (Odborný personál).....	68
RO	Instrucțiuni de utilizare (Personal de specialitate).....	74
HR	Upute za uporabu (Stručno osoblje).....	80
TR	Kullanma talimatı (Uzman personel).....	86
EL	Οδηγίες χρήσης (Τεχνικό προσωπικό).....	92
RU	Руководство по применению (Квалифицированный персонал).....	97
JA	取扱説明書 (有資格担当者).....	104
ZH	使用说明书 (专业人员).....	109

INFORMATION

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INFORMATION

Last update: 2021-04-14

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- Please keep this document for your records.

Explanation of symbols

▲ DANGER Warnings regarding directly impending risks of severe accident or injury.

▲ WARNING Warnings regarding possible risks of severe accident or injury.

▲ CAUTION Warnings regarding possible risks of accident or injury.

NOTICE Warnings regarding possible technical damage.

1 Intended use

The battery charger is to be used **exclusively** for charging the batteries of Ottobock products for the lower limb.

1.1 Contraindications

All conditions which contradict or go beyond the specifications listed in the section on „Safety“ and „Indications for use“.

2 Safety unstructions

INFORMATION

Advise your patients of the information in this section.

NOTICE

Non-observance of safety instructions. Failure to follow the safety instructions mentioned below can lead to faulty control or malfunction of the charger.

NOTICE

Penetration of dirt and humidity. Penetration of dirt and humidity can lead to faulty control or malfunction of the charger. This makes impossible to ensure an appropriate charging function. Do not let foreign particles or liquids get into the battery charger.

NOTICE

Mechanical overloading. Exterior mechanical impact or stress such as shocks and vibrations can lead to faulty control or malfunction of the battery charger. This makes impossible to ensure an appropriate charging function.

Do not expose the charger to mechanical vibrations or shocks.

Check the charger for any visible damage before every use.

NOTICE

Thermal overloading. Operating the battery charger at the temperature range below 0 °C can lead to faulty control or malfunction of the charger. This makes impossible to ensure an appropriate charging function.

Avoid the temperatures below 0 °C (see also Section "6 Recommendations").

NOTICE

Manipulation of system components. Any changes or modifications you make on the system components can lead to faulty control or malfunction of the charger. This makes impossible to ensure an appropriate charging function.

The charger components may be opened or repaired only by authorized Ottobock Service technicians.

NOTICE

Distance to HF communication devices is too short (e.g. mobile phones, Bluetooth devices, WiFi devices). If the distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi

devices) is too small, interference with internal data communications can result in unexpected prosthesis behaviour. Therefore, keeping the following minimum distances to these HF communication devices is recommended:

- Mobile phone GSM 850/GSM 900: 0.99 m
- Mobile phone GSM 1800/GSM 1900/UMTS: 0.7 m
- DECT cordless phones incl. base station: 0.35 m
- WiFi (routers, access points,...): 0.22 m
- Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

NOTICE

Operating near sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors). Avoid staying near high-tension power lines, transmitters, transformers or other sources of strong electromagnetic radiation (such as security systems for goods in department stores), as this can lead to malfunction of the product.

3 Scope of delivery

- 1 pc. 4E50-2 Battery charger for C-Leg
- 1 pc. 757L16-* AC-adapter
- 1 pc. 4X81-* Charging Cable
- 1 pc. Instruction Manual

4 Function

The battery charger is used for charging the C-Leg / C-Leg compact / C-Brace and Meridium. Charging is performed automatically after the charging plug has been connected to the charging receptacle of the component. LEDs indicate the readiness of the battery charger and the current battery charge level.

5 Battery procedures

A high quality Lithium-Ion battery provides the power necessary. Under normal circumstances this is equal to one full day's use of the electronic knee joint. Please charge the battery after each day's use of the electronic knee joint.

CAUTION

Incorrect behaviour while charging the joint. If the patient walks while the battery charger is connected to the prosthesis, he can get caught on the cable and fall.

Take off the prosthesis prior to charging.

CAUTION

Incorrect behaviour while charging the joint. If the joint is not flexed completely during charging, the charging plug can get damaged and the joint can no longer be charged. This can lead to insufficient power supply to the electronics of the joint and result in undefined conditions. This can cause the patient to fall.

During charging, keep the prosthesis fully flexed.

The 757L16-* AC-adapter for the 4E50-2 Battery Charger may be used with a power supply voltage range of 100 V to 240 V and a power supply frequency of 50 Hz to 60 Hz, making it suitable for use in most countries of the world.

The following procedure is recommended for charging the electronic knee joint's battery: Connect the individual components to the charging unit as shown in fig. 6.

- Connect the 757L16-* AC-adapter to the 4E50-2 Charger (Fig. 1, item B). Plug the 4X81-* Charging Cable into the 4E50-2 charger (Fig. 1, item A, Fig. 2). Plug the AC-adapter into the wall socket. The green LED lights up to indicate the Charger is ready for operation.

INFORMATION

When used daily, the complete charging unit may remain plugged into the wall socket.

- b) Lower the electronic knee joint's protective cover (if present) to expose the plug connections in front of the knee (Fig. 3).
- c) Remove the protective plug (Fig. 4) or move the slide to expose the charging receptacle of the electronic knee joint. Gently insert the charging cable plug (with notch facing forward) into the charging receptacle (Fig. 5). **Do not force it in!** The green and the yellow LED lights of the battery charger will light, verifying that the battery is charging.
- d) When the yellow LED goes out, the battery is fully charged. When the charging cable is unplugged, a short self-test is carried out. As soon as the self-test has been finished, the electronic knee joint is ready to function.

INFORMATION

The battery charger and power supply unit will heat up during the charging process!

CAUTION

Penetration of dirt and humidity. Penetration of dirt and humidity into the charging receptacle of the knee joint can lead to short-circuits in the electronics and malfunction of the knee joint. This can cause the patient to fall.

Replace the protective plug after the recharging (see fig. 3)!

Charging Time and Battery Capacity:

<i>Charging time</i>	<i>Capacity</i>
4.5 hours	45 hours of use (full battery)
1.5 hours	approx. 15 hours of use
20 minutes	approx. 5 hours of use

The indicated charging time and capacity values may vary, because, with increased capacity, charging time and patient use may increase.

Please note the following LED signals when using the 4E50-2 Battery Charger:

LED	Function
Green LED on	Battery Charger is ready for operation
Yellow LED on	Electronic knee joint is being charged
Yellow LED blinks	Battery at approx. 50% of full charge
Yellow LED off	Full battery, charging procedure over

6 Recommendations

You may leave the electronic knee joint connected to the battery charger even when the yellow LED is off. Overcharging or damage to the battery is not possible. Partial charging does not decrease battery shelf life (no memory effect). We recommend charging the product every day when used on a daily basis. Charging is only possible at temperatures over 0 °C (32 °F). If the temperature is lower than this, charging is interrupted to avoid battery damage.

Optional accessory – 4X74 Adapter for charging via 12-Volt auto/boat cigarette lighter receptacle.

7 Cleaning

Clean the product with a damp cloth (fresh water) when needed.

Dry the product with a lint-free cloth and allow it to air dry fully.

8 Disposal

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

9 Technical data

Battery charger

Reference number	4E50*
Storage and transport in original packaging	-25 °C/-13 °F to +70 °C/+158 °F
Storage and transport without packaging	-25 °C/-13 °F to +70 °C/+158 °F max. 93 % relative humidity, non-condensing
Operation	0 °C/+32 °F to +40 °C/+104 °F max. 93 % relative humidity, non-condensing
Input voltage	12 V
Service life	8 years

Power supply

Reference number	757L16*
Storage and transport in original packaging	-40 °C/-40 °F to +70 °C/+158 °F
Storage and transport without packaging	-40 °C/-40 °F to +70 °C/+158 °F 10 % to 93 % relative humidity, non-condensing
Operation	0 °C/+32 °F to +40 °C/+104 °F max. 90 % relative humidity, non-condensing
Input voltage	100 V~ to 240 V~
Mains frequency	50 Hz to 60 Hz
Output voltage	12 V

10 Symbols Used

Declaration of conformity according to the applicable European directives

Lot number (PPPP YYYY WW)

PPPP – plant

YYYY – year of manufacture

WW – week of manufacture

11 Liability

The manufacturer's warranty applies only if the product has been used under the conditions and for the purposes described. The manufacturer recommends that the product be used and maintained according to the instructions for use.

The 4E50-2 Battery Charger may be opened or repaired only by authorized Ottobock personnel.

12 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

13 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

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