




# ottobock.





## 10V40 MyoWrist 2Act




DE	Gebrauchsanweisung .....	7
EN	Instructions for use .....	17
FR	Instructions d'utilisation .....	26
IT	Istruzioni per l'uso .....	36
ES	Instrucciones de uso .....	45
PT	Manual de utilização .....	54
NL	Gebruiksaanwijzing .....	63
SV	Bruksanvisning .....	73
DA	Brugsanvisning .....	81
NO	Bruksanvisning .....	89
PL	Instrukcja użytkowania .....	97
HU	Használati utasítás .....	106
CS	Návod k použití .....	115
RO	Instrucțiuni de utilizare .....	124
TR	Kullanma talimatı .....	134
EL	Οδηγίες χρήσης .....	143
RU	Руководство по применению .....	153
JA	取扱説明書 .....	163
ZH	使用说明书 .....	171

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Deutsch

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---

## Explanation of symbols

**⚠ WARNING** Warnings regarding possible risks of severe accident or injury.

**⚠ CAUTION** Warnings regarding possible risks of accident or injury.

**NOTICE** Warnings regarding possible technical damage.

**INFORMATION** Additional information on the fitting/use.

### INFORMATION

Last update: 2016-03-18

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please contact the manufacturer if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- Please keep this document in a safe place.

The scope of delivery is shown on the front page.

## 1 Single components

### 1.1 Single components

For available single components refer to the catalogue.

### 1.2 Accessories

9S266 Chassis with Quick-Disconnect Wrist Unit

9E397=\* Coaxial Bushing

## 2 Description

### 2.1 Indications for use

The 10V40 MyoWrist 2Act is to be used **solely** for the exoprosthetic fitting of the upper extremities in combination with the 8E41 Ottobock Electric Hand.

## 2.2 Field of application

The 10V40 MyoWrist 2Act is used as a flexion wrist unit with a 8E41 Ottobock System Electric Hand with M12x1.5 Threaded Stud. The 10V40 MyoWrist 2Act supports the physiologically correct posture when working with the prosthesis and reduces the need for unnatural compensatory movements, especially when performing longer monotonous works.

## 2.3 Safety instructions

Please make your patients aware of the following safety instructions.

### WARNING

**Risk of accident when driving a motor vehicle.** An upper extremity amputee's ability to drive a vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities.

All persons are required to observe their country's national and state driving laws when operating motor vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorized test centre.

For maximum safety and convenience, Ottobock recommends that, at the very least, a specialist evaluate the need for any adaptations to the car (such as by installing a steering fork, automatic shift). Risk-free driving should be ensured even when the prosthesis is not functioning.

### WARNING

**Risk of accident when driving a motor vehicle.** Before operating a motor vehicle, the quick-disconnect wrist unit and flexion angle of the 10V40 MyoWrist 2Act must be so positioned that a slight turn of it, which can occur during steering, will not disconnect the hand from the prosthesis!

Switch the electric hand off. Driving with the electric hand turned ON may present a risk if the hand inadvertently moves due to unintentional muscle contraction or other causes.

**⚠ CAUTION**

**Risk of injury as a result of excessive strain.** The Ottobock prosthetic components described here have been developed for activities of daily living and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's safety!

Should prosthetic components be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the components inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Service Centre.

**NOTICE**

**Risk of corrosion.** Prosthetic components must not be exposed to corrosive elements such as fresh water, salt water and acids.

Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

**NOTICE**

**Damage caused by inappropriate environmental conditions.** Prosthetic components must not be subjected to intense smoke, dust, vibrations, shocks or high temperatures. Do not allow debris or liquids to get into the prosthetic components. Failure to follow these instructions can lead to malfunctioning and damage to the prosthesis.

**NOTICE**

**Damage caused by improper cleaning agents.** Do not use any aggressive cleaning agents, since they could cause damage to bearings, seals and plastic parts.

## INFORMATION

Please note that when carrying heavy objects with full flexion/extension, a torque moment can result causing the hand to rotate unintentionally.

Please be aware that the hand can flex unintentionally, if the flexion switch is pressed by accidentally.

### 2.4 Function

The 10V40 MyoWrist 2Act is suitable for hand sizes 7, 7¼, 7¾ and 8¼. It can be used on both sides. Individual flexion and extension of this wrist unit can be locked in 5 different positions (in increments of 20°).

The 10V40 MyoWrist 2Act is unlocked by pressing the release button (fig. 1). An integrated ratchet mechanism holds the hand until it is locked in the desired position and thus facilitates exact positioning. Flexion and extension are performed passively. Stop pressing the release button to relock the joint.

## 3 Technical data

<b>Weight [g]</b>	55
<b>Total height [mm]</b>	26
<b>Locking positions</b>	5 locking positions in increments of 20°
<b>Neutral or zero position</b>	0°
<b>Extension position</b>	-20° / -40°
<b>Flexion position</b>	20° / 40°
<b>Outside diameter [mm]</b>	36
<b>Maximum load capacity</b>	32.2 Nm $\hat{=}$ 27 kg

## 4 Handling

### 4.1 Assembly

Any version of the 8E41 Ottobock System Electric Hand is suitable for assembly.

To prepare the assembly of the 10V40 MyoWrist 2Act, the inner hand must be separated from the mechanism. For this purpose, use of the 711M1 Hand Holding Tool with 711M2 Mounting Plate is recommended (fig. 2).

To remove the hand cable of the 8E41 Ottobock System Electric Hand, first loosen the 9E83 Lock Nut using the 709Z7 Special Wrench (fig. 3).

#### **4.1.1 Disassembly of the hand cable**

##### *SensorHand Speed and MyoHand VariPlus Speed*

To loosen the hand cable in the mechanism, the control electronics must be removed.

1. Loosen the plug of the hand cable from the system electronics (Fig. 4).
2. Cut the hand cable in two (Fig. 5).
3. Pull the hand cable along with the 9E70 Cable Bushing (will be needed later) out of the threaded stud of the hand and dispose of the two cable residues.
4. Carefully remove the cover of the flat plugs (fig. 6).
5. Carefully loosen the plug connections (fig. 7).
6. Loosen the screw connections of the system electronics (fig. 8 + 9).
7. Remove the system electronics (fig. 10).

##### *DMC plus and Digital Twin System Electric Hands (except hand size 7)*

To loosen the hand cable in the mechanism, the control electronics must be removed.

1. Completely open the hand. Disassemble the 9E78 Bow Spring (fig. 11) using a screwdriver. Use a piece of leather, for example, to protect the mechanism from damage (fig. 11).
2. Unscrew the attachment screws of the control electronics (fig. 12).
3. Loosen the plug connections and remove the control electronics (fig. 13). Unsolder the hand cable from the plug using the 711M12 Soldering Jig (fig. 14). Pull out the exposed cable along with the 9E70 Cable Bushing (will be needed later) through the chassis on the side of the plug.

## *DMC plus and Digital Twin System Electric Hand – hand size 7*

To loosen the hand cable in the mechanism, the control electronics must be removed.

1. Loosen the attachment screw of the control electronics (fig. 15).
2. Carefully unplug the hand cable (fig. 16).
3. Apply a screwdriver behind the plug connector and turn with light pressure (fig. 17). This will press the control electronics out of the holder.
4. Unsolder the hand cable from the plug using the 711M12 Soldering Jig (fig. 14). Pull out the exposed cable along with the 9E70 Cable Bushing (will be needed later) through the chassis on the side of the plug.

### **4.1.2 Assembly of 10V40 MyoWrist 2Act**

#### **NOTICE**

**Risk of damage as a result of using an improper adhesive.** Using an unapproved adhesive can lead to malfunction/damage when torque moments are transferred and cause the hand to rotate unintentionally. To ensure a solid bond, **636K14 Loctite 601 Adhesive** is to be used.

#### **NOTICE**

**Risk of defects as a result of improperly prepared bonding surfaces.** To ensure a solid bond, the respective areas must be degreased with acetone before applying the adhesive.

Once you have degreased the areas, do not touch them anymore and immediately apply the adhesive to the surfaces after the acetone has dried.

#### **INFORMATION**

Follow all processing and safety instructions of the adhesive manufacturer.

Assembly of the 10V40 MyoWrist 2Act, of the 9S266 Chassis with Quick-Disconnect Wrist Unit and of the 9E397 Coaxial Bushing is the same with all 8E41 Ottobock System Electric Hands. The only difference for the dif-

ferent versions lies in the final connection of the hand cable. The assembly is here described on the example of a 8E41=8 SensorHand Speed.

1. To mount the 10V40 MyoWrist 2Act, pull the 9E70 Cable Bushing off the cable and shorten it on one side (fig. 18). The longer component will be used.
2. Degrease the threaded connectors of the 8E41 Ottobock System Electric Hand, of the 10V40 MyoWrist 2Act and of the 9S266 Chassis with Quick-Disconnect Wrist Unit with acetone, let them air dry and no longer touch them.
3. Lead the cable of the 9E397 Coaxial Bushing through the 9S266 Chassis with Quick-Disconnect Wrist Unit, the 10V40 MyoWrist 2Act, the modified 9E70 Cable Bushing and the 8E41 Ottobock System Electric Hand (fig. 19).
4. Apply a little bit of **636K14 Loctite 601 Adhesive** to the lower third of the thread of the 8E41 System Electric Hand (fig. 20). Take the 8E41 Ottobock System Electric Hand and screw in the 10V40 MyoWrist 2Act (fig. 21). To avoid twisting the cable, turn the 10V40 MyoWrist 2Act and not the System Electric Hand.

#### INFORMATION

Please be aware that the cable can become twisted in this process. The cable must run straight through the wrist unit in order to prevent it from being damaged.

5. Fully screw the 10V40 MyoWrist 2Act in. Turn back by maximally 1 turn to position the 10V40 MyoWrist 2Act. For good accessibility of the release switch, it is recommended to position the switch anteriorly, starting from the neutral or zero position (fig. 22).
6. Apply a little bit of **636K14 Loctite 601 Adhesive** to the thread of the 9S266 Chassis with Quick-Disconnect Wrist Unit and screw it into the 10V40 MyoWrist 2Act hand-tight (fig. 23).
7. Pull through the 9E397 Coaxial Bushing upward. Position the 9E397 Coaxial Bushing in the guide holes of the 9S266 Chassis with Quick-Disconnect Wrist Unit and tighten (fig. 24).
8. After fastening the 9E397 Coaxial Bushing, maximally flex the 10V40 MyoWrist 2Act once in both directions (fig. 25).

### 4.1.3 Connection of the hand cable

#### INFORMATION

Make sure that the cable is not stretched tight, since the flexion movement of the wrist will require a certain amount of cable length. For testing, maximal extension and flexion movements must be performed with the System Electric Hand.

#### *SensorHand Speed and MyoHand VariPlus Speed*

1. Shorten the cable to approx. 3 cm after it has exited the bore hole of the tube chassis of the 8E41=8 SensorHand Speed (fig. 26). Separate the ends of the individual strands of the cable, remove 2 mm of the insulation and tin (fig. 27).
2. Provide the strands with shrink tubing and solder them to the new, supplied plug (fig. 28). Start with the red strand on both ends.
3. Slide the shrink tubing over the solder joints and shrink.
4. Mount the control electronics and attach the plug. Perform a function check. If necessary, reposition the 10V40 MyoWrist 2Act in relation to the System Electric Hand.
5. Make sure to properly place the cable (fig. 29 + 30). If necessary, reposition the 10V40 MyoWrist 2Act once again in relation to the 8E41 Ottobock System Electric Hand (fig. 32).

#### *DMC plus and Digital Twin System Hands*

1. Shorten the cable to approx. 4 cm after it has exited the bore hole of the tube chassis of the 8E41 Ottobock System Electric Hand (fig. 26). Separate the ends of the individual strands of the cable, remove 2 mm of the insulation and tin (fig. 27).
2. Use the 711M12 Soldering Jig to solder on the plug. The marking in Fig. 31 shows the position of the red strand of the hand cable.
3. Reassemble the control electronics with the screws and attach the plug (fig. 13-11). Reassemble the 9E78 Bow Spring. Perform a function check. If necessary, reposition the 10V40 MyoWrist 2Act once again in relation to the 8E41 Ottobock System Electric Hand (fig. 32).

## *DMC plus and Digital Twin System Hands, hand size 7*

1. Shorten the cable to approx. 5 cm after it has exited the bore hole of the tube chassis of the System Hand (fig. 26). Separate the ends of the individual strands of the cable, remove 2 mm of the insulation and tin (fig. 27).
2. Use the 711M12 Soldering Jig to solder on the plug. The marking in Fig. 31 shows the position of the red strand of the hand cable.
3. Press in the control electronics, tighten the screws (fig. 17-15) and attach the plug. Perform a function check. If necessary, reposition the 10V40 MyoWrist 2Act once again in relation to the 8E41 Ottobock System Electric Hand (fig. 32).

## **5 Legal information**

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### **5.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### **5.2 CE conformity**

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH