







## 12R6=\*, 12R7=\*, 13R1=\*

<b>DE</b>	Gebrauchsanweisung .....	4
<b>EN</b>	Instructions for use .....	12
<b>FR</b>	Instructions d'utilisation .....	20
<b>IT</b>	Istruzioni per l'uso .....	29
<b>ES</b>	Instrucciones de uso .....	38
<b>PT</b>	Manual de utilização .....	47
<b>NL</b>	Gebruiksaanwijzing .....	56
<b>SV</b>	Bruksanvisning .....	65
<b>DA</b>	Brugsanvisning .....	73
<b>NO</b>	Bruksanvisning .....	81
<b>FI</b>	Käyttöohje .....	89
<b>PL</b>	Instrukcja użytkowania .....	97
<b>HU</b>	Használati utasítás .....	106
<b>CS</b>	Návod k použití .....	115
<b>RO</b>	Instrucțiuni de utilizare .....	123
<b>HR</b>	Upute za uporabu .....	131
<b>SL</b>	Navodila za uporabo .....	139
<b>SK</b>	Návod na používanie .....	147
<b>BG</b>	Инструкция за употреба .....	156
<b>TR</b>	Kullanma talimatı .....	165
<b>EL</b>	Οδηγίες χρήσης .....	173
<b>RU</b>	Руководство по применению .....	182
<b>JA</b>	取扱説明書 .....	192
<b>ZH</b>	使用说明书 .....	200
<b>KO</b>	사용 설명서 .....	206

## INFORMATION

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### INFORMATION

Date of last update: 2020-04-08

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.

► Please keep this document for your records.

## 1.1 Construction and Function

The 12R6=\* and 12R7=\* modular kits and the 13R1=\* elbow joint modular kit are used for the prosthetic fitting of patients with amputations from the level of the upper arm to the shoulder girdle. The focus is on a lightweight prosthetic fitting. The modular kits contain structural components for the alignment of an arm prosthesis. The modular elbow joint is included in the modular kits. It can be adjusted to 7 locking positions.

Reference number	Designation	Fitting
12R6=*	Transhumeral modular kit	Upper arm amputation
12R7=*	Shoulder disarticulation modular kit	Shoulder girdle amputation, shoulder joint disarticulation
13R1=*	Modular elbow joint	Component of the modular prosthesis

## 1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

## 2 Intended use

### 2.1 Indications for use

The product is to be used solely for upper limb prosthetic fittings.

### 2.2 Area of application

The product serves to create passive arm prostheses.

### 2.3 Environmental conditions

Allowable environmental conditions
Temperature range for use: -10 °C to +45 °C
Relative humidity: max. 90%, non-condensing
Moisture: fresh water and salt water as splashing water, sweat
Unallowable environmental conditions
Mechanical vibrations or impacts
Acids
Dust, sand, highly hygroscopic particles (e.g. talcum), smoke

## 3 Safety

### 3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 3.2 General safety instructions

 **CAUTION**

#### **Excessive strain on the product**

Risk of injury due to breakage of load-bearing components

- ▶ Use the product according to the specified area of application (see Page 13).

 **CAUTION**

#### **Unallowable combination of prosthetic components**

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

 **CAUTION**

#### **Mechanical damage to the product**

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

 **CAUTION**

#### **Reuse on another patient**

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Only use the product for a single patient.

**NOTICE****Use under unallowable environmental conditions**

Damage to product due to unallowable environmental conditions

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop, etc.).

**Signs of changes in or loss of functionality during use**

Instability or lack of joint function indicate a loss of function.

**4 Scope of delivery**

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

12R6=*				
Fig.	Item	Quantity	Designation	Reference number
–	–	1(■)	Instructions for use	647G471
–	–	1(■)	Modular kit (Modular arm components, modular elbow joint)	–
–	–	1(■)	Adapter	13R6

12R7=*				
Fig.	Item	Quantity	Designation	Reference number
–	–	1(■)	Instructions for use	647G471
–	–	1(■)	Modular kit (Modular arm components, modular elbow joint)	–
–	–	1(■)	Adapter	13R7

13R1=*				
Fig.	Item	Quantity	Designation	Reference number
–	–	1(■)	Instructions for use	647G471
–	–	1(■)	Modular elbow joint	–

## 5 Preparing the product for use

### CAUTION

#### **Incorrect alignment or assembly**

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

### 5.1 Determining length and position

The position and length of the prosthesis are determined on the patient. The prosthesis is assembled temporarily for this purpose, and the patient tries it on. During this trial fitting, the dimensions of the prosthesis are adjusted and then corrected.

Length range available		
Reference number	Forearm range [mm]	Upper arm range [mm]
12R6=*	90 to 294	67 (with 13R6), 75 to 170 (without 13R6)
12R7=*	90 to 294	59 (with 13R7), 57 to 296 (without 13R7)

- 1) Temporarily assemble the prosthesis (prosthetic shoulder joint/prosthetic socket up to terminal device).
- 2) Perform the trial fitting on the patient.

#### **Unilateral fitting**

- ▶ Determine the dimensions of the prosthesis on the basis of the dimensions of the contralateral side. Use the pivot point of the elbow joint and the distal end of the middle finger as a reference.

#### **Bilateral fitting**

- ▶ Determine the dimensions of the prostheses to achieve the maximum function and flexibility for the patient.

### 5.2 Upper arm fitting

The 12R6=\* modular kit is used for prosthetic upper arm fittings of patients with various residual limb lengths. The 13R6 adapter can be attached for long residual limbs.

#### **Optional: replacing the adapter**

##### > **Tools required:**

Torque wrench (adjustable to **0.5 Nm**)

- 1) Unscrew the set screw that secures the proximal tube to the modular elbow joint. Pull the tube off the modular elbow joint.

- 2) Unscrew the cap screw of the adapter on the proximal end of the tube and remove the adapter. The tube is no longer needed.
- 3) Unscrew the 3 truss head screws of the adapter and remove the lamination ring. The adapter is no longer needed.
- 4) Place the new adapter on the lamination ring and screw in the 3 truss head screws (torque: **4 Nm**).
- 5) Fit the new adapter completely on the modular elbow joint and tighten the set screw (torque **0.5 Nm**).

### **Trial fitting**

#### > **Tools required:**

Torque wrench (adjustable to **0.5 Nm**)

- > The inner socket and moulded hard foam part for the outer socket have been prepared.
- 1) Assemble the inner socket and the moulded hard foam part.
- 2) Connect all prosthesis components to the moulded hard foam part.
- 3) **Optional:** adjust the tube length to the patient's measurements (see Page 19). Adjust the length of the moulded hard foam part to the patient's measurements.
- 4) **If the measurements are correct:** mark the position of the lamination ring on the moulded hard foam part.

### **Definitive mounting**

#### > **Tools required:**

Torque wrench (adjustable to **0.5 Nm**)

- > The inner socket and moulded hard foam part for the outer socket have been prepared.
- 1) Unscrew the 3 truss head screws on the proximal end of the arm component and remove the lamination ring.
- 2) Insulate the inside of the lamination ring.
- 3) Transfer the inner contour of the lamination ring to the marked position on the moulded hard foam part with lightweight putty.
- 4) Remove the lamination ring after the lightweight putty hardens.
- 5) Grind the shape of the moulded hard foam part.
- 6) Pull on a PVA bag and tie off distally.
- 7) Position the lamination ring on the PVA bag.
- 8) Laminate the outer socket. In the process, secure the connection to the lamination ring with carbon fibre cloth.
- 9) Grind away or cut away the distal surface of the lamination ring after the laminating resin hardens. Remove the moulded hard foam part from the outer socket.

- 10) Attach the adapter of the arm component to the lamination ring.  
**INFORMATION: Check the dimensions of the definitive prosthesis with a trial fitting and correct them if necessary.**

### **Adjusting the humeral rotation**

- > **Tools required:**
  - Torque wrench (adjustable to 0.5 Nm)
- ▶ Adjust the torque of the 3 truss head screws to achieve the desired resistance for the humeral rotation (torque **1.5 Nm to 4 Nm**).

### **5.3 Shoulder fitting**

The 12R7=\* modular kit is used for the prosthetic fitting of patients with a shoulder girdle amputation or shoulder disarticulation. The attached tube is for connecting to the 12S7 prosthetic shoulder joint. The 13R7 adapter is used for the connection to the 12S4 prosthetic shoulder joint.

#### **Optional: replacing the adapter**

- > **Tools required:**
  - Torque wrench (adjustable to **0.5 Nm**)
- 1) Unscrew the set screw that secures the proximal tube to the modular elbow joint.
- 2) Pull the tube off the modular elbow joint; it is no longer needed.
- 3) Fit the adapter completely on the modular elbow joint and tighten the set screw (torque **0.5 Nm**).

#### **Trial fitting**

- > **Required tools and materials:**
  - Torque wrench (adjustable to **0.5 Nm**)
- 1) Connect all prosthesis components to each other.  
**INFORMATION: The 10R2\* adapter is needed for assembly of the 12S7 prosthetic shoulder joint (not included in scope of delivery).**
- 2) **Optional:** adjust the tube length to the patient's measurements (see Page 19). **Only for 12S4:** adjust the length of the upper arm bar to the patient's measurements.

#### **Definitive mounting**

- > **Required tools and materials:**
  - Torque wrench (adjustable to **0.5 Nm**), 636K13 Loctite® 241
- 1) **Only for 12S4:** Unscrew the truss head screws of the adapter. Apply Loctite to the truss head screws and tighten them (torque: **4 Nm**).
- 2) **Only for 12S7:** fit the 10R2\* adapter completely on the tube and tighten the cap screw (torque: **4 Nm**).

- 3) Connect all prosthesis components to each other.

**INFORMATION: Check the dimensions of the definitive prosthesis with a trial fitting and correct them if necessary.**

#### 5.4 Optional: shortening the tube

##### CAUTION

##### **Incorrect processing of tube**

Fall due to damage to the tube

- ▶ Do not clamp the tube into a vice.
- ▶ To shorten the tube, use only a tube cutter or a cutting device.

- > **Required materials:** 710D4 torque wrench, 719R3 tube cutter or 704Y14\* cutting device, 718R1 tube deburrer, degreasing cleaning agent (e.g. 634A58 isopropyl alcohol)
- 1) Shorten the tube according to the patient's measurements.
- 2) Deburr the inside and outside of the cut edge with the tube deburrer.

#### 5.5 Modular elbow joint

The 13R1=\* modular elbow joint has 7 locking adjustment steps. It covers a flexion range from 0° to 141°. The first locking position is at 15°, and the other positions are adjusted in 21° increments. The modular elbow joint is available in left and right versions. This makes it possible to always position the lock lever medially.

- > **Tools required:**
  - Torque wrench (adjustable to **0.5 Nm**)
- 1) Fit the connection components on the modular elbow joint.
- 2) **INFORMATION: Position the lock lever medially.**
  - Tighten the set screws (torque: **0.5 Nm**).

#### 6 Cleaning

- 1) Clean the product with a damp, soft cloth.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

#### 7 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

## 8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## 9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

### 9.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

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