
4WR95=1, 4WR95=2

DE Gebrauchsanweisung	3
EN Instructions for use	8
FR Instructions d'utilisation.....	14
IT Istruzioni per l'uso	20
ES Instrucciones de uso	26
PT Manual de utilização.....	32
NL Gebruiksaanwijzing.....	38
SV Bruksanvisning	44
DA Brugsanvisning	49
NO Bruksanvisning	54
FI Käyttöohje.....	60
PL Instrukcja użytkowania.....	65
HU Használati utasítás	71

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INFORMATION

Last update: 2015-10-08

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

1.1 Construction and Function

The 4WR95=1 and 4WR95=2 lamination anchors are laminated into a prosthetic socket. They serve to connect with the distal prosthetic components. The 4WR95=1 lamination anchor is equipped with a pyramid receiver and a more angled anchor arm. The 4WR95=2 lamination anchor features a pyramid adapter.

The product is intended for the fabrication of a waterproof modular prosthesis.

1.2 Combination possibilities

This prosthetic component is part of the Ottobock modular system. As a result of its design, it may be combined with other prosthetic components of the modular system. You can find possible restrictions in this section. Contact the manufacturer with questions.

- The product may only be combined with Ottobock prosthetic components that are approved for waterproof additional prostheses.

2 Application

2.1 Indications for use

The product is to be used solely for lower limb prosthetic fittings.

2.2 Area of application

Approved for a body weight of **up to 150 kg**.

- The product may only be used for waterproof **additional prostheses**.

The set screws of stainless steel can optionally be replaced by titanium set screws, which increase the protection against corrosion.

4WR95=1

The product must be used in TF prostheses only.

2.3 Environmental conditions

Allowable environmental conditions

Temperature range for use: -10 °C to +60 °C

Storage and transport conditions: -10° C to +60° C, relative humidity: no restrictions

Humidity: fresh water, salt water – cleaning and drying is required after contact and use.

Allowable immersion depth: 1 m

Unallowable environmental conditions

Sauna, extreme watersports

Mechanical vibrations or impacts

Unallowable environmental conditions

Perspiration, urine, acids or bases

Dust, sand, highly hygroscopic particles (e.g. talcum)

2.4 Service Life

In general, all modular adapters are tested by the manufacturer for 3 million load cycles. Depending on the patient's level of activity, this corresponds to a service life of 3 to 5 years.

3 Safety

3.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION

Exceeding the service life

Risk of injury due to change in or loss of functionality and damage to the product

- ▶ Ensure that the approved service life is not exceeded (see section "Service life").

CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).

- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

⚠ CAUTION

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved according to the section "Combination possibilities".
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

⚠ CAUTION

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

- ▶ Use the prosthetic components according to the classification (see section "Area of application").

NOTICE

Use under unallowable environmental conditions

Damage to product due to unallowable environmental conditions

- ▶ Do not expose the product to unallowable environmental conditions (see the section "Environmental Conditions").
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of delivery

The following single components and accessories are included in the scope of delivery in the specified quantities and are available for reordering as single components (■), single components with minimum order quantities (▲) or single-component packs (●):

4WR95=1, 4WR95=2				
Fig.	Item	Quantity	Designation	Reference number
-	-	1(■)	Instructions for use	-
-	-	1(■)	Lamination anchor	-
-	-	1(■)	Cap screw	501T24=M5x25
-	-	1(■)	Rounded washer	507U16=5.2-NI-RO
For 4WR95=1				
-	-	4(▲)	Set screw	506G3=M8x12-"-NIRO"

5 Preparation for use

⚠ CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

⚠ CAUTION

Improper assembly of the screw connections

Fall due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified installation torque values.
- ▶ Observe the instructions for securing the screw connections and the use of the correct length.

CAUTION

Wrong processing of the anchor arms of the lamination anchor

Fall due to breakage of the anchor arms

- ▶ Only bend anchor arms made of steel.
- ▶ Avoid excessive or too frequent bending of the anchor arms.
- ▶ Use the 711S4* bending irons for bending.

INFORMATION

- ▶ Only use the titanium cap screw at the end of the finishing process for the prosthesis.
- ▶ Apply Loctite® to the titanium cap screw.

5.1 Mounting of the screwed insert

> Required tools and materials:

710D4 torque wrench, cap screw

- 1) Screw in the screwed insert completely into the thread of the lamination anchor.
- 2) Tighten the cap screw from the threaded clamping using the torque wrench (mounting torque: **10 Nm**).

5.2 Laminating a Transfemoral Socket

> Required tools and materials:

711S4* bending iron, 99B81 PVA bag, 623T3 Perlon stockinette, 616B1 carbon fibre cloth strap, 616G6 Dacron felt, 616G12 carbon fibre cloth, 616G15 carbon fibre stockinette, 617H119 Orthocryl lamination resin 80:20 PRO

- 1) Mount a screwed insert according to section "Mounting of the screwed insert" (see Page 11).
- 2) Pull a soaked PVA bag over the model.
- 3) Apply one layer of Dacron felt and thin out the fibres along the edge.
- 4) Pull two layers of Perlon stockinette over the model.
- 5) Place three layers of carbon fibre cloth (e.g. **15 cm x15 cm**), off-set on the distal end of the model.
- 6) Pull two layers of Perlon stockinette over the model.
- 7) Pull a soaked PVA bag over the model.
- 8) Perform the casting process with Orthocryl.
- 9) Once the lamination resin has cured, remove the PVA bag.
- 10) **INFORMATION: Align an anchor arm that has been angled at the factory in the posterior direction.**
Bring the lamination anchor arms in the correct **a - p** and **m - l** position. Never align the threaded clamping in the anterior or posterior direction.
- 11) Place the lamination anchor at the socket end as required for proper alignment and set.
- 12) **Optional:** Adapt the anchor arms of the lamination anchor to the model with a bending iron.
- 13) Apply the lamination dummy according to the instructions in the section "Installing the lamination dummy" (see Page 12).
- 14) Stabilise and fix the lamination anchor with a sealing resin/talcum powder mixture.
- 15) Pull one layer of Perlon stockinette over the model.

- 16) Place two layers of carbon fibre cloth (e.g. **15 cm x 15 cm**) offset over the lamination anchor arms.
- 17) Pull one layer of Perlon stockinette over the model.
- 18) Pull two layers of carbon fibre stockinette over the model.
- 19) Pull two layers of Perlon stockinette over the model and tie off below the cap screw of the lamination anchor.
- 20) Pull a soaked PVA bag over the model.
- 21) Perform the casting process with Orthocryl.
- 22) After the lamination resin has cured, remove the lamination dummy.

5.3 Laminating a Transtibial Socket

> Required tools and materials:

711S4* bending iron, 99B81 PVA bag, 623T3 Perlon stockinette, 616B1 carbon fibre cloth strap, 699B1 fibreglass roving, 616G12 carbon fibre cloth, 616G15 carbon fibre stockinette, 617H119 Orthocryl lamination resin 80:20 PRO

- 1) Mount a screwed insert according to section "Mounting of the screwed insert" (see Page 11).
- 2) Pull a soaked PVA bag over the model. The PVA bag should be sealed by heat welding at the distal end.
- 3) Pull two layers of Perlon stockinette over the model.
- 4) Place three layers of carbon fibre cloth (e.g. **15 cm x15 cm**), offset on the distal end of the model.
- 5) Pull two layers of Perlon stockinette over the model.
- 6) Bring the lamination anchor arms in the correct **a - p** and **m - l** position. Never align the threaded clamping in the anterior or posterior direction.
- 7) Place the lamination anchor at the socket end as required for proper alignment and set.
- 8) **Optional:** Adapt the anchor arms of the lamination anchor to the model with a bending iron.
- 9) Apply the lamination dummy according to the instructions in the section "Installing the lamination dummy" (see Page 12).

- 10) Stabilise and fix the lamination anchor with a sealing resin/talcum powder mixture.
- 11) Pull the fibreglass roving through the holes in the anchor arms and allow it to hang in a loop. If there are no holes, allow the fibreglass roving to hang around the anchor arms in a loop.
- 12) Place three layers of carbon fibre cloth (e.g. **15 cm x 15 cm**) offset over the lamination anchor arms.
- 13) Pull one layer of Perlon stockinette over the model.
- 14) Pull two layers of carbon fibre stockinette over the model.
- 15) Pull two layers of Perlon stockinette over the model and tie off below the cap screw of the lamination anchor.
- 16) Pull on a soaked PVA bag.
- 17) Perform the casting process with Orthocryl.
- 18) After the lamination resin has cured, remove the lamination dummy.

5.4 Installing the lamination dummy

> Required tools and materials:

4X46 lamination dummy

- 1) Tighten the cap screw of the threaded clamping to a point where the screwed insert can still be screwed in and out.
- 2) Unscrew the screwed insert from the lamination anchor.
- 3) Screw in the lamination dummy.

5.5 Installing the distal prosthetic components

Mounting of a screwed insert

The threaded connector is used for installing a connecting component with an external thread (e.g. pyramid adapter or prosthetic knee joint).

> Required tools: 710D4 torque wrench, 636K13 Loctite®

- 1) Screw in the screwed insert (e.g. pyramid adapter) completely into the thread of the lamination anchor.
- 2) The screwed insert may be unscrewed by a **maximum of 1/4 turn** out of the thread in order to align the insert.

- 3) Secure the cap screw of the threaded clamping using Loctite® and tighten using the torque wrench (mounting torque: **10 Nm**).

Connecting the pyramid adapter and pyramid receiver

The pyramid adapter is fixed with the set screws of the pyramid receiver.

> Required tools and materials:

710D4 torque wrench, 636K13 Loctite®

1) Trial fitting:

Screw in the set screws.

Use the torque wrench to tighten the set screws (**10 Nm**).

2) Definitive mounting:

Use Loctite® to secure the set screws.

Screw in the set screws.

Pre-tighten the set screws with the torque wrench (**10 Nm**) and then tighten them (**15 Nm**).

- 3) Replace any set screws that are protruding or recessed too much with suitable ones (see selection table).

Selection table for set screws		
Reference number	Reference number (titanium)	Length (mm)
506G3=M8X12-"Niro"	506G5=M8X12	12
506G3=M8X14-"Niro"	506G5=M8X14	14
506G3=M8X16-"Niro"	506G5=M8X16	16

Alignment

The set screws in the pyramid receiver can be used to make static adjustments during alignment, trial fittings and after the prosthesis is finished.

Replacement and disassembly

The set position of the prosthetic component can be maintained during replacement or disassembly. In order to do this, unscrew the two set screws that are screwed in the furthest and located next to each other.

6 Cleaning

- 1) Rinse the product with clear fresh water.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

7 Maintenance

⚠ CAUTION

Failure to follow the maintenance instructions

Risk of injuries due to changes in or loss of functionality and damage to the product

▶ Observe the following maintenance instructions.

- ▶ The prosthetic components should be inspected after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

8 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage

caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

10 Technical data

Reference number	4WR95=1	4WR95=2
Weight [g]	165	165
System height [mm]	44	2
Material	Stainless steel	

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Otto Bock HealthCare GmbH