

ottobock.

2R10, 2R33, 2R51

DE	Gebrauchsanweisung	4
EN	Instructions for use	7
FR	Instructions d'utilisation	10
IT	Istruzioni per l'uso	14
ES	Instrucciones de uso	17
PT	Manual de utilização.....	21
NL	Gebruiksaanwijzing.....	24
SV	Bruksanvisning.....	28
DA	Brugsanvisning.....	31
NO	Bruksanvisning.....	34
FI	Käyttöohje.....	38
PL	Instrukcja użytkowania.....	41
HU	Használati utasítás	44
CS	Návod k použití.....	48
RO	Instrucțiuni de utilizare.....	51
TR	Kullanma talimatı.....	55
EL	Οδηγίες χρήσης	58
RU	Руководство по применению	62
JA	取扱説明書.....	66


2R33 2R51 2R10


①

②

2S88=22-23

2S88=24-25

2S88=26-27

3

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English

Last update: 2015-06-10

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

Note: The scope of delivery is shown on the front page.

1 Single components (Fig. 1)

1.1 Single component ■

–

1.2 Single-component pack ●

- 2D5** consisting of:
- 1 lower joint bushing (1)
 - 1 plantar flexion bumper, hard (2)
 - 1 plantar flexion bumper, medium (2)
 - 1 plantar flexion bumper, soft (2)
 - 1 shell

- 1 washer (3)
- 2 hexagon nuts (4)

1.3 Minimum order quantity ▲

–

2 Description

2.1 Indications for use

The structural parts shown on the front page are to be **solely** used for the prosthetic fitting of the lower limb.

2.2 Field of application

Field of application according to the **Ottobock MOBIS Mobility System:**

Approved for a **patient weight of up to 100 kg (220 lbs)**

Attention!

Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water, and acids. Using this medical product under such conditions will render all claims against Otto Bock HealthCare null and void.

Please also inform your patients!

2.3 Function

The Single Axis Foot Adapters 2R33, 2R51 and 2R10 are used for connection of the Single-Axis Feet 1G9 and 1H* with an adapter pyramid receiver.

3 Technical data

	2R33=*		2R51=*		2R10=*	
Size range	22-25 cm	26-30 cm	22-25 cm	26-27 cm	22-25 cm	26-30 cm
Weight	200 g	210 g	230 g	235 g	325 g	340 g
System height	already considered in the foot					
Material	Titanium	Titanium	Aluminum	Aluminum	Steel	Steel

Patient weight max.	100 kg		100 kg		100 kg	
Mobility grade	–	–	–	–	–	–

4 Handling

4.1 Assembly

In comparison to the 2R10 and 2R33 adapters, the dorsal stop is changed with 2R51 so as to enable an elastic dorsal movement. A Pedilan stop from the 2S88 Dorsal Stop Set should be used for assembly of the adapter with single-axis foot with foam cover (1H40, 1H38) (Fig. 2).

The following dorsal stops are available (Fig. 3):

2S88=22–23 Dorsal Stop Set –Pedilan–, 2-part, 1 pcs. each soft and hard, for 2R51=22–25 in conjunction with the single-axis foot with foam cover, size 22 and 23 cm

2S88=24–25 idem, for 2R51=22–25 in conjunction with the single-axis foot with foam cover, size 24 and 25 cm

2S88=26–27 idem, for 2R51=26-27 in conjunction with the single-axis foot with foam cover, size 26 and 27 cm

The dorsal set stop must be ordered separately as an option.

The initial sluggishness of the joint has no adverse effect on the function.

4.2 Adjustments and final assembly

Attention!

Hex nut assembly torque: **6 Nm**

4.3 Maintenance instructions

Note:

As a basic principle, all Ottobock modular adapters are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

We recommend carrying out regular safety checks once a year.

5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

5.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

5.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

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Otto Bock HealthCare GmbH