

# ottobock.

## 17B95/17B96

<b>DE</b>	Gebrauchsanweisung .....	3
<b>EN</b>	Instructions for Use .....	6
<b>FR</b>	Instructions d'utilisation .....	9
<b>IT</b>	Istruzioni per l'uso .....	12
<b>ES</b>	Instrucciones de uso .....	15
<b>PT</b>	Manual de utilização.....	18
<b>NL</b>	Gebruiksaanwijzing.....	21
<b>SV</b>	Bruksanvisning .....	24
<b>DA</b>	Brugsanvisning.....	26
<b>NO</b>	Bruksanvisning .....	29
<b>PL</b>	Instrukcja użytkowania .....	32
<b>CS</b>	Návod k použití.....	35
<b>TR</b>	Kullanma talimatı.....	38
<b>EL</b>	Οδηγίες χρήσης .....	40
<b>RU</b>	Руководство по применению .....	43

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Date of the last update: 2013-03-26

- Please read this document carefully.
- Follow the safety instructions.

## 1 Application

The 17B95/17B96 System Knee Joints with Wedge Lock and Remote Cable-Release must be used in pairs **solely** for the orthotic fitting of the lower limb.

### Attention!

The 17B95=16 System Knee Joints with 16 mm connection width can be used only with a loading up to 50 kg (110 lbs.)

### Note

These System Knee Joints are designed to be used together with free-motion or dorsiflexion assist ankle joints. Please note: using an ankle joint with a dorsal stop significantly increases the wear on the Ottobock Orthotic Components.

## 2 Indication

Partial paralysis or complete paralysis of the lower limb musculature. Specific applications are to be determined by the physician.

## 3 Construction

### 3.1 Components (fig. 1)

- |                  |                     |
|------------------|---------------------|
| 1) 30X62         | Upper Joint Section |
| 2) 30U90         | Lower Joint Section |
| 3) 513D19=3.8x14 | Compression Spring  |
| 4) 17Y88         | Wedge Lock          |
| 5) 30Z22         | Locking Rocker      |
| 6) 17Y93=*       | Cotter Bolt         |
| and 501S32=*     | Joint Screw         |

The System Knee Joints with Remote Cable-Release must be used with the 513D19=3.8x14 Compression Spring which holds the 17Y88 Wedge Lock against the 30U90 Lower Joint Section. The 17Y88 Wedge Lock is released via the 30Z22 Lock Rocker by pulling on the 21A18=2 Perlon Cable.

## 4 Processing

Please consider the higher structural loading on the Ottobock System Joints due to different clinical parameters, such as forefoot lever length or angular limb deformities which must be accommodated in the orthotic alignment.

Insert the 17B38, 17B39, 17B40, 17B41, 17B84, 17B85 System Side Bars or the 17Y105, 17Y107 System Lamination Bars as connection elements into the Knee Joint connection channels. The use of 636W28 Ottobock Special Adhesive is required to stabilize the bars in the connection channels.

To do this, degrease the attachment surfaces using 634A1 Ottobock Thinner and Solvent. Mix the contents of the two tubes containing special adhesive and hardener (included), and apply to the joint's connection channels. Insert the System Bars and fasten in the 501T7=7.5x9xM5 Phillips Head Screws. Allow three to four hours for the adhesive to dry prior to use; final bonding strength will be achieved after about 16 hours.

The 513D19=3.8x14 Compression Spring (3) must always be installed in the **upper joint section** (1), even during trial fittings.

The Ottobock System of Joint Bars and System Joints are made of steel that is both rust-resistant and acid-resistant. We recommend final buffing and polishing for increased protection against corrosion and 618T40 Sintering Powder for coating.

### Note

Do not apply heat in excess of 300°C (572° F), or the material's temper will be destroyed.

Hardened joint parts should not be heat-treated.

### Attention!

Before delivering the orthosis to the patient, make sure that both locks are parallel and engaging safely. Instruct the patient on the proper handling of the orthosis.

## 5 Maintenance

Use only 633F7 Special Lubricant. To ensure proper function, make sure no dirt particles enter the joint mechanism. Use appropriate protection devices (e.g. joint covers), as necessary.

### **Attention!**

Visible wear of the lock mechanism may affect the patient's safety. Replace the wedge lock and/or the lower joint section, as necessary.

### **Note**

It is recommended that the function and condition of this orthosis be checked every six months. The 513D19=3.8x16 Compression Springs should be replaced at this time, or at the very latest once each year.

## 6 Liability

The manufacturer's warranty applies only if the product has been used under the conditions and for the purposes described. The manufacturer recommends that the product be used and maintained according to the instructions for use.

## 7 CE conformity

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH