

1E91 Runner / 1E93 Runner Junior

1E91 Runner / 1E93 Runner Junior

The Runner and Runner Junior are characterized by a durable and lightweight carbon spring providing the athletic adult or child with a powerful drive and stable turning characteristics. The innovative sliding and rotatable pyramid allows for adjustment of the carbon spring dynamics and stiffness of the Runner and Runner Junior feet. The feet are adaptable according to individual requirements. This has no effect on the basic alignment.

1 Pyramid adapter, rotatable

The adapter's pyramid is continuously rotatable to facilitate the individual adjustment of the prosthesis. For the 1E91 Runner, the 4R216 four-hole adapter is an alternative option (see fig. 2)

- Anterior positioning of the adapter for relaxed running and comfort due to increased spring deflection. (see fig. 3)
- Posterior positioning of the adapter for powerful running due to an increased forward propulsion. (see fig. 4)

5 All-terrain sole and spike pad

Soles for the Runner include 2Z540 (all-terrain) and 2Z541 (spike). Sole for the Runner Junior is 2Z543 (all-terrain).



©2016 Otto Bock HealthCare LP · 13590B 6/16

Stiffness chart

| Body weight | Stiffness version | |
|--------------------------------|------------------------------|---------------|
| 15 to 20 kg / 33 to 44 lbs | SPR-1 | - |
| 20 to 25 kg / 44 to 55 lbs | SPR-2 | - |
| 25 to 30 kg / 55 to 66 lbs | SPR-3 | - |
| 30 to 37 kg / 66 to 81 lbs | SPR-4 | - |
| 37 to 45 kg / 81 to 99 lbs | SPR-5 | - |
| | Long-distance running | Sprint |
| 40 to 50 kg / 88 to 110 lbs | SPR-1 | SPR-2 |
| 50 to 60 kg / 110 to 132 lbs | SPR-2 | SPR-3 |
| 60 to 72 kg / 132 to 158 lbs | SPR-3 | SPR-4 |
| 72 to 86 kg / 158 to 189 lbs | SPR-4 | SPR-5 |
| 86 to 104 kg / 189 to 229 lbs | SPR-5 | SPR-6 |
| 104 to 125 kg / 229 to 275 lbs | SPR-6 | - |

1E93

1E91



Technical data

| | |
|---|--|
| Max. body weight | Runner Junior – 99 lbs / 45 kg Runner – 275 lbs / 125 kg |
| Weight | 6 – 7.4 oz |
| Structural Height (Under standing load) | Runner Junior – 6 1/4 in (16 cm) Runner – 10 5/8 in (27 cm) |

Ordering example

| Article number | = | Stiffness | - | S - N |
|----------------|---|-----------|---|-------|
| 1E91 | = | SPR-2 | - | S - N |
| 1E93 | = | SPR-5 | - | S - N |

* Note - Proximal adapter and sole must be ordered separately. Not included with carbon foot spring.

Uses



Running

1E91 Runner 1E93 Runner Junior

Order Form



Shipping Options:

- Next Day Ground
 2-Day Other _____

Account Information

Date _____ Account Number _____

Bill To _____

Phone Number _____ Fax Number _____

Email Address _____ Ordered By _____

P.O. Number _____

Ship To

Name _____ Phone Number _____

Address _____ City _____ State/Zip Code _____

1 4R111
Lamination Anchor with
Pyramid Receiver, rotating
Quantity: _____

6 3S80=1 Sport (max 100kg)
Knee Joint
Quantity: _____

7 4R72*
Double Adapter

Size:
 32 45 60 75
Quantity: _____


10 4R218
Runner pyramid adapter,
rotatable
Quantity: _____

4R224
Runner junior pyramid
adapter, rotatable
Quantity: _____

4R216
Runner 4-hole adapter
(not pictured)
Quantity: _____

11 1E91 Runner
Quantity: _____

1E93 Runner Junior
Quantity: _____
Stiffness: _____



2 4R116
Lamination Anchor with
Pyramid Adapter, rotating
Quantity: _____

3 4R72*
Double Adapter
Size:
 32 45 60 75
Quantity: _____

4 2R58 (long)
Tube Adapter
Quantity: _____
or
2R57 (short)
Tube Adapter
Quantity: _____

5 4R82
Tube Clamp Adapter
Quantity: _____

8 2R58 (long)
Tube Adapter
Quantity: _____
or
2R57 (short)
Tube Adapter
Quantity: _____

9 4R82
Tube Clamp Adapter
Quantity: _____

12 2Z540 All-terrain Runner sole
Quantity: _____
2Z541 Spike Runner sole
Quantity: _____
2Z543 Runner Junior sole
Quantity: _____