8E600=*^{*}

DE Gebrauchsanweisung (Fachpersonal)	3
EN Instructions for use (qualified personnel)	19
FR Instructions d'utilisation (Personnel spécialisé)	36
IT Istruzioni per l'uso (Personale tecnico specializzato)	53
ES Instrucciones de uso (Personal técnico especializado)	69
PT Manual de utilização (Pessoal técnico)	85
NL Gebruiksaanwijzing (Vakmensen)	101
SV Bruksanvisning (Fackpersonal)	117
DA Brugsanvisning (Faguddannet personale)	132
NO Bruksanvisning (Fagpersonell)	148
FI Käyttöohje (Ammattihenkilöstö)	163
CS Návod k použití (Odborný personál)	179
RU Руководство по применению (Квалифицированный персонал)	195
JA 取扱説明書 (有資格担当者)	212

INFORMATION

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INFORMATION

Date of last update: 2021-04-15

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.


The 8E600=R/L AxonHook is referred to as the product/Axon-Bus terminal device below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

2 Product description

2.1 Design

The product consists of the following components:

1. Hook tip L/R
2. Cap nut cover
3. Gearbox with connected flexion wrist unit
4. Hook - wrist unit cuff

2.2 Function

The product is a myoelectrically controlled Axon-Bus gripping component in the modular Ottobock Axon-Bus prosthetic system.

The product consists of two hooks ("hook tips"). One fixed and one movable hook perform the gripping movements. Objects can be grasped by opening and closing the hook tips. A passive flexion wrist unit is permanently fixed to the product. This wrist unit offers a natural pattern of movement in both the "flexible mode" and the "rigid mode".

Power is supplied to the product by the battery which comes integrated in the socket.

Control settings for the product can only be carried out when all Axon-Bus components are connected. For further information, see the instructions for use provided with the Axon-Bus components.

Depending on the Axon-Bus prosthesis components currently used, specific standard programmes are available. With the 560X500 AxonSoft adjustment software, these standard programmes can be adapted/modified for the specific needs of each respective patient.

The product is equipped with the "neutral position" function. If there is no active muscle signal within a defined period of time, the product moves to the neutral position. The two hook tips are closed. The process is stopped if the hook tips encounter resistance or an obstacle. The period of time can be configured in the adjustment software.

2.3 Combination possibilities

This product is to be combined exclusively with components in the Axon-Bus prosthetic system. Components of the Ottobock MyoBock system or components from other manufacturers cannot be used with this product.

3 Intended use

3.1 Indications for use

The 8E600=* AxonHook is to be used **exclusively** for upper limb exoprosthetic fittings.

3.2 Conditions of use

The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports which involve excessive strain and/or shocks to the wrist joint (pushups, downhill racing, mountain biking, ...) or extreme sports (free climbing, paragliding, etc.). Furthermore, the product should not be used to operate motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment.

The product is intended **exclusively** for use on **one** patient. Use of the product by another person is not approved by the manufacturer.

Permissible ambient conditions are described in the technical data (see page 30).

3.3 Indications

- Amputation level transradial and transhumeral
- For unilateral or bilateral amputation
- Dismelia of the forearm or upper arm
- The user must be able to understand usage instructions and safety notices and put them into practice.
- The user must fulfil the physical and mental requirements for perceiving optical/acoustic signals and/or mechanical vibrations

3.4 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

3.5 Qualification

The fitting of a patient with the product may only be carried out by O&P professionals who have been authorised with the corresponding Ottobock training.

4 Safety

4.1 Explanation of warning symbols

 Warning regarding possible risks of accident or injury.

 Warning regarding possible technical damage.

4.2 Structure of the safety instructions

The heading describes the source and/or the type of hazard

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in case of failure to observe the hazard
- > E.g.: Consequence 2 in case of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

4.3 General safety instructions

Non-observance of safety instructions

Personal injury/damage to the product due to using the product in certain situations.

- ▶ Observe the safety instructions and the stated precautions in this accompanying document.

Risk of pinching between the hook tips

Injury due to pinching of body parts.

- ▶ Ensure that no body parts are between the hook tips when using or rotating the product.

⚠ CAUTION

Independent manipulation of the product

Injury due to product damage or malfunction.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

⚠ CAUTION

Signs of wear and tear on the product components

Injury due to faulty control or malfunction of the product

- ▶ In the interest of preventing patient injuries and to maintain product quality, we recommend performing regular service.
- ▶ Please see the service pass for more information about the service intervals. This is included in the purchased extended warranty packages.

⚠ CAUTION

Operator errors during the adjustment process with the adjustment software

Injury due to unexpected product behaviour.

- ▶ Participation in an Ottobock product training course is mandatory prior to using the product. During this product training course you will receive a password giving you access to the adjustment software. Additional product training courses may become necessary to qualify for software updates.
- ▶ Do not share your unlock PIN.
- ▶ Use the online help which is integrated into the software.

⚠ CAUTION

Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)

Falling due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Therefore, keeping a minimum distance of 30 cm to HF communication devices is recommended.

⚠ CAUTION

Proximity to sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)

Injury due to unexpected behaviour of the product caused by interference with internal data communication.

- ▶ Avoid remaining in the vicinity of visible or concealed theft prevention systems at the entrance/exit of stores, metal detectors/body scanners for persons (e.g. in airports) or other sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ When walking through theft prevention systems, body scanners or metal detectors, watch for unexpected behaviour of the product.

⚠ CAUTION

Remaining in areas outside the allowable temperature range

Injury due to faulty control or malfunction of the product.

- ▶ Avoid remaining in areas outside the allowable temperature range (see page 30).

⚠ CAUTION

Penetration of product with dirt and humidity

Injury due to unexpected product behaviour or malfunction.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.

⚠ CAUTION

Overloading due to unusual activities

Injury due to unexpected product behaviour as a result of a malfunction.

- ▶ The product was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist joint (pushups, downhill mountain biking, ...) or extreme sports (free climbing, paragliding, etc.).
- ▶ Careful handling of the product and its components not only increases their service life but, above all, ensures the patient's personal safety!
- ▶ If the product and its components have been subjected to extreme loads (e.g. due to a fall, etc.), then the product must be inspected for damage immediately. If necessary, forward the product to an authorised Ottobock Service Centre.

⚠ CAUTION

Storage of the product in the closed state

Injury due to faulty control or malfunction of the product as a result of damage to the sensors and mechanism.

- ▶ Store the product only in the neutral position or in the open state.

⚠ CAUTION

Changing Axon-Bus gripping components when switched on

Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- ▶ Power down the Axon-Bus prosthetic system by pressing the button in the charging receptacle before changing any Axon-Bus components (e.g. Axon-Bus gripping component).

⚠ CAUTION

Unintentional unlocking of the Axon-Bus gripping component

Injury due to releasing the Axon-Bus gripping component from the forearm (e.g. while carrying objects).

- ▶ The two release buttons should only be activated in order to change the Axon-Bus gripping component, deliberately and with due consideration of the respective situation.

⚠ CAUTION

Grasping objects with different gripping forces

Injury due to unexpected product behaviour.

- ▶ Note that the gripping force can vary depending on the characteristics (soft/hard) of the grasped object.

⚠ CAUTION

Incorrect electrode settings

Injury due to unexpected product behaviour.

- ▶ The electrodes are to be placed on intact skin only and with as much electrode-skin contact as possible. In the case of strong interference from electronic devices, the position of the electrodes should be checked and changed if necessary. If the interference cannot be eliminated or if you do not achieve the expected results by adjustment or selection of the appropriate control programme, please contact the Ottobock branch responsible for your country.
- ▶ Set the electrode gain as low as possible in order to reduce interference from powerful electromagnetic radiation (e.g. visible or concealed theft prevention systems at the entrance/exit of stores), metal detectors/body scanners for persons (e.g. in airports) or other sources of strong electromagnetic interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean).
- ▶ Use only the following products for cleaning/disinfecting the inner socket:
Cleaning: 453H10=1-N Ottobock DermaClean
Disinfection: Colourless, conventional medical disinfectants

INFORMATION

If you do not achieve the results expected by adjusting and selecting the appropriate control mode, please contact Ottobock Myo-Service.

5 Scope of Delivery and Accessories

5.1 Scope of delivery

- 1 pc. 8E600=* AxonHook
- 1 pc. Instructions for use (qualified personnel)
- 1 pc. Instructions for use (user)

5.2 Accessories

The following components are not included in the scope of delivery and must be ordered separately:

- 757B501 AxonEnergy Integral
- "560X500=* AxonSoft" adjustment software
- 757L500 AxonCharge Integral
- 13E500 AxonMaster (firmware 2.0.70 or higher)
- 9S503 AxonRotation (active rotation unit)

or

- 9S501 AxonRotation Adapter (passive rotation unit)

6 Charging the battery

The following points must be observed when charging the battery:

- The capacity of a fully charged battery is sufficient for one full day.
- We recommend charging the product once a day when used by the patient on a daily basis.
- The battery should be charged for at least 3 hours prior to initial use.

6.1 Connecting the power supply and battery charger

- 1) Slide the country-specific plug adapter onto the power supply until it locks into place.
- 2) Connect the round plug of the power supply to the receptacle on the battery charger so that the plug locks into place.

INFORMATION: Ensure correct polarity (guide lug). Do not use force when connecting the cable plug to the battery charger.

- 3) Plug the power supply unit into the wall socket.

6.2 Connect battery charger to the product

- 1) Connect the charging plug to the charging receptacle of the product.
→ The correct connection of the battery charger to the product is indicated by feedback (see page 33).
- 2) The charging process starts.
→ The Axon-Bus prosthetic system is switched off automatically.
- 3) Disconnect the product after the charging process is complete.

6.3 Display of the current charge level

Power supply and battery charger operating states

During the charging process the current battery charge level is indicated by the LEDs on the battery charger lighting up in sequence (Status signals).

Displaying the charge level on the product

The battery charge level can be queried at any time.

- 1) With the prosthesis switched on, press the charging receptacle button and hold for less than one second.
- 2) The LED display on the charging receptacle provides information on the current battery charge level (Status signals).

7 Preparing the product for use

7.1 Reducing the length of the wrist unit cuff

The length of the Hook's wrist unit cuff can be reduced as required for adaptation to the patient.

- ▶ Use scissors to cut along the marked grooves (see fig. 7).

7.2 Settings

Adjustments to the product can be made via Bluetooth data transmission using the 560X500=* AxonSoft software. For this purpose, a wireless connection must be established between the 13E500 AxonMaster and the PC with the help of the 60X5 BionicLink PC.

For further information, see the instructions for use provided with the 13E500 AxonMaster and the 560X500=* AxonSoft software.

INFORMATION

Cybersecurity

- ▶ Keep your operating system up to date and always install any available security updates.
- ▶ Protect your computer from unauthorised access (e.g., by using virus scans, password protection etc.).
- ▶ Do not use unsecured networks.
- ▶ Please contact the manufacturer if you suspect cybersecurity problems.

8 Use

8.1 Changing the Axon-Bus gripping component

Attaching the Axon-Bus gripping component to the socket

- 1) Slide the terminal device on to the socket (rotation anchored in the socket) until it locks in place audibly.
- 2) Verify proper attachment by pulling on the terminal device.

Separating the Axon-Bus gripping component from the socket

- 1) To switch off the Axon-Bus prosthetic system, press the charging receptacle button and hold for more than one second.
- 2) Press both unlock buttons on the flexion wrist unit.

INFORMATION: For safety reasons the Axon-Bus gripping component cannot be removed if only one unlock button is pressed.

- 3) Separate the Axon-Bus gripping component from the socket.

8.2 Switching the product on and off

CAUTION

Storage of the product in the closed state

Injury due to faulty control or malfunction of the product as a result of damage to the sensors and mechanism.

- ▶ Store the product only in the neutral position or in the open state.

- ▶ 1) Push and hold the charging receptacle button until a confirmation sound is heard (at least 1 second).
→ The prosthesis and prosthetic components are switched on.
Repeat this process to switch the prosthesis and prosthetic components off.

8.3 Opening the Axon-Bus gripping component in an emergency

This safety function allows the Axon-Bus gripping component to be opened regardless of the control signals present.

- 1) With the Axon-Bus prosthetic system switched on, press the charging receptacle button and hold for approximately three seconds until the Axon-Bus gripping component begins to open.
→ A pulsating beep sounds as it opens.
- 2) Releasing the button immediately stops the process of opening the Axon-Bus gripping component and turns the entire Axon-Bus prosthetic system off.

8.4 Proper handling of the Axon-Bus gripping component

Hand flexion, hand extension procedure

This function enables the flexion or extension of the Axon-Bus gripping component.

Flexible mode

- 1) Push the unlock lever of the Axon-Bus gripping component (to the stop) until it engages.
→ The wrist unit can be moved without engaging at the ratchet positions.
- 2) When the unlock lever is pushed again, the wrist unit engages at the next available position.

Rigid mode

- 1) Push the unlock lever of the Axon-Bus gripping component lightly (not to the stop) and move the wrist unit to the desired position.
- 2) When the unlock lever is released, the wrist unit engages at the next available position.

Pronation, supination procedure

⚠ CAUTION

Risk of pinching between the hook tips

Injury due to pinching of body parts.

- ▶ Ensure that no body parts are between the hook tips when using or rotating the product.

INFORMATION

The following steps apply only to the product "9S501 AxonRotation adapter" (mechanical rotation). When using the product "9S503 AxonRotation", pronation and supination are controlled by the associated muscle signals.

The wrist unit supports pronation and supination.

- ▶ Grasp the gearbox and turn the wrist unit to the desired position.

9 Cleaning

NOTICE

Improper product care

Damage to the product due to the use of incorrect cleaning agents.

- ▶ Only clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean).

- 1) Clean the product with a damp cloth and mild soap (e.g. 453H10=1-N Ottobock DermaClean) when needed.
Make sure that no liquids get into the product and product components.
- 2) Dry the product with a lint-free cloth and allow it to air dry fully.

10 Maintenance

Performing regular maintenance (service inspections) every 24 months is recommended to prevent injuries and maintain the quality of the product.

The grace period is no more than one month before or three months after maintenance is due.

In general, all products are subject to compliance with the maintenance intervals during the warranty period. This is the only way to maintain full warranty cover.

Additional services such as repairs may be provided in the course of maintenance. These additional services may be provided free of charge or can be billable according to an advance cost estimate, depending on the extent and validity of the warranty.

The following components must always be sent in for maintenance and repairs:

The product, battery charger and power supply. The shipping container for the loaner unit you receive must be reused for sending back the components requiring inspection.

10.1 Replacing hook tip

Required Tools

- Crowfoot spanner
- 9S305-1 Hook tip

Removing the hook tip

- 1) Slide the cap nut cover over the hook tip and remove it.
- 2) Loosen the cap nut with the spanner (see fig. 1).
- 3) Unscrew the cap nut (see fig. 2) and pull off the hook tip with the cap nut.

Mounting the hook tip

- 1) Place the hook tip on the connector of the gearbox (see fig. 4).

INFORMATION: Mount the hook tip so that the groove of the hook tip lies parallel over the pin of the connector (see fig. 3).

- 2) Screw on the hook tip.
- 3) Tighten the hook tip with the spanner (see fig. 5).

INFORMATION: The hook tip must be tightened using a torque of 17 Nm.

- 4) Pull the cap nut cover over the hook tip and position it over the cap nut (see fig. 6).

10.2 Initialising the Axon-Bus gripping component

The Axon-Bus gripping component is delivered in an initialised state.

⚠ CAUTION

Risk of pinching between the hook tips

Injury due to pinching of body parts.

- ▶ Ensure that no body parts are between the hook tips when using or rotating the product.

- > If the two hook tips no longer close completely, the following steps must be carried out for initialisation:
 - 1) Remove the AxonHook from the socket while opening the hook tips (see page 26).
 - 2) Switch off the Axon-Bus prosthetic system (see page 26).
 - 3) Reattach the AxonHook to the socket.
 - 4) Switch the Axon-Bus prosthetic system back on.
 - 5) Open the hook tips completely.
 - The hook tip is opened at low speed.
 - 6) Close the hook tips completely.
 - The hook tip is closed at low speed.
- The AxonHook has been initialised.
The complete opening width and maximum speed is available once again.

11 Legal information

11.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for

damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

11.2 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

11.3 CE conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

This product meets the requirements of the 2014/53/EU directive.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

11.4 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and
- 2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation.

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party:

Otto Bock Health Care, LP
3820 West Great Lakes Drive
Salt Lake City, Utah 84120-7205 USA

Phone + 1-801-956-2400

Fax + 1-801-956-2401

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website

<http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.

5470 Harvester Road

L7L 5N5 Burlington, Ontario

Canada

Phone + 1-800-665-3327

Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.

12 Technical data

Ambient conditions	
Storage (with and without packaging)	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing
Transport (with and without packaging)	-20 °C/-4 °F to +60 °C/+140 °F Max. 90% relative humidity, non-condensing
Operation	-5 °C/+23 °F to +45 °C/+113 °F Max. 95% relative humidity, non-condensing
Charging the battery	+5 °C/+41 °F to +40 °C/+104 °F Max. 85% relative humidity, non-condensing

General information	
Reference number	8E600=*
Min. opening width	130 mm/5.12 inch
Min. closing time	0.75 sec
Max. closing time	10 sec
Weight	440 g/15.52 oz
Expected lifetime given compliance with the recommended maintenance intervals	500,000 gripping cycles/5 years

The following gripping forces and load limits apply only when the battery of the Axon-Bus prosthetic system is fully charged and at room temperature.

Maximum gripping forces	
Gripping force at the hook tips	approx. 110 N
General wrist unit	
Flexion	75°/4 ratchet positions
Extension	45°/3 ratchet positions
Prosthesis battery	
Battery type	Li-Ion
Battery service life	2 years
Charging time until battery is fully charged	4 hours
Behaviour of the product during the charging process	The product is non-functional
Operating time of the product with fully charged battery	1 day with average use
Load limits	
Maximum lateral load on the hooks with locked wrist unit (e.g. holding an object/horizontal forearm)	16 kg/35.2 lb
Maximum load on the hooks in the closed hook position (e.g. carrying bags/vertical forearm)	50 kg/110.2 lb
Maximum vertical weight on the hooks (e.g. bracing)	60 kg/132.3 lb

13 Appendices

13.1 Symbols Used

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

Serial number (YYYY WW NNN)
 YYYY – year of manufacture
 WW – week of manufacture
 NNN – sequential number

Article number

Manufacturer

IP33

Protection against penetration of solid foreign objects with a diameter > 2.5 mm, protection against water spray

MD

Medical device

13.2 Operating states/error signals

The product indicates operating states audibly and visually.

13.2.1 Signals for operating states

Battery charger connected/disconnected


Beep signal	Event
1 x short	Battery charger connected
2 x short	Charging begins
1 x long	Charging ends

Switching the Axon-Bus prosthetic system on/off


Charging receptacle	Beeper	Event
	2 x short	Switch on the Axon-Bus prosthetic system
	1 x long	Switch off the Axon-Bus prosthetic system

13.2.2 Warnings/error signals

Error during use





Charging receptacle	Beeper	Explanation/correction
	1x long	Error in the Axon-Bus prosthetic system <ul style="list-style-type: none">Switch the Axon-Bus prosthetic system off/onOpen/close the Axon-Bus prosthetic systemConsult prosthetist if error persists
	-	Axon-Bus prosthetic system outside operating temperature range <ul style="list-style-type: none">Allow the Axon-Bus prosthetic system to warm up or cool down (Technical data)
-	3 x short	Axon-Bus prosthetic system switches itself off <ul style="list-style-type: none">Charge battery
-	Pulsating	Opening the Axon-Bus gripping component in an emergency <ul style="list-style-type: none">System is switched off, system can be used without restriction after being switched on again

Error while charging the product

Battery charger	Explanation/correction
	Axon-Bus prosthetic system is not charging <ul style="list-style-type: none"> • Disconnect and then reconnect the charging plug. • Clean charger contacts.
	Charger is overheated or overcooled <ul style="list-style-type: none"> • Disconnect charging plug and allow charger to cool down/warm up.
	Battery charger or power supply unit is defective <ul style="list-style-type: none"> • Consult prosthetist

13.2.3 Status signals

Battery charger is connected

Battery charger	Explanation
	Power supply and battery charger operational
	Battery is charging
	Battery is charging, battery 25% charged
	Battery is charging, battery 50% charged
	Battery is charging, battery 75% charged
	Battery is fully charged, battery 100% charged

Battery charge level

Charging receptacle	Event
	Battery is fully charged Battery charge level is 100%
	Battery charge level is over 50%
	Battery charge level is less than 50%

13.3 Glossary

The term "Axon" stands for **Adaptive exchange of neuroplacement data**. The Axon-Bus is an Ottobock innovation for the field of exoprosthetics: a data transmission system, derived from safety-related bus systems in the aviation and automobile industries. For the user this means enhanced safety and reliability because of a considerably reduced sensitivity to electromagnetic interference in comparison with conventional systems.

13.4 Directives and manufacturer's declaration

13.4.1 Electromagnetic environment

This product is designed for operation in the following electromagnetic environments:

- Operation in a professional healthcare facility (e.g. hospital, etc.)
- Operation in areas of home healthcare (e.g. use at home, use outdoors)

Observe the safety notices in the section "Information on proximity to certain areas" (Information on Proximity to Certain Areas).

Electromagnetic emissions

Interference measurements	Compliance	Electromagnetic environment directive
HF emissions according to CISPR 11	Group 1/class B	The product uses HF energy exclusively for its internal functioning. Its HF emissions are therefore very low, and interference with neighbouring electronic devices is unlikely.
Harmonics according to IEC 61000-3-2	Not applicable – power below 75 W	–
Voltage fluctuations/flicker according to IEC 61000-3-3	Product meets the requirements of the standard.	–

Electromagnetic interference immunity

Phenomenon	EMC basic standard or test procedure	Interference immunity test level
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air,
High-frequency electromagnetic fields	IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1 kHz
Magnetic fields with rated power frequencies	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz
Electrical fast transients/bursts	IEC 61000-4-4	± 2 kV 100 kHz repetition rate
Surges Line against line	IEC 61000-4-5	± 0.5 kV, ± 1 kV
Conducted interference induced by high-frequency fields	IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM and amateur frequency bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz

Phenomenon	EMC basic standard or test procedure	Interference immunity test level
Voltage drops	IEC 61000-4-11	0% U _T ; 1/2 period At 0, 45, 90, 135, 180, 225, 270 and 315 degrees
		0% U _T ; 1 period and 70% U _T ; 25/30 periods Single phase: at 0 degrees
Voltage interruptions	IEC 61000-4-11	0% U _T ; 250/300 periods

Interference resistance against wireless communication devices

Test frequency [MHz]	Frequency band [MHz]	Radio service	Modulation	Maximum power [W]	Distance [m]	Interference immunity test level [V/m]
385	380 to 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 to 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	1.8	0.3	28
710	704 to 787	LTE band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 to 960	GSM 800/90-0, TETRA 800, iDEN 820, CDMA 850, GSM 800/90-0, LTE band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1,720	1,700 to 1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1,845						
1,970						
2,450	2,400 to 2,570	Bluetooth WLAN 802.1-1 b/g/n, RFID 2450 LTE band 7	Pulse modulation 217 Hz	2	0.3	28
5,240	5,100 to 5,800	WLAN 802.1-1 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5,500						
5,785						

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- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

— Consult the dealer or an experienced radio/ TV technician for help.

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<http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.
5470 Harvester Road
L7L 5N5 Burlington, Ontario
Canada
Phone + 1-800-665-3327

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