4R220=1 Dynamic Vacuum System



DE	Gebrauchsanweisung	6
EN	Instructions for use	15
FR	Instructions d'utilisation	24
IT	Istruzioni per l'uso	34
ES	Instrucciones de uso	43
PT	Manual de utilização	53
NL	Gebruiksaanwijzing	63
SV	Bruksanvisning	72
DA	Brugsanvisning	81
NO	Bruksanvisning	90
FI	Käyttöohje	98
PL	Instrukcja użytkowania	107
HU	Használati utasítás	117
CS	Návod k použití	126
RO	Instrucțiuni de utilizare	135
HR	Upute za uporabu	144
SL	Navodila za uporabo	153
SK	Návod na používanie	162
BG	Инструкция за употреба	171
TR	Kullanma talimatı	180
EL	Οδηγίες χρήσης	189
RU	Руководство по применению	199
JA	取扱説明書	210
ZH	使用说明书	219
KO	사용 설명서	226

INFORMATION

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INFORMATION

Date of last update: 2020-06-04

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The 4R220* Dynamic Vacuum System (DVS) enables an active vacuum prosthesis. It consists of a component with an integrated piston pump, a liner and a sealing sleeve.

The Dynamic Vacuum System uses the pistoning movements between the residual limb and prosthetic socket to generate the vacuum. The piston of the pump is equipped with magnets. A metallic counter-piece is installed on the liner. Air is suctioned out of the space between the liner and prosthetic socket into the cylinder during the swing phase, and this air is pressed out during the stance phase. The result is a permanent vacuum, with a level that adapts to the user's activity level.

1.2 Combination possibilities

INFORMATION

To ensure proper product function, only use the product with the components specified in this section.

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

Limited combination options for Ottobock components

Allowable liner: 6Y94

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

The product is embedded in the prosthetic socket and is no structural component. The maximum body weight is determined by the prosthetic socket's laminate and the prosthetic components used.

2.3 Environmental conditions

Allowable environmental conditions

Temperature range for use: -10 °C to +60°C

Allowable relative humidity 0 % to 90 %, non-condensing

Unallowable environmental conditions

Mechanical vibrations or impacts

Perspiration, urine, fresh water, salt water, acids

Dust, sand, highly hygroscopic particles (e.g. talcum)

2.4 Service life

Depending on the patient's level of activity, the service life of the product is 3 to 5 years.

3 Safety

3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

3.2 General safety instructions

 **CAUTION**

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

 **CAUTION**

Use under unallowable environmental conditions

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

 **CAUTION**

Exceeding the service life and reuse on another patient

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded.
- ▶ Only use the product on a single patient.

⚠ CAUTION**Mechanical damage to the product**

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

⚠ CAUTION**Product contains a magnet**

Influence on or damage to devices or items due to a strong magnetic field

- ▶ Keep the product away from devices and items that are sensitive to magnetic fields.

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of delivery

4R220=1 Dynamic Vacuum System		
Quantity	Designation	Reference number
1	Instructions for use	-
1	Cylinder body	4X324=4
1	Piston	-
1	Fixation ring	4X339
1	Plaster screw	503S3
1	Valve bushing	-
1	Duckbill valve	-
1	Piston dummy	-
1	Silicone dummy	-
2	Valve dummy (1 short, 1 long)	-
1	Mounting wrench	4X338
1	Polylub GLY 801 special grease	633F30=2

Spare parts/accessories (not included in the scope of delivery)	
Designation	Reference number
Single-component pack (1 piston, 1 fixation ring)	4X320
Single-component pack (1 valve bushing, 2 duckbill valves)	4X322
Dummy set (1 piston dummy, 1 silicone dummy, 1 short valve dummy, 1 long valve dummy, 1 plaster screw)	4X326=1

5 Preparing the product for use

CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

INFORMATION

- ▶ Not all of the materials may be available in your country. In this case, please contact the local branch of the manufacturer to obtain information on alternative materials.

5.1 Preparing for socket fabrication

INFORMATION

The Dynamic Vacuum System can be used with a total surface weight-bearing socket or a specific weight-bearing socket. Observe the following instructions:

- ▶ Mould the residual limb length with liner unchanged while taking the plaster cast.
- ▶ Do not shorten the plaster model during modelling.

- > **Required materials:** piston dummy, silicone dummy, short or long valve dummy, 99B25 nylon stockinette, 99B81=70X19X5 PVA bag, string, 633W8 wax, 636K8* Plastaband
- 1) Distally flatten the plaster model so that the piston dummy can be positioned in the longitudinal residual limb axis (see fig. 2). Position the piston dummy on the plaster model and loosely secure it with a plaster screw.
 - 2) Pull a nylon stockinette over the plaster model.
 - 3) Tie off the nylon stockinette between the piston dummy and plaster model and trim off the excess nylon stockinette (see fig. 3).
 - 4) Turn the plaster screw all the way in (see fig. 4).

5) **INFORMATION: No PVA bag is required for fabricating a check socket.**

Soak the shorter PVA bag and pull it over the plaster model.

- 6) Use string to tie off the PVA bag in the undercut of the piston dummy (see fig. 5).
- 7) Cut off the excess PVA bag with a scalpel.
- 8) Position the silicone dummy on the piston dummy. Check that no excess PVA bag is protruding during this process (see fig. 6).
- 9) **INFORMATION: Do not close the opening in the cylinder body from the inside. Screw in the valve dummy after positioning the cylinder body.**
Position the cylinder body on the piston dummy. Perform the mediolateral alignment of the exhaust opening taking into account the prosthetic alignment.
- 10) Screw the valve dummy into the exhaust opening until the O-ring closes the opening (see fig. 7).
- 11) Close the slot in the valve dummy with Plastaband.

5.2 Fabricating the prosthetic socket

5.2.1 Optional: fabricating the check socket

- > **Recommended materials:** ThermoLyn (e.g. 616T52* ThermoLyn rigid, 616T83* ThermoLyn clear), 616S5* BetaSil
- 1) Prepare for fabrication of the prosthetic socket (see Page 19).
 - 2) Attach two strips of Plastaband lengthwise over the cylinder on both the anterior and the posterior side. This makes it easier to remove the cylinder after the trial fitting (see fig. 8).
 - 3) Attach a strip of Plastaband circularly over the proximal edge of the cylinder to seal the prosthetic socket (see fig. 8).
 - 4) Complete the vacuum forming process.
 - 5) Finish the prosthetic socket (see Page 21).
 - 6) **If the prosthetic socket leaks:** seal the transition from the cylinder body to the inside of the prosthetic socket with BetaSil.

5.2.2 Laminating the prosthetic socket

- > **Required materials:** 623T3=8 or 623T3=10 Perlon stockinette, 81A1=8 or 81A1=10 nylon stockinette, 616B1=25x* carbon fibre cloth strap, 99B81=100X19X5 PVA bag, 617H119 Orthocryl lamination resin 80:20 PRO, string, 627B40 polyethylene adhesive tape
- 1) Prepare for fabrication of the prosthetic socket (see Page 19).
 - 2) Cut a piece of Perlon stockinette (twice the length of the plaster model).

- 3) Pull a Perlon stockinette over the plaster model to the edge. Tie off the second half of the Perlon stockinette and fold it over the plaster model.
- 4) **INFORMATION: Tie this off tightly so the prosthetic socket is airtight after the lamination process.**
Wrap string around the valve dummy and tie the Perlon stockinette off tightly in the undercut of the cylinder body (see fig. 9).
- 5) Expose the valve dummy using a scalpel (see fig. 10).
- 6) Position a layer of carbon fibre cloth strap from the medial condyle over the cylinder body to the lateral condyle (see fig. 11).
- 7) Position a layer of carbon fibre cloth strap from the MPT (mid-patella-tendon) point over the cylinder body to the opening for the back of the knee (see fig. 12).
- 8) Tie off the carbon fibre cloth strap tightly in the undercut of the cylinder body (see fig. 13).
- 9) Position a layer of carbon fibre cloth strap circularly at the height of the MPT point (see fig. 13).
- 10) Expose the valve dummy.
- 11) Cut a piece of Perlon stockinette (twice the length of the plaster model).
- 12) Pull a Perlon stockinette over the plaster model to the edge. Tie off the second half of the Perlon stockinette and fold it over the plaster model (see fig. 14).
- 13) Soak a PVA bag and pull it over the plaster model.
- 14) Carry out the lamination process with Orthocryl.
- 15) **Once the lamination resin is sufficiently distributed:** shift excess lamination resin in the area of the cylinder body to the distal end by wrapping with polyethylene adhesive tape. Wrap tightly so the prosthetic socket becomes airtight.
- 16) Let the lamination resin set and roughen the laminate with wet sanding paper.
- 17) Select a socket attachment block/lamination anchor. If the prosthetic socket will be foamed with rigid foam, remove the short valve dummy and exchange it for the long valve dummy (see fig. 15, see fig. 16).
- 18) Carry out the second lamination process as described in the instructions for use of the selected socket attachment block/lamination anchor.
- 19) Finish the prosthetic socket (see Page 21).

5.3 Finishing the prosthetic socket

- > **Required materials:** mounting wrench, 633F30=2 Polyub GLY 801
- 1) Mark and trim the contour of the prosthetic socket.
 - 2) Sand free the valve dummy and remove the Plastaband.

- 3) Unscrew the valve dummy and retain it. After fabricating the check socket, the valve dummy is still required for the fabrication of the final prosthetic socket.
- 4) Remove the prosthetic socket from the plaster model.
- 5) Remove the piston dummy and silicone dummy from the plaster model and insert into the cylinder so it is closed (see fig. 17).
- 6) Close the opening in the piston dummy with Plastaband.
- 7) Screw in the valve dummy (see fig. 18).
- 8) Sand the contour of the prosthetic socket.
- 9) Unscrew the valve dummy.
- 10) Remove the piston dummy and silicone dummy from the prosthetic socket.
- 11) Wipe out the cylinder space with a lint-free cloth.
- 12) Grease the cylinder space and O-rings of the valve bushing with the supplied Polylub GLY 801.
- 13) **NOTICE! Do not allow the duckbill valve to come into contact with grease.**
Insert the duckbill valve into the valve bushing with the pointed end first (see fig. 19).
- 14) Screw the valve bushing into the prosthetic socket (see fig. 20, torque value: **3 Nm**).
- 15) Insert the piston into the cylinder to the stop (see fig. 21).
- 16) Insert the fixation ring into the thread of the cylinder and tighten it with the mounting wrench (see fig. 22).

6 Use

CAUTION

Using the product without a sealing sleeve

Risk of injury due to unexpected loosening of the connection

- ▶ Always use the product with a sealing sleeve so that a vacuum can be generated.

6.1 Donning the prosthesis

- 1) Apply the liner to the residual limb so that the distal end is aligned with the longitudinal axis of the residual limb. Unroll the liner onto the residual limb so that it has no wrinkles, without shifting the soft tissue.
- 2) Slide into the prosthetic socket with the liner. The magnets in the piston of the Dynamic Vacuum System connect to the metallic counter-piece in the liner.

- 3) Unroll the sealing sleeve over the prosthetic socket onto the thigh.

6.2 Doffing the prosthesis

- 1) Unroll the sealing sleeve from the thigh onto the prosthetic socket.
- 2) Doff the residual limb and liner out of the prosthetic socket.

7 Cleaning

CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

Dirt in the prosthetic socket or on the liner can impair the function of the Dynamic Vacuum System.

- 1) Rinse the product with clear fresh water.
- 2) Dry the product with a soft cloth.
- 3) Allow to air dry in order to remove residual moisture.

8 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.
- ▶ To lubricate the cylinder or piston seal, use **Polylub GLY 801**.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

9.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

10 Technical data

Reference number	4R220=1
Weight [g]	110
System height [mm]	27
Material	Aluminium

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