

17A0100

DE	Gebrauchsanweisung	4
EN	Instructions for use	12
FR	Instructions d'utilisation	20
IT	Istruzioni per l'uso	29
ES	Instrucciones de uso	37
PT	Manual de utilização	46
NL	Gebruiksaanwijzing	55
SV	Bruksanvisning	63
DA	Brugsanvisning	71
NO	Bruksanvisning	79
FI	Käyttöohje	87
IS	Notkunarleiðbeiningar	95
TR	Kullanma talimatı	103
JA	取扱説明書	111
ZH	使用说明书	118

2

INFORMATION

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INFORMATION

Last update: 2018-06-20

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.
- ▶ Please keep this document in a safe place.

These instructions for use provide important information for the processing of the 17AO100 ankle joint.

2 Product description

The 17AO100 ankle joint is a joint with a foot-lifting effect. It can be adapted to the individual requirements of the patient by combining the "adjustable stop" and "compression spring stop" modules.

2.1 Available sizes

The joint is available in one size 17AO100=22-T.

2.2 Components/design

Scope of delivery (see fig. 1)			
Item	Quantity [Piece (s)]	Article designation (wrench size)	Reference number
		Adjustable stop	
1	1	Set screw (T30)	
2	1	Stop bumper	
3	1	Stop sleeve with stop pin	
		Ankle joint	
4	2	Flat head screws (T30)	
5	1	Base body	
6	1	Joint screw (T30)	
7	2	Axial washers	
8	1	Foot stirrup	
9	1	Lamination dummy cover	
10	1	Foot stirrup bearing bushing	
11	1	Joint nut	
16	1	Lamination bar	
		Compression spring stop	
12	1	Spring guide with ball	
13	1	Compression spring	
14	1	Spring buffer	
15	1	Set screw (T30)	

Not included in the scope of delivery (see fig. 2)			
Item	Quantity [Piece (s)]	Article designation (wrench size)	Reference number
		Lamination dummy set	17LD1N=22
1	1	Shoulder screw (T30)	
2	1	Lamination dummy – lower part	
The cover for the lamination dummy (see fig. 1, item 9) is not included in the lamination dummy set.			

Service sets and spare parts (see fig. 1, see fig. 2)			
Item	Quantity [Piece (s)]	Article designation	Reference number
6	1	Joint screw	17AO100S=22-1
11	1	Joint nut	
12	1	Spring guide with ball	17AO100S=22-2
13	1	Compression spring	
14	1	Spring buffer	
2	1	Stop bumper	17AO100S=22-3
3	1	Stop sleeve with stop pin	
7	1	Axial disc set Sizes: 0.85, 0.90, 0.95, 1.0, 1.05, 1.10 mm	17LA33N=20
7	1	Set screw	30Y450=22
8	1	Foot stirrup	17SF100=OS-22
10	1	Bearing bushing	
16	1	Lamination bar	17S100=LB-20

3 Intended use

3.1 Indications for use

The ankle joint is intended **exclusively** for orthotic fittings of the lower limbs on one patient and for use only with the 17KO1=* C-Brace knee joint.

3.2 Area of application

Approved for a body weight of **max. 110 kg** with unilateral use.

3.3 Indications

Partial or total paralysis of the leg muscles.
Indications must be determined by the physician.

3.4 Qualification of the O&P professional

Adapting the product to the patient and configuring settings may only be performed by qualified personnel authorised after completing the corresponding training.

4 Safety

4.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

4.2 Safety instructions

 **CAUTION**

Excessive strain on load-bearing components

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

 **CAUTION**

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

 **CAUTION**

Incorrect alignment, assembly or adjustment

Risk of injury due to change in or loss of functionality

- ▶ Assembly, adjustment and maintenance operations may only be completed by qualified personnel.
- ▶ Observe the alignment, assembly and adjustment instructions.

⚠ CAUTION

Excessive strain due to use on more than one patient

Risk of injury and loss of functionality as well as damage to the product

- ▶ Use the product on only one patient.
- ▶ Observe the maintenance recommendations.

NOTICE

Thermal overloading of the product

Damage due to improper thermal treatment

- ▶ Do not carry out any heat treatment at temperatures above **300 °C (570 °F)**.
- ▶ Prior to thermal treatment, remove all temperature-critical components (such as plastic parts).

NOTICE

Thermal overloading of product components

Loss of function due to improper thermal treatment

- ▶ Prior to the thermal treatment of the product, remove all temperature-critical components (e.g. plastic parts).

NOTICE

Exposure of the product to unsuitable environmental conditions

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C or above +40 °C (e.g. sauna, excessive sunlight, drying on a radiator).

5 Preparation for use

The components of the ankle joint are pre-assembled on delivery. The components must be disassembled to adapt the joint to the patient. The joint assembly is laid out for the right side. If necessary, convert the positions of the "adjustable stop" and "compression spring stop" modules as well as the foot stirrup (see fig. 1).

5.1 Processing

Prerequisite: The patient-specific plaster model has been fabricated.

Bending the lamination bar

To avoid a risk of breakage, do not bend or otherwise modify the insertion zone on the bar. The allowable bending zone of the bar begins at the slotted hole (see fig. 1, item 16, hatched area).

- 1) Clamp bending irons in a vice.
- 2) Bend the bar with bending irons, noting the allowable bending zone.

Bending the foot stirrup

To avoid a risk of breakage, do not bend or otherwise modify the zone up to the point where the material tapers. The allowable bending zone of the foot stirrup starts where the material tapers (see fig. 1, item 8, hatched area).

- 1) Clamp bending irons in a vice.
- 2) Bend the foot stirrup with bending irons starting where the material tapers.

Laminating the orthosis

See the technical information for the C-Brace® orthotronic mobility system for further processing details: 646T5325 (prepreg technique) and 646T5326 (lamination resin technique).

The blasted zone on the lamination bar and foot stirrup identifies the lamination area to obtain a good form fit when laminating the orthosis. **Laminate the orthosis up to the insertion zone of the lamination bar and foot stirrup (see fig. 1).**

Required materials: 17LD1N=22 lamination dummy set, lamination dummy cover, 17SF100=OS-22 foot stirrup, 17S100=LB-22-T lamination bar, 501S136=M6x1-2-1 flat head screws, 636K8 Plastaband, 743Y56=4 alignment adapter

- 1) Mount the lamination bar on the lower section of the lamination dummy with the two flat head screws using a **T 30** Torx wrench.
- 2) Slide the lamination dummy cover onto the foot stirrup in order to press the bearing bushing out of the foot stirrup.
- 3) Save the bearing bushing for later use.
- 4) Place the foot stirrup in the lamination dummy cover.
- 5) Set the cover onto the lower section of the lamination dummy and secure it with the shoulder screw. Due to the design of the cover, it is brought to the neutral position.
- 6) Attach the lamination dummy to the plaster model by inserting the shoulder screw through the cover, foot stirrup and lower section and screwing it to the alignment adapter with a **T 30** Torx wrench.
- 7) Seal the screws of the bar and dummy with Plastaband.
- 8) Laminate the orthosis to fit the patient.

- 9) After the laminate hardens, remove the laminated blank from the plaster model. To do so, expose the lamination dummy with a knife and remove it.
- 10) Remove the Plastaband from the screws.
- 11) Remove the two flat head screws of the lamination bar.
- 12) Loosen the lamination dummy screws and remove the cover.
- 13) Expose the removal edges with the oscillating saw.
- 14) Remove the components from the plaster.

5.2 Assembly

Mounting the foot stirrup

- 1) Check the fit of the bearing bushing in the foot stirrup. If it is missing, insert a bearing bushing into the foot stirrup.
- 2) Position two axial washers on the foot stirrup, using one axial washer per side.
- 3) Slide the foot stirrup together with the two axial washers into the base body of the joint. To do so, the stops of the axial washers must be positioned in the receivers on the underside of the joint.
- 4) Make sure that the axial washers prevent direct contact between the metal basic piece and the foot stirrup.
- 5) Slide the split-pin bolt through the base body, axial washers and foot stirrup from the rear.
- 6) Install the joint screw from the front with a **T 30** wrench.
- 7) Before handing over the finished orthosis to the patient, secure the screw connection with Loctite 241 and a tightening torque **8 Nm**.

Mounting the bar

- 1) Insert the lamination bar into the base body of the joint.
- 2) Secure it with the two flat head screws and a **T 30** wrench, securing the screw connections with Loctite 241 and a tightening torque of **8 Nm**.

5.3 Adaptation

The joint is adapted to the needs of the patient through the individual combination of the "adjustable stop" and "compression spring stop" modules.

Prerequisite: The foot stirrup is mounted in the joint.

Installing the adjustable stop

- 1) Insert the stop sleeve with the integrated stop pin into the base body. When installing the stop, ensure that the tapered side of the stop is resting on the foot stirrup.
- 2) Then insert the stop bumper into the base body.
- 3) Install the set screw in the base body with a **T 30** wrench.

Installing the compression spring stop

- 1) The ball must be positioned in the receiver of the spring guide housing when installing the spring that lifts the foot. First slide the ball together with the spring guide housing into the base body (see fig. 1, item 12).
- 2) Position the compression spring and spring pin in the base body.
- 3) Install the set screw in the base body with a **T 30** wrench.

Patient trial fitting

- 1) Perform a trial fitting with the patient.
- 2) After final adjustment of the joint, apply Loctite 241 to the threads of the set screws and install them into the base body with a **T 30** wrench.

6 Cleaning

Promptly clean the joint after contact with water containing salt, chlorine or soap, or if it gets dirty.

- 1) Clean the joint with a damp cloth when needed.
- 2) Dry it with a lint-free cloth and allow it to air dry fully. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

7 Maintenance

INFORMATION

The product may be exposed to increased loads by the patient.

- ▶ Shorten the maintenance intervals according to the expected loads.

The manufacturer requires at least a semi-annual inspection of the product to verify functionality and check for wear.

As part of maintenance, inspect the product for wear and proper functionality.

Spare parts are listed in the section "Components/design".

Corresponding service sets can be ordered for servicing the product.

Do not use lubricants.

Mounting the axial washers

If the preinstalled axial washers become worn, they have to be replaced with new washers of the same thickness.

- ▶ **Optional:** Optimise the movement of the joint by replacing the axial washers step-by-step: Use thicker washers if there is play in the medial/lateral direction, thinner washers if there is resistance in plantar flexion/dorsal extension.

8 Disposal

Dispose of the product in accordance with national regulations.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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