
50P22 Rhizo Arexa

DE Gebrauchsanweisung	5	FI Käyttöohje.....	33
EN Instructions for use	7	PL Instrukcja użytkowania.....	36
FR Instructions d'utilisation.....	10	HU Használati utasítás	39
IT Istruzioni per l'uso	13	CS Návod k použití.....	42
ES Instrucciones de uso	16	SK Návod na používanie	44
PT Manual de utilização.....	19	ET Kasutusjuhis	47
NL Gebruiksaanwijzing.....	22	RU Руководство по применению	50
SV Bruksanvisning	25	JA 取扱説明書.....	53
DA Brugsanvisning	28	ZH 使用说明书.....	56
NO Bruksanvisning.....	30		

1


50P22





2

Right hand

Left hand

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INFORMATION

Date of last update: 2022-08-22

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.

- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50P22 Rhizo Arexa thumb brace.

2 Intended use

2.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the upper limbs and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

2.2 Indications

- Tendinopathy
- Ligament injury, ligament strain in the thumb
- Ligament rupture (partial rupture)
- Sprains of the saddle or basal joint of the thumb
- Subluxation/dislocation
- Instability of the saddle or basal joint of the thumb
- Paralysis of the opponens pollicis

The indication must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

None known.

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases / injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the device will be applied; sensory and circulatory disorders in the body region being treated and in the area of the fingers.

2.4 Lifetime

The product is designed for a maximum lifetime of **1 year**.

2.5 Mechanism of Action

The brace stabilizes the carpometacarpal and metacarpophalangeal joints.

3 Safety

3.1 Explanation of warning symbols

 **CAUTION**

Warning regarding possible risks of accident or injury.

 **NOTICE**

Warning regarding possible technical damage.

3.2 General safety instructions

 **CAUTION**

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

 **CAUTION**

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage

- ▶ Keep the product away from open flames, embers and other sources of heat.

 **NOTICE**

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

NOTICE**Use of a worn or damaged product**

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

NOTICE**Improper molding or application**

Damage to the product due to overloading of the material and improper fit of the product due to breakage of load-bearing components

- ▶ The product may only be molded by qualified personnel.
- ▶ Do not make any improper changes to the product.
- ▶ Always apply the product according to the information in the instructions.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Selecting the Size

- 1) Measure the wrist circumference.

- 2) Determine the brace size (see size chart).

4.2 Fitting and Application

CAUTION**Incorrect or excessively tight application**

Skin irritation due to overheating, local pressure points due to tight application to the body

- ▶ Ensure that the product is applied properly and fits correctly.
- ▶ When in doubt, do not continue to use the product if skin irritation occurs.

Adapting the side

The brace is pre-configured for the right hand upon delivery. For use on the left hand, remove the wrist strap, rotate it by **180°** and fasten it to the brace again (see fig. 2).

Adapting the splints

Two metal splints are located in the outer pockets of the brace. The two metal splints are anatomically pre-shaped and can be customized to fit the user (see fig. 3).

- 1) Pull the splints out of the pockets on the outside of the brace (see fig. 4).
- 2) Shape the splints, paying attention to molding them around the thumb.
- 3) Place the shaped splints back in the outer pockets.

Application

- > The wrist strap is open.
 - > The thumb strap is loosened.
- 1) Pull the brace over the thumb and wrist.
 - 2) Guide the wrist strap around the wrist and fasten it with the hook-and-loop closure (see fig. 5).
 - 3) Use the two hook-and-loop closures to fasten the thumb strap (see fig. 6).

→ The brace fits snugly without pinching. (Use the hook-and-loop closures to optimize the fit of the brace if necessary.)

4.3 Removal

- 1) Open the hook-and-loop closure of the brace.
- 2) Take off the brace.

5 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

- 1) Remove all metal splints.
- 2) Fasten all hook-and-loop closures.
- 3) Recommendation: Use a laundry bag or net.
- 4) Wash the brace in warm water at **40 °C** (104 °F) with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 5) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).
- 6) Reinsert the metal splints.

6 Disposal

Dispose of the product in accordance with national regulations.

7 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

7.1 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

7.2 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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