







## 50R300N Dyneva

<b>DE</b> Gebrauchsanweisung .....	5
<b>EN</b> Instructions for use .....	11
<b>FR</b> Instructions d'utilisation .....	17
<b>IT</b> Istruzioni per l'uso .....	24
<b>ES</b> Instrucciones de uso .....	31
<b>PT</b> Manual de utilização .....	37
<b>NL</b> Gebruiksaanwijzing .....	44
<b>SV</b> Bruksanvisning .....	50
<b>DA</b> Brugsanvisning .....	56
<b>NO</b> Bruksanvisning .....	62
<b>FI</b> Käyttöohje .....	68
<b>PL</b> Instrukcja użytkowania .....	74
<b>HU</b> Használati utasítás .....	81
<b>CS</b> Návod k použití .....	87
<b>SK</b> Návod na používanie .....	93
<b>TR</b> Kullanma talimatı .....	99
<b>ZH</b> 使用说明书 .....	105
<b>KO</b> 사용 설명서 .....	111

	Größe/Size	Abstand/Distance	
		cm	inch
	XS	a ~ 27; b ~ 30	a ~ 10.6; b ~ 11.8
	S	a ~ 30; b ~ 33	a ~ 11.8; b ~ 13.0
	M	a ~ 33; b ~ 36	a ~ 13.0; b ~ 14.2
	L	a ~ 36; b ~ 39	a ~ 14.2; b ~ 15.4
	XL	a ~ 39; b ~ 42	a ~ 15.4; b ~ 16.5
	XXL	a ~ 42; b ~ 45	a ~ 16.5; b ~ 17.7

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
### INFORMATION

Date of last update: 2020-04-24

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use contain important information for fitting and applying the 50R300 Dyneva lumbar spine flexion brace.

## 2 Product description

Item	Component	Item	Component
1	Connecting rods	7	Sacral pad
2	Y-hook-and-loop closure tension belt	8	Spring adapter with abdominal strap element
3	Tension belt	9	Cover caps
4	Abdominal strap closure	10	Upper tube frame
5	Abdominal strap segments	11	Back pads
6	Lower tube frame		

## 3 Intended use

### 3.1 Indications for use

The brace is intended **exclusively** for orthotic fittings of the back and **exclusively** for contact with intact skin.

The brace must be used in accordance with the indications.

### 3.2 Indications

With pain syndromes of the lumbar spine caused by:

- Lumbar spinal canal stenosis
- Lumbar facet joint syndrome
- Lumbar radicular syndrome (anterior thigh pain, lumboschialgia)
- Lumbar instability (spondylolisthesis, spondylolysis)
- Lumbar spondylarthrosis
- Intervertebral disc prolapse and protrusion in the lumbar spine area
- General degenerative changes with reduced intervertebral disc height

The indication must be determined by the physician.

### 3.3 Contraindications

#### 3.3.1 Absolute Contraindications

None known.

### 3.3.2 Relative Contraindications

Report to and discuss with the physician any of the following: contractures of the spine, restricted chest breathing ability, skin diseases/injuries, inflammation, raised scars with swelling, reddening, and warming in the area where the brace is fitted, tumors, lymphatic congestion – including unclear swelling of soft tissue distal to the medical aid, and dysesthesia in the legs.

### 3.4 Mechanism of Action

The three-point principle and the spring mechanism result in an effective and statistically significant reduction of the muscle-force-induced compression of the movement segments and the small vertebral joints in the lumbar spine region. Spine deformation occurs in both the sagittal and the frontal planes. This applies for moderate as well as maximum loads. The brace effectively, significantly and systematically reduces the load on the lumbar spine.

## 4 Safety

### 4.1 Explanation of warning symbols

**CAUTION** Warning regarding possible risks of accident or injury.

**NOTICE** Warning regarding possible technical damage.

### 4.2 General safety instructions

#### **CAUTION**

##### **Excessive strain on load-bearing components**

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

#### **CAUTION**

##### **Reuse on other persons and improper cleaning**

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

#### **CAUTION**

##### **Contact with heat, embers or fire**

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

#### **NOTICE**

##### **Contact with oils, salves, lotions or other products that contain oils or acids**

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

#### **NOTICE**

##### **Exposure of the product to unsuitable environmental conditions**

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C (14 °F) or above +50 °C (122 °F) (e.g. sauna, excessive sunlight, drying on a radiator).

#### NOTICE

#### Contact with saltwater or water containing chlorine/soap, fresh water or abrasive substances (e.g. sand)

Damage and premature product wear


- ▶ Following any contact with the substances identified above, promptly clean the product in accordance with the section "Cleaning".

## 5 Handling

#### INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

### 5.1 Selecting the Size

- > Required tools: 743S1=60 diameter or 743E10 measuring tool
- 1) Measure the width of the torso **4 cm** below the shoulder blades (dimension A).
- 2) Measure the width of the torso at the level of the iliac crest (dimension B).
- 3) Determine the size of the tube frames using the measured dimensions and the size chart (see Page 3).
- 4) **INFORMATION: Determine the size of the brace (XS/S/M/L/XL/XXL) on the basis of the width of the torso below the shoulder blades. If necessary, order the lower tube frame one size smaller or bigger and then replace it (reference number 29R414=XS/S/M/L/XL/XXL).**

### 5.2 Adaptation

#### NOTICE

#### Improper changes

Damage to the product due to overloading of the material and improper fit of the product due to breakage of load-bearing components

- ▶ Do not make any improper changes to the product.

The Dyneva is assembled on delivery and must be adapted to the requirements of the patient.

- 1) Unfasten the hook-and-loop closure of the abdominal strap.
- 2) Place the brace frame on the patient's back in such a way that the sacral pad is positioned on the sacrum, and the back pads are positioned approximately **4 cm** below the shoulder blades (see fig. 1).
- 3) Optional: Adapt the fit by bending the tube frames (see fig. 2) or replacing the lower tube frame (see the option "Replacing the lower tube frame").
- 4) Wrap around the abdominal strap (see fig. 3) and determine the required number of abdominal strap segments, paying attention to the minimum and maximum overlapping of the strap ends (see fig. 4).
- 5) Remove unnecessary abdominal strap segments by pulling out the connecting rods (see fig. 4).
- 6) Reconnect the abdominal strap segments to be used with a connecting rod.
- 7) Apply the brace frame on the back.
- 8) Fasten the abdominal strap, ensuring that it is positioned evenly. The strap must not cut into the abdomen.
- 9) Attach the Y-hook-and-loop closures for the tension belts to the abdominal belt and determine the optimal length of the tension belts (see fig. 5).  
**NOTICE! In doing so, ensure that the angle of the springs does not exceed 20° when the tension belts are taut (see fig. 6).**
- 10) Shorten the length of the tension belts with scissors (see fig. 5).
- 11) Attach the Y-hook-and-loop closures to the strap ends and fasten them to the abdominal strap with hook-and-loop.

**Optional: Adjusting the height position of the back pads**

- 1) **INFORMATION: To adjust the height of the back pads, slide the upper tube frame inside the spring adapters. The maximum pull-out height is indicated on the tube frame (see fig. 7).**
- 2) Loosen the upper screw of the cover on the upper tube frame.
- 3) Slide the tube frame in the holder, noting the maximum pull-out height.
- 4) Tighten the screw.

**Optional: Horizontally adjusting the back pads**

- 1) **INFORMATION: The back pads can be inserted in one of the three holes for horizontal positioning. Install all cover caps after positioning (see fig. 7).**
- 2) Remove the back pad cover caps.
- 3) Loosen the back pads from the brace frame.
- 4) Insert the back pads into one of the three holes.
- 5) Reinstall the cover caps.

**Optional: Replacing the lower tube frame**

- 1) **INFORMATION: Optimize the fit of the lower tube frame by replacing the tube frame. The replacement size may only be one size smaller or bigger (reference number 29R414=XS/S/M/L/XL/XXL).**
- 2) Loosen the screws on the lower tube frame.
- 3) Replace the tube frame.
- 4) Tighten the screws.

**Optional: Adjusting the angle of the abdominal strap**

- 1) Optimize the fit of the abdominal strap by adapting the abdominal strap angle (see fig. 8).
- 2) Loosen the screws of the covers for the upper tube frame.
- 3) Pull off the spring adapter covers.
- 4) Adjust the hook-and-loop position of the abdominal strap on the tube frame.
- 5) Replace the covers.
- 6) Tighten all screws.

## Optional: Attaching lateral pads to the tube frame

- 1) To improve wearer comfort, attach and adapt the lateral pads, which can be ordered separately, to the ends of the tube frames. One set of lateral pads (29R390) is required for the upper and lower tube frames.
- 2) Loosen the screws of the tube frame.
- 3) Remove the tube frame.
- 4) Slide the lateral pads onto the tube frame.
- 5) Insert the tube frame into the spring adapters and tighten the screws.

### 5.3 Application

#### **⚠ CAUTION**

##### **Direct skin contact with the product**

Skin irritation due to chafing or sweat formation

- ▶ Do not wear this product directly on the skin.

#### **⚠ CAUTION**

##### **Incorrect or excessively tight application**

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

#### **NOTICE**

##### **Use of a worn or damaged product**

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.
- ▶ Do not continue using a product that is no longer functional, or that is worn or damaged.

- 1) Place the brace frame in such a way that the sacral pad is positioned on the sacrum, and the two back pads are positioned below the shoulder blades.
- 2) Fasten the abdominal strap.

- 3) Evenly tighten the tension belts and attach them to the positions identified during the fitting by means of hook-and-loop closures.

## 6 Cleaning

### NOTICE

#### Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the brace regularly:

- 1) Fasten all hook-and-loop closures.
- 2) Hand-wash the straps and pads in warm water at **30 °C (86 °F)** using ph-neutral detergent. Rinse thoroughly.
- 3) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove, or radiator).
- 4) Clean contamination off the brace frame with a moist cloth and allow to air dry.

## 7 Disposal

Dispose of the product in accordance with national regulations.

## 8 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

### 8.1 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use after this chapter.

### 8.2 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

### 8.3 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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