## 17H100

<b>DE</b> Gebrauchsanweisung .....	3
<b>EN</b> Instructions for use .....	21
<b>FR</b> Instructions d'utilisation.....	39
<b>IT</b> Istruzioni per l'uso .....	57
<b>ES</b> Instrucciones de uso.....	75
<b>NL</b> Gebruiksaanwijzing .....	93
<b>SV</b> Bruksanvisning.....	111
<b>NO</b> Bruksanvisning.....	129
<b>FI</b> Käyttöohje .....	147
<b>TR</b> Kullanma talimatı .....	165

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

Diagnose: \_\_\_\_\_

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## Table of Contents

<b>1 Description</b>	21	<b>4 Maintenance</b>	35
<b>2 Components and construction</b>	24	<b>5 Disposal</b>	35
<b>3 Mounting instructions and Trial Fitting</b>	30	<b>6 Legal information</b>	35

### INFORMATION

Last update: 2020-04-09

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document in a safe place.

### Explanation of symbols

**⚠ CAUTION** Warnings regarding possible risks of accident or injury.

**NOTICE** Warnings regarding possible technical damage.

**INFORMATION** Additional information on the fitting/use.

### INFORMATION

For a correct fitting, the following information must be observed strictly in order to ensure safety and function.

## 1 Description

### 1.1 Indications for use

The 17H100 RGO Hip Joint System is intended exclusively for the orthopaedic fitting of the lower extremities with reciprocating gait orthoses.

### 1.2 Indications and contraindications

The Reciprocating Hip Joint System (RGO) was specifically designed for patients with a body weight up to **65 kg**. It can be used with patients suffering from spina bifida or myopathies of other genesis (e.g. traumatic paraplegia) with a lesion level from Th5 to L2.

Restrictions exist for patients with ICP, motor perception disorders, deficits of movement of the upper extremity, insufficient muscle strength in the shoulder girdle, loss of trunk mobility in the frontal

or sagittal plane, severe deformities of the skeletal system (e.g. scoliosis, dislocations), hip flexion contracture (> 20°), knee joint flexion contracture (> 15°), and torsional deformities of the legs.

### 1.3 Safety instructions

#### ⚠ CAUTION

#### Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

**⚠ CAUTION**

**Risk of falls due to improper use.** Explain the function and effects of the RGO hip joint system to your patients and to those persons who must be familiarized with handling of the joint system.

**⚠ CAUTION**

**Risk of injury as a result of improper use.** No modifications may be made to the RGO hip joint system, other than the tasks described in these Instructions for Use.

Disassembly and repair of the sitting and walking joints may only be performed by Ottobock staff.

**⚠ CAUTION****Incorrect alignment or assembly**

Injuries due to changes in or loss of functionality

- ▶ The product may only be installed by trained, qualified personnel.
- ▶ Observe the alignment and assembly instructions.
- ▶ **Inform the patient.**

**⚠ CAUTION****Mechanical damage to the product**

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product.
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service etc.).
- ▶ **Inform the patient.**

**⚠ CAUTION****Excessive strain on load-bearing components**

Injuries due to changes in or loss of functionality

- ▶ Only use the product for the defined area of application.
- ▶ If the product has been exposed to extreme strain (e.g. due to falling), take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service etc.).
- ▶ **Inform the patient.**

**⚠ CAUTION****Improper installation**

Patient injuries due to falling, loss of functionality

- ▶ The product may only be installed by trained, qualified personnel.

**NOTICE****Exposure of the product to unsuitable environmental conditions**

Damage, brittleness or destruction due to improper handling

- ▶ Avoid storage in condensing ambient humidity.
- ▶ Avoid contact with abrasive substances (e.g. sand, dust).
- ▶ Do not expose the product to temperatures below -10 °C or above +60 °C (e.g. sauna, excessive sunlight, drying on a radiator).
- ▶ **Inform the patient.**

Fig. 1  
Initial position. The legs are positioned parallel to each other in the frontal plane.

Fig. 2  
At heel strike with rotated pelvis the legs remain oriented in walking direction.

## 1.4 Function

The connection between the leg bars and the pelvic module is especially important for fittings with reciprocating gait orthoses.

The RGO hip joint system was developed specifically for a largely physiological rotation of the pelvis when walking. A bi-axial axis design allows for a pelvic rotation of 15° while maintaining the foot position in the walking direction while walking. This leads to greater freedom of movement and improved mobility for the user (see Fig. 1/2).

The following special features distinguish the Otobock RGO Hip Joint System

- Easy to use modular system.
- Replacing individual components allows the orthosis to "grow" with the user.
- Lightweight alloy pelvic tube with high torsional strength allows for modular connection of joints.
- The bi-axial axis system allows for a pelvic rotation of 15° in the orthosis without changing the walking direction.
- Smooth and maintenance-free power transmission during walking through the push-pull cable system.
- Sitting joint lock release with pre-release of the locking mechanism.
- Safety button for re-locking in case of accidental pre-release of the sitting joint.
- Choice of thermoplastic or laminate materials for the integration of an individually fabricated pelvic module.
- Torso bar is easily disassembled.
- Torso bar hip flexion adjusts by up to 10°.
- Prefabricated lightweight alloy thigh bars with 0 mm, 5 mm and 10 mm offset available.
- The thigh bar is easy to connect to and disconnect from the hip joint.

Pelvic width, clear width in mm	Pelvic tube Article number	Joint size
200	17Y140=200	17H100=2
210	17Y140=210	17H100=2
220	17Y140=220	17H100=2
230	17Y140=230	17H100=2
240	17Y140=240	17H100=2
250	17Y140=250	17H100=2
260	17Y140=260	17H100=2
270	17Y130=270	17H100=1
280	17Y130=280	17H100=1
290	17Y130=290	17H100=1
300	17Y130=300	17H100=1
310	17Y130=310	17H100=1
320	17Y130=320	17H100=1
330	17Y130=330	17H100=1
340	17Y130=340	17H100=1
350	17Y130=350	17H100=0
360	17Y130=360	17H100=0
370	17Y130=370	17H100=0
380	17Y130=380	17H100=0

Table 1

### 1.5 Size selection

For children ages 4 – 8 with a body weight up to 34 kg: 17H100=2.

For patients 9 and older up to a body weight of 65 kg: 17H100=1 or 17H100=0.

Pelvic tubes (Fig. 3, item 2) from 200 mm to 380 mm in 1 cm increments (Table 1) are available to adjust the pelvic width.

When measuring, take the padding and wall thickness of the pelvic module being fabricated into account (see last page of the order form).

**INFORMATION**

We recommend using size classification according to the UL Joint System for fabrication of the leg orthosis components.

Refer to 646D197=D for information.

## 2 Components and construction


Fig. 3

### 2.1 Components (Fig. 3 and 4)

1. **17H100 Joint System** consisting of:

- 1.1. Hip joint (Lightweight alloy, Niro)
- 1.2. Push-pull cable
- 1.3. 17Y132/142 Torso bar (lightweight alloy)
- 1.4. 17Y131/141=10 Thigh bar (lightweight alloy)
- 1.5. 501S32=M4x10x10 Truss head screw (Niro)
- 1.6. 636W28 OttoBock Special Adhesive
- 1.7. 616Z9 Shrink Tubing

2. **17Y130/140 Pelvic Tube** (lightweight alloy)


Fig. 4


Fig. 5

### 2.1.1 Service parts for 17H100=1/0

17Y130=*	Pelvic tube
17Y132=L	Torso bar, left
17Y132=R	Torso bar, right
17Y131=0	Thigh bar 0 mm offset
17Y131=5	Thigh bar, 5 mm offset
17Y131=10	Thigh bar, 10 mm offset
501S32=M4x10x9.5	Truss head screw
501S32=M4x10x9.5	Truss head screw
636W28	Ottobock Special Adhesive
616Z9=25.4	Shrink Tubing

### 2.1.2 Service parts for 17H100=2

17Y140=*	Pelvic Tube
17Y142=L	Torso bar, left
17Y142=R	Torso bar, right
17Y143=0	Thigh bar 0 mm offset
17Y143=5	Thigh bar, 5 mm offset
17Y143=10	Thigh bar, 10 mm offset
501S32=M4x10x10	Truss head screw
636W28	Ottobock Special Adhesive
616Z9=19	Shrink Tubing

#### INFORMATION

The component kit contains 2 thigh bars (Fig. 4, item 4) with 10 mm offset (Fig. 5). In addition, thigh bars with and without a 5 mm offset are available.


Fig. 6


Fig. 7

## 2.2 Construction and effects

When walking with the orthosis, the movements are transferred reciprocally; for this purpose the two hip joints are functionally coupled to each other via a push-pull-cable (Fig. 3, item 1.2). The hip joints consisting of an upper joint section (a), middle joint section (b), and the lower joint section (c), which together function as a complex joint mechanism and are adjusted accordingly. The two hip joints are connected by a lightweight alloy pelvic tube (Fig. 3. item 2). The selected materials and construction details ensure strength and safety.

The RGO Hip Joint System is delivered mostly preassembled.

### INFORMATION

To ensure safe function of the joint mechanism, the hip joints must not be disassembled. The push and pull cable may not be adjusted or loosened under any circumstances.

The biaxial joint system consists of a sitting joint **S** and a walking joint **G**. The walking joint is locked when the patient is sitting. The walking joint is inclined by 35° in its axis, resulting in a 15° inward and 15° outward rotation. This results in 15° of pelvic rotation when the legs are positioned for walking straight ahead. The joint system can be adjusted to an actual pelvic width of approx 200 mm to approx. 380 mm by selecting corresponding pelvic tubes. The torso bars (Fig. 3, item 1.3) can be inclined by up to 10° ventrally to take an abducted hip position into account (see page 11, Fig. 18). The individually shaped pelvic module (fabricated of lamination resin or thermoplastic material, for example) is mounted between the torso bars. The components of the leg orthosis are attached rigidly to the Hip Joint System with the corresponding thigh bars.

### INFORMATION

The pelvic tube is important for stability and must be mounted for the following working steps (see chapter 2.3, Fig. 8). Position the sitting joint (S) at the anatomical hip rotation point (see chapter 2.4, page 10).

Fig. 8

### ***2.3 Connecting the joint components to the pelvic tube***

Insert the pelvic tube all the way into the tube adapter of the upper joint section (required length = 25 mm).

Fig. 9

The joint components must be positioned at 90° to the pelvic tube.

Fig. 10

Fix the pelvic tube in the tube adapter by tightening with a 5 mm Allen wrench.

Fig. 11

Fig. 12

Fig. 13

## ***2.4 Joint mechanism functionality***

When the release button is pressed, the sitting joint lock mechanism will be pre-released with an audible click.

If the mechanism has been pre-released by mistake, it can be re-locked by pressing the mechanism's safety button.

Complete release of the sitting joints is only possible in the neutral position (thigh bars are parallel). The lower joint sections are brought to the parallel position from a forward or backward stepping position, with the walking joint being blocked and the sitting joint being released in the process.


Fig. 14

Now the joint can be flexed up to the stop. During flexion, the safety bolt automatically returns to its initial position with an audible click, so that the sitting joint will be locked again during subsequent extension.

Fig. 15

To verify the correct width of the pelvic tube, hold the joint system over the fabricated pelvic module and check the clear width (allow room for padding, etc.).

Fig. 16

The pelvic tube may be shortened if necessary (possibly after trial fitting).

Simply cut the tube to the desired length.

Deburr the end of the tube.

#### **INFORMATION**

Important information for alignment of the orthosis: During alignment of the orthosis, the mechanical axis of the sitting joint (S) must be lined up with the anatomical hip rotation point.

This will ensure the most appropriate physiological movement due to the oblique axis of the walking joint (G).

### 3 Mounting instructions and Trial Fitting

Fig. 17

Fig. 18

#### 3.1 Adjusting and disassembling the Torso Bar

Upon delivery, the torso bars are assembled straight. A hip flexion up to 10° can be set. To do so, loosen the 4 mm Allen head screws.

We recommend using this adjustment only during trial fitting.

The torso bar can be removed from the upper joint section to make adjusting easier. Remove the 4 mm Allen head screws (Fig. 18) remove the bolt from the inner side of the joint using a pin wrench. The 711S5=7x9 bending irons are recommended for bending the thigh bars.

#### NOTICE

**Damage while working the surface of the lightweight alloy:** No heat treatment for forming. Do not bend the lightweight alloy bars using sharp-edged bending irons. Avoid a bending radius smaller than **10 mm**. Remove surface damage such as scoring or burrs which result from forming through grinding, brushing, or polishing (e.g. 749Z12=160x16 sanding brush), otherwise a fracture risk develops at this location through stress concentration.

Fig. 19

For assembly, re-insert the bolt into the torso bar and secure with 636K13 Loctite.

Fig. 20

To use the maximum flexion position, the finished torso bars must be parallel at the stop.

For hip flexion adjustment, see page 11, Fig. 17.

Fig. 21

### **3.2 Adjusting and mounting of the thigh bars**

The thigh bars are available pre-formed in order to obtain a high stability for the bar components. The offset at the bar head end is 10 mm (see Fig. 5). **Forming with bending irons is intended only for correct abduction and adduction positioning as well as for rotation of the leg orthosis components.**

The 711S5=7x9 bending irons are recommended for bending the thigh bars.


Fig. 22

**INFORMATION**

Please observe the instructions for working light-weight alloy (see bottom of page 11).

After bending and alignment of the orthosis, the thigh bars are fixed to the lower joint section tube adapter with the supplied screws (501S32=M4x10x10) for trial fitting.

Fig. 23

**3.3 Trial Fitting of the orthosis**

For the recommended use of the unilateral joint bar system, plug connectors for the thigh bars which make length adaptation much easier are available (Fig. 23).

Depending on the design of the orthosis the thigh bars can either be used without leg shells in the thigh area, or you may attach leg shells to the thigh bars, preferably in the area of the connectors. A fixing clamp (different versions would be possible: e.g. metal; carbon fibre) for additional reinforcement of the thigh shells is recommended to provide rotational stability.

**INFORMATION**

The walking orthosis must be aligned so that the two thigh bars are mounted parallel to one another on the leg orthosis components (especially in case of unequal flexion contractures in the hip or knee joints).

Otherwise the functionality of the release mechanism in the sitting joint is impaired, since the sitting joint can only be unlocked when both thigh bars are parallel to one another!

Fig. 24

With the orthosis applied, the sitting joint can be pre-released in all positions by pressing the release button (Fig. 11). Under this **load-bearing** condition, the sitting joint only releases when the legs are parallel and **the hip is also over-extended briefly**.

In case of incorrect alignment see chapter 3.3.

Fig. 25

Released sitting joint while seating the patient.

To finish the orthosis, the thigh bars must be glued in place with the included 636W28 Special Adhesive (**observe processing instructions**). If the bar must be replaced, heat the adhesive with a hot air gun so that the bar can be pulled from the lower joint section.

#### INFORMATION

Processing instructions for 636W28 Ottobock Special Adhesive:

Degrease the surfaces to be glued with 634A1 Thinner or Solvent, thoroughly mix the contents of the two tubes of special adhesive and hardener, and apply the mixture to the surfaces to be glued. The adhesive takes 3 to 4 hours to set. Final bonding strength will be reached after 16 hours.

Fig. 26


Fig. 27

### 3.4 Gluing the Pelvic Tube

#### ⚠ CAUTION

**Risk of falling.** For safety reasons, the pelvic tube must be secured in the tube adapter of the joint upper section with 636W28 Ottobock Special Adhesive before final delivery of the orthosis to the patient. Observe processing instructions for the Special Adhesive (see page 10, beside Fig. 25).

Cut the pelvic tube if required (Fig. 16).

Degrease the surfaces to be glued (pelvic tube and tube adapter) and roughen with fine sandpaper. Next, apply Special Adhesive to both surfaces and insert the pelvic tube into the tube adapter up to the stop (25 mm length). Tighten the 5 mm Allen head screw so that the position of the pelvic tube **can still be changed** easily. The resulting thin layer of adhesive between the pelvic tube and the tube adapter is necessary for a durable bond.

Verify the orthosis' alignment in the assembled state and let the adhesive set.

## 4 Maintenance

In the interest of the patient, the joint system must be inspected at least once every 6 months or sooner in specific cases.

### **⚠ CAUTION**

#### **Locking mechanism wear**

Injuries caused by loss of functionality

- ▶ Check the locking mechanism for wear.
- ▶ In case of wear, replace the corresponding components.

## 5 Disposal

Dispose of the product in accordance with national regulations.

## 6 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### **6.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### **6.2 CE conformity**

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

## Order Form

**The Reciprocating Hip Joint System (RGO) was developed specifically for patients using the hip joint system 17H100**

*Please fax this order form to your Ottobock service facility.*

Company: \_\_\_\_\_ Technician: \_\_\_\_\_

Tel: \_\_\_\_\_ Fax: \_\_\_\_\_ Date: \_\_\_\_\_

### General instructions:



The RGO Hip Joint System is suitable for patients up to 65 kg. The system is available in three different versions. 17H100=2 for ages 4-8, up to a body weight of 34 kg. 17H100=1/ 0 for ages 9 and up, up to a body weight of 65 kg. Pelvic tubes are available in 1 cm increments for clear pelvic widths of 200 mm to 380 mm.

The pelvic width **must be measured** between the two trochanters in both the sitting and prone position.

When selecting the pelvic tube, also take into account the wall thickness of the pelvic module to be fabricated including the padding!

**Please checkmark the desired combination!**

Combination	17H100=2 up to 34 kg	17H100=1 up to 65 kg	17H100=0 up to 65 kg
17Y140=210	<input type="radio"/>		
17Y140=220	<input type="radio"/>		
17Y140=230	<input type="radio"/>		
17Y140=240	<input type="radio"/>		
17Y140=250	<input type="radio"/>		
17Y140=260	<input type="radio"/>		
17Y130=270		<input type="radio"/>	
17Y130=280		<input type="radio"/>	
17Y130=290		<input type="radio"/>	
17Y130=300		<input type="radio"/>	
17Y130=310		<input type="radio"/>	
17Y130=320		<input type="radio"/>	
17Y130=330		<input type="radio"/>	
17Y130=340		<input type="radio"/>	
17Y130=350			<input type="radio"/>
17Y130=360			<input type="radio"/>
17Y130=370			<input type="radio"/>
17Y130=380			<input type="radio"/>

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