

4R110

4R110




DE	Gebrauchsanweisung	4
EN	Instructions for use	9
FR	Instructions d'utilisation	14
IT	Istruzioni per l'uso	19
ES	Instrucciones de uso	24
PT	Manual de utilização.....	29
NL	Gebruiksaanwijzing.....	34
DA	Brugsanvisning.....	39
PL	Instrukcja użytkowania	44
HU	Használati utasítás	49
CS	Návod k použití.....	54
EL	Οδηγίες χρήσης	59

4R110

	■	●	▲
DE	Einzelteile	Einzelteile-Pack	Mindestmenge
EN	Single components	Single-Component Pack	Minimum order quantity
FR	Pièces à l'unité	Kit de pièces de rechange	Quantité minimum
IT	Singoli componenti	Imballo dei singoli componenti	Quantità minima
ES	Componentes	Kit componentes	Cantidad mínima
PT	Componentes	Conjunto de peças de reposto	Quantidade mínima
NL	Onderdelen	Onderdelenpakket	Minimikvantitet
SV	Separata delar	Service-Set	Minimum oantal
DA	Komponenter	Komponentpakke	Mindste mængde
NO	Enkle komponenter	Enkeltkomponentpakke	Minste ordremengde
PL	Podzespoły	Zestaw naprawczy	Ilość minimalna
HU	Alkatrészek	Alkatrész-csomag	Minimum mennyiség
CS	Jednotlivé díly	Sada jednotlivých dílů	Minimální množství pro objednání
RO	Componente individuale	Set componente individuale	Cantitate minimă
HR	Pojedinačni dijelovi	Pakiranje pojedinačnih dijelova	Minimalna količina
TR	Münferit parçalar	Münferit parça ambalajı	Minimum miktar
EL	Μεμονωμένα εξαρτήματα	Σετ μεμονωμένων εξαρτημάτων	Ελάχιστη ποσότητα
RU	Отдельные детали	Комплект деталей	Минимальное количество
JA	単一成分パック	単一成分パック	最低発注量
ZH	单个部件	维修组件	最低起订量

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Date of last update: 2020-11-30

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

1 Scope of Delivery

4R110 Lamination Anchor with Pyramid Receiver for Children
647H252 Instructions for Use

2 Single Components (Fig. 1)

2.1 Single Component ■

–

2.2 Single-Component Pack ●

–

2.3 Minimum Order Quantity ▲

- (2) **506G3=M6×10** Set screws
- (3) **506G3=M6×12** Set screws
- (4) **501Z2=M4×20** Head screw

3 Safety Instructions

- Please read the relevant material safety data sheet of the lamination resin and hardener before laminating.
- Use torque wrench for tightening the adjustment screws and tube clamp and observe the torque values!

4 Description

4.1 Intended Use

The 4R110 Lamination Anchor is suitable for both transtibial and transfemoral prostheses and is to be integrated into the prosthetic socket during lamination. During trial fitting the adapter can be screwed onto the thermoplastic test socket. The enclosed protective cap protects the mechanical components during the lamination process and allows for an exact fit on the distal end of the laminated socket.

4.2 Field of Application

The 4R110 Lamination Anchor is to be used **solely** for the exoprosthetic fitting of unilateral or bilateral lower limb amputations for children with body weight up to **45 kg (99 lbs)**.

Attention!

Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water, and acids. Using this medical product under such conditions will render all claims against Otto Bock Healthcare Products GmbH null and void.

Please also inform your patients!

These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

5 Technical Data

	4R110
Weight	55 g
System Height	35 mm
Material	Aluminium
Patient weight max.	45 kg (99 lbs)
Mobility grade	–

6 Symbols used

Medical device

7 Processing Instructions

7.1 Laminating

Attention!

Please read the relevant material safety data sheet of the lamination resin and hardener before laminating.

- Seal the clamp screw, adjustment ring and clamp shell using 633G6 Special Lubricant (→, Fig. 2). Carefully attach the cover cap.

Attention!

Do not modify any part of the lamination dummy! It protects the mechanical components during the lamination process and allows for an exact fit on the distal end of the laminated socket. The pegs on the cover must not be removed. They serve as spacers.

Attention!

Position the lamination anchor so that the clamp shell is positioned medially when viewed from the top (Fig. 3). If incorrectly positioned, the support for the lamination anchor is not provided. This would result in an overloading of the connection adapter, and could damage the adapter.

Attention!

When adjusting the adapter make sure that the adjustment ring is screwed in and tighten the clamping screw to **3 Nm**.

- The layup is based on the expected loading and consists of nyloglas stockinette with additional glass or carbon fiber reinforcement as necessary. The individual stockinette layers are tied off in the groove between the lamination cap and lamination anchor (Fig. 4).
- After the resin hardens, break off the cover cap and remove any excess plastic.

7.2 Adjustments and Final Assembly

Attention!

The fitting of a patient with the lamination anchor may only be carried out by prosthetists.

Attention

The rotation adjustment must not be used for length adjustment. The adjustment range is limited to about 45° by safety marking (center punch) on the top of the adjustment ring. Turn the adjustment ring to the desired position and clamp with the clamping screw and torque **3 Nm**. Use 710D12 Torque Wrench.

Attention!

In the case of pronounced adjustment, replace the set screws that are deeply inserted with 506G3=M6×12 Set Screws enclosed.

Static corrections are possible at any time during alignment and trial walking as well as after completion of the prosthesis using the four 506G3 Set Screws. To exchange a module or during disassembly the previous position may be retained by loosening only the two most deeply screwed in, adjacent set screws.

Replace set screws that appear too long (to prevent foam cover damage) or too short (stability) with the appropriate screw length.

Adjustment screws are available in the following sizes:

<i>Article number</i>	<i>Length</i>
-----------------------	---------------

506G3=M6×10	10 mm
--------------------	-------

506G3=M6×12	12 mm
--------------------	-------

Use 710D12 Torque Wrench to tighten the set screws.

Tighten the 506G3 Set Screws clockwise to **7 Nm**.

To finish the prosthesis: Secure the 506G3 Set Screws with 636K13 Loctite.

8 Cleaning and Care

Clean the product with a damp cloth and mild soap (e.g. Ottobock 453H10=1 Derma Clean) when needed.

Ensure that no liquid penetrates into the system component(s).
Dry the product with a lint-free cloth and allow it to air dry fully.

9 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this product may only be used in combination with components that were authorized by the manufacturer (see instructions for use and catalogs). The manufacturer does not assume liability for damage caused by component combinations which it did not authorize.

10 Trademarks

All denotations within this accompanying document are subject to the provisions of the respective applicable trademark laws and the rights of the respective owners, with no restrictions.

All brands, trade names or company names may be registered trademarks and are subject to the rights of the respective owners.

Should trademarks in this accompanying document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

11 CE Conformity

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with European requirements for medical devices. The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

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