
3R32, 3R23

DE Gebrauchsanweisung	5
EN Instructions for use	11
FR Instructions d'utilisation	17
IT Istruzioni per l'uso	23
ES Instrucciones de uso	29
PT Manual de utilização	35
NL Gebruiksaanwijzing	41
SV Bruksanvisning	47
DA Brugsanvisning	53
PL Instrukcja użytkowania	58
HU Használati utasítás	64
CS Návod k použití	69
RO Instrucțiuni de utilizare	75
HR Upute za uporabu	82
TR Kullanma talimatı	87
EL Οδηγίες χρήσης	93
RU Руководство по применению	100
JA 取扱説明書	106
ZH 使用说明书	112





● Einzelteile-Pack 4D9
 Single-Component Pack 4D9
 Kit de pièces de rechange 4D9
 Imballo dei singoli componenti 4D9
 Kit componentes 4D9
 Onderdelenpakket 4D9
 4D9 Service-Set
 Conjunto de peças de reposto 4D9
 Alkatrész-csomag 4D9-hez

▲ Mindestmenge
 Minimum order quantity
 Quantité minimum
 Quantità minima
 Cantidad mínima
 Minimikvantitet
 Minimum cantal
 Quantidade minima
 Minimum mennyiség

■ Einzelteile
 Single components
 Pièces à l'unité
 Singoli componenti
 Componentes
 Separata delar
 Onderdelen
 Componentes
 Alkatrészek

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Last update: 2020-10-07

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

1 Scope of delivery (Fig.1)

1.1 4D9 Single Component Pack (Fig. 1 ●)

Consisting of: 2 stops (15), 4 truss head screws (17), 4 two-hole nuts (18), 4 slotted bushings (9), 4 Belleville spring washers (7), 3 set screws (12, 12a), 2 lock rings (13), 2 lock nuts (14)

1.2 Minimum order quantity ▲

502R1=M5x16 Two-hole Nut (18), 501T1=M5x16 Truss Head Screw (17), 506G3=M8x12-V Set Screw (12), 506G3=M8x10 Set Screw (12a), 502S19=M6 Lock Nut (14), 4B82 Slotted Bushing (9), 513T4=18.8x10.2 Belleville Spring Washer (7), 507S16=10x1 Lock Ring (13), 4Z43=H Stop (15), 507U16=6.4 Spherical Washer (10), 501T15=M6x40 Socket Head Cap Screw (11), 4B86 Axle Bushing (8)

1.3 Single component ■

4G70 Lamination Anchor (1), 4F46 Lock, complete (19)

2 Description

2.1 Indications for use

The polycentric Modular Knee Joints 3R32 (TITAN) and 3R23 (INOX) are to be used **exclusively** for the prosthetic fitting of the lower extremities.

2.2 Field of application

Field of application according to the **Ottobock MOBIS Mobility System**:

Recommended for amputees **with Mobility Grade 1** (indoor walkers).

Approved for a patient weight of up to **125 kg**.

Given the lock, the joints are suitable in particular for patients who require a high degree of safety.

Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water, acids and other liquids. Using this medical device under such conditions will render all claims against Otto Bock HealthCare null and void.

Please inform the patient.

Reuse on another patient

Fall due to loss of functionality as well as damage to the product

Only use the product for a single patient.

Inform the patient.

2.3 Lifetime

This prosthetic component has been tested by the manufacturer for 3 million load cycles according to ISO 10328. The maximum lifetime is 5 years.

2.4 Construction and function

The upper joint section (3) with lamination anchor and lower joint section (2) are connected to one another by axle brackets (4+5) and linkage bars (6). Due to the way in which the joints move, the stretch of the cosmetic foam cover decreases in the sitting position. The lamination anchor (1) is attached to the residual limb bedding and constitutes the adjustable connection piece to the coupling core. The centre of rotation (instantaneous centre of rotation) therewith changes its position in relation to the flexion. In an extended position, the centre of rotation lies at the height of the femoral condyles and moves downward and forward with increasing flexion.

The lock made of plastic is attached between the axle brackets (4+5) with the anchor bracket (26) and secures the joint in extended position. To initiate knee flexion, e.g. for sitting down, the joint is released with the lock control handle (22). In extended position, the lock engages automatically. The lock cable should run along the lateral prosthesis side.

The lock (19) can be replaced with the extension assist of the Joint 3R21 or 3R30 of the same production series.

3 Technical data

Article number	3R23	3R32
Proximal connection	Lamination anchor	
Distal connection	Adjustment pyramid	
Knee flexion angle	110°	
Weight	880 g	655 g
System height	99 mm	
System height from alignment reference point	17 mm	
System height until alignment reference point	82 mm	
Max. patient weight	125 kg / 275 lbs	
Mobility level	1	

4 Handling

4.1 Alignment

The three-dimensional arrangement of the prosthetic socket and the modular components affects the static and dynamic functions of the prosthesis. Function of the joint is influenced by the position of the axes. In extension, the instantaneous centre of rotation is located above the pyramid adapter and behind the alignment reference line, thus providing knee stability in the stance phase. The advantages of the 3R23 and 3R32 Knee Joints can only be made optimal use of in case of a correct alignment.

The optimal residual limb position must be anticipated when positioning the socket connector. Plumb lines in the frontal and sagittal planes (drawn from the hip joint's centre of rotation and marked during plaster cast taking and trial fitting of the test socket) will facilitate correct positioning of the lamination anchor or socket adapter.

To align the prosthesis please proceed in two steps:

1. First make the bench alignment using an alignment tool such as 743L200 L.A.S.A.R. Assembly. For precise clamping of the knee joint in the 743L200 L.A.S.A.R. Assembly we recommend using the special devices 743Y89=1 and 2.
2. Then use 743L100 L.A.S.A.R. Posture for static alignment optimisation.

4.1.1 Bench alignment with alignment tool (the following steps refer to Fig. 6)

- ① Position the middle of the foot **30 mm anterior** to the alignment reference line.
- ② Add 5 mm to the required heel height of the foot. Set correct outward rotation of the foot.
- ③ Clamp the knee joint using the appropriate adapter inserts. For bench alignment, the **alignment reference line should run 0 mm in front of the upper anterior axis (alignment reference point)**. At that point the pyramid base should be horizontal. Pay attention to the knee ground distance and outward rotation of the knee (adapter inserts provide for a rotation of approx. 5°). Recommended positioning of the alignment reference point: 20 mm above the medial tibial plateau.
- ④ Connect the foot to the modular knee joint using a tube adapter.
- ⑤ Mark the centre of the socket proximally and distally on the lateral side. Draw a line through both marks from socket brim to the distal end of the socket.
- ⑥ Now position the socket such that the alignment reference line passes through the proximal centre mark of the socket. Set the socket flexion to somewhere between 3° and 5°; however, the individual situation (e.g. hip joint contractures) must be taken into account and, if necessary, more flexion should be provided. Also pay attention to the ischial tuberosity to ground distance.
- ⑦ Connect the socket and modular knee joint using a corresponding lamination anchor (see Section 4.3.1.2).

4.1.2 Static alignment optimization using L.A.S.A.R. Posture 743L100

(the following steps refer to Fig. 7)

The bench alignment can be substantially improved using L.A.S.A.R. Posture. In order to ensure appropriate stability combined with easy swing phase initiation, please proceed as follows:

- ① To make the load line visible, the transfemoral amputee stands on the L.A.S.A.R. Posture with the prosthetic side on the force plate and with the other leg on the height compensation panel. The prosthesis side should be sufficiently loaded (> 35 % of the body weight).
- ② Now adapt the alignment by only **adjusting the plantar flexion** of the foot. The **load line** (laser line) should be approx. **35 mm anterior to the front knee axis** (see Fig. 7).
- ③ After step 2, dynamic optimisation can take place between parallel bars.

4.2 Combination possibilities

INFORMATION

In a prosthesis, all prosthetic components have to meet the patient's requirements regarding the amputation level, body weight, activity level, environmental conditions and field of application.

4.3 Adjustments and final assembly

4.3.1 Socket fabrication

4.3.1.1 Lamination prior to trial fitting

Take a piece of 623T3 Perlon Stockinette twice the length of the plaster model and pull half of it over the isolated socket. Twist the other half of the stockinette at the distal end and pull it over the socket as well. To strengthen the modular knee joint, the socket must be incrementally reinforced with a 616G13 Fiberglass Stockinette. To do so, pull a first layer of fiberglass stockinette two-thirds over the socket length. Tie off and pull remaining stockinette halfway back over the socket. Place 2 layers of 616G12 Carbon Fiber Cloth on the distal area such that the to-be-mounted 4G70 Adapter will have a carbon fiber basis with a 3 cm rim. Pull over two layers of 623T3 Perlon Stockinette. The socket will be laminated twice with two different resins. For the first casting, the socket is laminated up to two-thirds of its length with 617H19 Orthocryl Lamination Resin. Once the casting has hardened, apply 2 more layers of 623T3 Perlon Stockinette. For the second casting, the proximal part of the socket is laminated with 617H17 Orthocryl Flexible Resin.

The lamination anchor serves to restrict flexion by stopping at the axle bushing of the lower joint section (2) (Fig. 2). This stop function must also be taken into consideration after the final lamination. To ensure a proper flexion stop and to prevent damage, the set screw (12), which may be too long for certain adjustment positions, should be replaced by a shorter one (12a).

4.3.1.2 Applying the lamination anchor

The correctly fitted lamination anchor must be prepared before it is glued to the socket. To do so, fill the coupling opening with 636K8 Plastaband (Fig. 3). Clamp the socket and knee joint into the alignment apparatus.

Attention!

Remove any soft, flexible or porous materials from underneath the supporting surface of the lamination anchor. **Mix the putty exclusively from 617H21 Orthocryl Sealing Resin and 639A1 Talcum Powder. Glue on the lamination anchor.**

Secure with 627B2 Adhesive Tape for trial fitting. Test the stop function. If necessary, form the required contact surface with putty. If need be, glue a stop made of Pedilin onto the outer laminate.

4.3.1.3 Finalizing the socket

Once the trial fitting has been completed, fasten the lamination anchor to the socket with truss head screws (17) and two-hole nuts (18). Then laminate. Then proceed with reinforcement as follows: Pull one layer of 623T3 Perlon Stockinette over the entire socket and tie circularly on the top. Apply two layers of 616G12 Carbon Fiber Cloth over the arms of the 4G70 Adapter. Then pull the remaining half of the 623T3 Perlon Stockinette back over the socket. The carbon fiber layers are thus sandwiched between the stockinette layers. Again reinforce the socket in increments with 616G13 Fiberglass Stockinette (see Section 4.3.1.1) Apply two layers of 623T3 Perlon Stockinette. Proceed with lamination as with the first casting.

Attention!

Deviating from the processing instructions or the recommended materials for the lamination of the lamination anchor can result in a loosening or breaking of the adapter.

Attention!

Upon completion of assembly, tighten the set screws (12-506G3=M8×12-V, 12a-506G3=M8×10) with a torque wrench (710D1) to the required assembly torque. **Trial fitting:** Tighten the set screws to **10 Nm**. **Completion:** Secure the 506G3 Set Screws with 636K13 Loctite, pre-tighten to **10 Nm** and then tighten to **15 Nm**.

4.4 Adjustment of the bearing

The axis pins of the linkage bars, encased in bushings, are connected to the upper and lower joint section. The type of bearing and screw connection allows to increase or decrease the friction and to remove play resulting from ordinary wear and tear. For this, attach lock nut (14) using enclosed 709S18 Socket Wrench (16). Then tighten socket head cap screw (11) with 5-mm Allen wrench (Fig. 4). It is of utmost importance that an adjustment to remove play take place after a breaking-in period of 2 to 4 weeks. At every maintenance check, the adjustment must be verified and the 4Z43=H Stops replaced.

To lubricate the axle bushings (8) and slotted bushings (9) use 633G6 Ottobock Special Lubricant.

4.5 Adjusting the lock (Fig. 5)

The Modular Knee Joint 3R23 or 3R32 is delivered checked and adjusted for assembly. The lock easily and audibly engages automatically when the joint is in extension. It can be readjusted by turning the countersunk screw with hexagon socket (24) above the brass wedge (23). Turning the screw to the **right** = impedes release, reduces play, engages less easily. Turning the screw to the **left** = facilitates release, increases play, engages more easily.

Following any adjustment, ensure that the lock engages and that the joint is reliably locked. Secure the adjustment screw (24) with 636K13 Loctite. At every maintenance check, the lock should be verified and the 4Z43=H Stops (15) replaced. If the locks fail to engage securely, the stops must be reworked. Do not alter the plastic parts of the lock.

5 Cosmetic foam cover

Use Foam Cover 6R6 (3R6, 3S106) for Joint 3R23 and 3R32.

Attention!

Do not use talcum powder to eliminate noise in the cosmetic foam cover. Talcum draws out grease from the mechanical components. This causes considerable mechanical malfunctions and can result in blockage of the knee joint, in turn leading to the fall of the patient. Therefore, using talcum with this medical device will render any warranty claims against the manufacturer null and void.

Note

To reduce friction and to eliminate noise, apply 519L5 Silicone Spray directly onto the contact surfaces of the cosmetic foam cover.

6 Maintenance instructions

Attention!

Do not use any aggressive cleaning agents since they could cause damage to bearings, seals and plastic parts.

Do not disassemble the knee joint! If you have a problem, please send the joint to Ottobock.

Attention – Please inform your patients!

The knee joint function can be altered by environmental conditions and type of use. To avoid risk to the patient, the knee joint should not be used after noticeable functional changes have taken place. Noticeable functional changes may include, among others, poor response, incomplete extension, diminished swing phase control or stance phase stability and abnormal noises.

Action to be taken:

Have the prosthesis examined at a specialist workshop.

Ottobock recommends readjusting the knee joint settings once the patient has gotten used to the prosthesis. The duration of this breaking-in period may differ from patient to patient. Check the knee joint for wear and proper function at least once a year and make adjustments if necessary. When doing so, pay special attention to the resistance to movement and to the development of unusual noise. A complete flexion and extension must be ensured.

7 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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