## BionicLink 60X3

<b>DE</b> Gebrauchsanweisung (Fachpersonal) .....	3
<b>EN</b> Instructions for use (qualified personnel) .....	15
<b>FR</b> Instructions d'utilisation (Personnel spécialisé) .....	29
<b>IT</b> Istruzioni per l'uso (Personale tecnico specializzato) .....	43
<b>ES</b> Instrucciones de uso (Personal técnico especializado) .....	55
<b>PT</b> Manual de utilização (Pessoal técnico) .....	67
<b>NL</b> Gebruiksaanwijzing (Vakmensen) .....	79
<b>SV</b> Bruksanvisning (Fackpersonal) .....	91
<b>NO</b> Bruksanvisning (Fagpersonell) .....	103
<b>PL</b> Instrukcja użytkowania (Personel fachowy) .....	115
<b>CS</b> Návod k použití (Odborný personál) .....	127
<b>EL</b> Οδηγίες χρήσης (Τεχνικό προσωπικό) .....	139
<b>RU</b> Руководство по применению (Квалифицированный персонал) .....	153
<b>JA</b> 取扱説明書 (有資格担当者) .....	167
<b>ZH</b> 使用说明书 (专业人员) .....	179

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<b>1</b>	<b>Foreword .....</b>	<b>16</b>
<b>2</b>	<b>Product description .....</b>	<b>16</b>
2.1	Function .....	16
2.2	Combination possibilities.....	16
<b>3</b>	<b>Intended use .....</b>	<b>16</b>
3.1	Indications for use .....	16
3.2	Conditions of use .....	16
3.3	Contraindications .....	16
3.4	Qualification .....	16
<b>4</b>	<b>Safety .....</b>	<b>17</b>
4.1	Explanation of warning symbols.....	17
4.2	Structure of the safety instructions.....	17
4.3	General safety instructions.....	17
<b>5</b>	<b>Scope of Delivery and Accessories.....</b>	<b>19</b>
5.1	Scope of delivery.....	19
5.2	Accessories .....	19
<b>6</b>	<b>Preparing the product for use.....</b>	<b>19</b>
6.1	Preparing the BionicLink .....	19
6.2	Attaching the BionicLink.....	19
6.3	Removing the BionicLink .....	20
<b>7</b>	<b>Use .....</b>	<b>20</b>
<b>8</b>	<b>Legal information .....</b>	<b>20</b>
8.1	Liability .....	20
8.2	Trademarks .....	20
8.3	CE conformity .....	20
8.4	Local Legal Information .....	21
<b>9</b>	<b>Technical data .....</b>	<b>22</b>
<b>10</b>	<b>Appendices .....</b>	<b>22</b>
10.1	Symbols Used.....	22
10.2	Directives and manufacturer's declaration .....	23

# 1 Foreword

## INFORMATION

Date of last update: 2020-12-23

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

The product "60X3 BionicLink" is referred to as the product/BionicLink below.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

## 2 Product description

### 2.1 Function

The BionicLink provides wireless data communication between the knee joint and a computer with 60X5 BionicLink PC using Bluetooth technology.

Patient-specific adjustments can therefore be carried out without a cable connection, and the patient has more freedom of movement during the adjustment process. The BionicLink must be removed from the patient's prosthesis after adjustment.

The BionicLink is suitable for patients with a left-sided or right-sided fitting.

## INFORMATION

The 60X3 BionicLink does not require any batteries for operation, as it is powered by the battery of the prosthesis.

### 2.2 Combination possibilities

- 4X181-2 C-Leg wireless adapter cable

## 3 Intended use

### 3.1 Indications for use

The product is to be used **exclusively** for data communication for adjusting a C-Leg or C-Leg compact prosthesis system.

### 3.2 Conditions of use

Permissible ambient conditions are described in the technical data (see page 22).

### 3.3 Contraindications

- All conditions which contradict or go beyond the specifications listed in the section on "Safety" and "Indications for use".

### 3.4 Qualification

The fitting of a patient with the product may only be carried out by an O&P professional.

## 4 Safety

### 4.1 Explanation of warning symbols

**NOTICE**

Warning regarding possible technical damage.

### 4.2 Structure of the safety instructions

**NOTICE****The heading describes the source and/or the type of hazard**

The introduction describes the consequences in case of failure to observe the safety instructions. Consequences are presented as follows if more than one consequence is possible:

- > E.g.: Consequence 1 in the event of failure to observe the hazard
- > E.g.: Consequence 2 in the event of failure to observe the hazard
- ▶ This symbol identifies activities/actions that must be observed/carried out in order to avert the hazard.

### 4.3 General safety instructions

**NOTICE****Non-observance of safety notices**

Damage to the product due to usage of the product in specific situations.

- ▶ Observe the safety notices and the stated precautions in this accompanying document.

**NOTICE****Use by untrained personnel**

Limited function due to product malfunction or damage.

- ▶ The product may only be used by trained orthopaedic technicians.

**NOTICE****Penetration of dirt and humidity into the product**

Limited function due to product malfunction or damage.

- ▶ Ensure that neither solid particles nor liquids can penetrate into the product.
- ▶ If liquid has penetrated into the product, forward it to an authorised Ottobock Service Centre before using it again.

**NOTICE****Mechanical stress on the product**

Limited function due to product malfunction or damage.

- ▶ Do not subject the product to mechanical vibrations or impacts.
- ▶ Check the product for visible damage before each use.

**NOTICE****Independent user changes or modifications made to the product**

Limited function due to product malfunction or damage.

- ▶ Manipulations to the product other than the tasks described in these instructions for use are not permitted.
- ▶ The product and any damaged components may only be opened and repaired by authorised, qualified Ottobock personnel.

#### NOTICE

##### **Use of the product with a damaged wireless adapter cable**

Malfunction due to incorrect or interrupted data transfer.

- ▶ Avoid bending or squeezing the wireless adapter cable.
- ▶ Replace defective or damaged wireless adapter cables immediately.
- ▶ Do not make any unauthorised changes to the wireless adapter cable.

#### NOTICE

##### **Operating the product outside of the permissible temperature range**

Limited function due to product malfunction or damage.

- ▶ Avoid operation in areas outside the permissible temperature range (see page 22).

#### NOTICE

##### **Operating near sources of strong magnetic or electrical interference (e.g. theft prevention systems, metal detectors)**

Limited function due to product malfunction.

- ▶ Avoid operating the product in the vicinity of sources of strong magnetic and electrical interference (e.g. high-voltage lines, transmitters, transformer stations, computer tomographs, magnetic resonance tomographs, etc.).
- ▶ Use the product only in orthopaedic workshops, or during in-patient stays in rehabilitation facilities.

#### NOTICE

##### **Insufficient distance to HF communication devices (e.g. mobile phones, Bluetooth devices, WiFi devices)**

Product malfunction due to interference during data transfer.

- ▶ Therefore, keeping the following minimum distances to these HF communication devices is recommended:
  - Mobile phone GSM 850 / GSM 900: 0.99 m
  - Mobile phone GSM 1800 / GSM 1900 / UMTS: 0.7 m
  - DECT cordless phones incl. base station: 0.35 m
  - WiFi (routers, access points ...): 0.22 m
  - Bluetooth devices (third-party products not approved by Ottobock): 0.22 m

#### NOTICE

##### **Operating the product very close to other electronic devices**

Limited function due to product malfunction.

- ▶ Do not stack the product during operation nor operate it in the immediate vicinity of other electronic devices in operation.
- ▶ If simultaneous operation cannot be avoided, monitor the product and verify proper use in the existing setup.

#### NOTICE

##### **Use of unapproved accessories**

- > Limited function due to product malfunction as a result of reduced interference resistance.
- > Interference of other electronic devices due to increased emissions.
- ▶ Use the product only in combination with the accessories, signal converters and cables listed in the section "Combination possibilities" (see page 16).

- ▶ The use of USB extensions is not recommended, since it can lead to disturbances.

#### **INFORMATION**

The BionicLink must be removed from the patient's prosthesis after adjustment.

#### **INFORMATION**

The BionicLink cannot recognise any defects on the prosthesis and prosthetic components by itself. For this reason, attentive observation by the orthopaedic technician and knowledge about the function of the prosthetic components are indispensable.

#### **INFORMATION**

The BionicLink is powered by the battery in the prosthesis. For this reason, adjustments should only be made when the battery in the prosthesis is fully charged.

#### **INFORMATION**

##### **Cybersecurity**

- ▶ Keep your operating system up to date and always install any available security updates.
- ▶ Protect your computer from unauthorised access (e.g., by using virus scans, password protection etc.).
- ▶ Do not use unsecured networks.
- ▶ Please contact the manufacturer if you suspect cybersecurity problems.

## **5 Scope of Delivery and Accessories**

### **5.1 Scope of delivery**

- 1 pc. 60X3 BionicLink
- 1 pc. 4X181-2 C-Leg wireless adapter cable
- 1 pc. hook-and-loop strap
- 1 pc. Instructions for use (qualified personnel)

### **5.2 Accessories**

- "4X180=V2.6 C-Soft" or higher adjustment software
- 60X5 BionicLink PC

## **6 Preparing the product for use**

### **6.1 Preparing the BionicLink**

The following steps only need to be carried out during first-time use.

- 1) Insert the wireless adapter cable plug into the BionicLink.  
**INFORMATION: The plug does not need to be removed again.**
- 2) Insert the proper side of the hook-and-loop strap into the slots on the back of the BionicLink (see fig. 1).

### **6.2 Attaching the BionicLink**

- 1) Open the plug cover on the knee joint.
- 2) Attach the wireless adapter cable plug to the socket on the joint module (see fig. 4).  
**INFORMATION: Observe the plug direction.**  
→ The blue LED on the BionicLink has to light up. If it does not, charge the prosthesis battery or check the connecting cable.
- 3) Engage the BionicLink on the knee joint (see fig. 2)

- 4) Loop the free end of the hook-and-loop strap around the joint module, thread it through the plastic ring and pull it tight (see fig. 3).

#### **INFORMATION**

Apply the BionicLink to the knee joint medially for a fitting on the left side, laterally for a fitting on the right side.

### **6.3 Removing the BionicLink**

Remove the BionicLink and wireless adapter cable following adjustment.

1. Disconnect the wireless adapter cable from the knee joint (see fig. 4). Do not pull on the cable, but pull the plug on the plug housing.
2. Reinstall or close the plug cap.

## **7 Use**

- 1) Launch the adjustment software on the computer (C-Soft version 2.6 or later).
- 2) Start the connection process in C-Soft.
- 3) The adjustments can be carried out.

**INFORMATION: The prosthesis wearer must stand or sit comfortably and calmly during the adjustment.**

→ The blue LED on the BionicLink flashes during the data transfer process.

#### **INFORMATION**

The range of the Bluetooth® connection is around 10 metres, but this can vary depending on the local conditions. The prosthesis should not be taken outside of this range while it is connected.

#### **INFORMATION**

If the prosthesis is outside the transmitter's range or if the transmission path is shadowed (e.g. by steel or reinforced concrete), the connection may be interrupted. This will be indicated by the warning "**connection interrupted**" on the computer.

## **8 Legal information**

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### **8.1 Liability**

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### **8.2 Trademarks**

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

### **8.3 CE conformity**

Otto Bock Healthcare Products GmbH hereby declares that the product is in compliance with applicable European requirements for medical devices.

The product meets the requirements of the RoHS Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic devices.

This product meets the requirements of the 2014/53/EU directive.

The full text of the regulations and requirements is available at the following Internet address: <http://www.ottobock.com/conformity>

#### 8.4 Local Legal Information

Legal information that applies **exclusively** to specific countries is written in the official language of the respective country of use in this chapter.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference, and
- 2) This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**Caution: Exposure to Radio Frequency Radiation.**

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party:

Otto Bock Health Care, LP  
3820 West Great Lakes Drive  
Salt Lake City, Utah 84120-7205 USA  
Phone + 1-801-956-2400  
Fax + 1-801-956-2401

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of this device.

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**Caution: Exposure to Radio Frequency Radiation.**

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website <http://www.hc-sc.gc.ca/rpb>.

Responsible party:  
 Otto Bock Healthcare Canada Ltd.  
 5470 Harvester Road  
 L7L 5N5 Burlington, Ontario  
 Canada  
 Phone + 1-800-665-3327

**Caution: Federal law (USA) restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.**

**9 Technical data**

<b>Ambient conditions</b>	
Storage and transport in original packaging	±0 °C/+32 °F to +60 °C/+140 °F max. 93 % relative humidity, non-condensing
Storage and transport without packaging	±0 °C/+32 °F to +60 °C/+140 °F
Operation	±0 °C/+32 °F to +60 °C/+140 °F

<b>General</b>	
Reference number	60X3 BionicLink
Product service life	5 years
Range of the Bluetooth connection in buildings	approx. 10 m
Bluetooth class	Class 2
Operating voltage	3.6V $\approx$
Power supply	By the prosthesis battery
Recommended adjustment software	C-Soft, version 1.0 or later

<b>Data transfer</b>	
Wireless technology	Bluetooth 2.0
Range	approx. 10 m / 32.8 ft
Frequency range	2402 MHz to 2480 MHz
Modulation	GFSK
Data rate (over the air)	704 kbps
Maximum output power (EIRP)	-1.78 dBm (0.66 mW)

**10 Appendices**

**10.1 Symbols Used**

Manufacturer

Lot number (PPPP YYYY WW)  
 PPPP – plant  
 YYYY – year of manufacture  
 WW – week of manufacture

BT ID: YYWWNNN

Bluetooth ID number

Compliance with the requirements according to "FCC Part 15" (USA)

Declaration of conformity according to the applicable European directives

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Non-ionising radiation

Compliance with the requirements under the "Radiocommunications Act" (AUS)

Medical device

## 10.2 Directives and manufacturer's declaration

Tables according to IEC 60601-1-2 for the product 'BionicLink 60X3'.

**Table 1 – Electromagnetic emissions for all devices and systems**

<b>Directives and manufacturer's declaration – electromagnetic interference</b>		
The product 'BionicLink 60X3' is designed for operation in an electromagnetic environment as specified below. The customer or user of the product 'BionicLink 60X3' must ensure that the device is operated in such an environment.		
<b>Interference measurements</b>	<b>Compliance</b>	<b>Electromagnetic environment directive</b>
HF emissions according to CISPR 11	Group 1	The product 'BionicLink 60X3' uses HF energy exclusively for its internal functioning. Therefore its HF emissions are very low and interference with neighbouring electronic devices is unlikely.
HF emissions according to CISPR 11	Class B	The product 'BionicLink 60X3' is intended for use in other

Interference measurements	Compliance	Electromagnetic environment directive
Harmonics according to IEC 61000-3-2	N/A	than residential facilities and in facilities connected directly to a public low-voltage network that also supplies residential buildings.
Voltage fluctuations/flicker according to IEC 61000-3-3	N/A	

**Table 2 – Electromagnetic interference immunity for all devices and systems**

**Directives and manufacturer's declaration – electromagnetic interference immunity**

The product 'BionicLink 60X3' is designed for operation in an electromagnetic environment as specified below. The customer or user of the product 'BionicLink 60X3' must ensure that the device is operated in such an environment.

Interference immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment directive
Electrostatic discharge (ESD) according to IEC 61000-4-2	± 6 kV discharge on contact ± 8 kV discharge through air	± 6 kV discharge on contact ± 8 kV discharge through air	Floors should be wood, concrete or ceramic tile. If the floor covering is made of synthetic material, the relative humidity must be at least 30 %.
Electrical fast transient/bursts according to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines	N/A	The quality of the supply voltage should be equivalent to a typical commercial or hospital supply.
Surges according to IEC 61000-4-5	± 1 kV phase-to-phase voltage ± 2 kV even voltage	N/A	The quality of the supply voltage should be equivalent to a typical commercial or hospital supply.
Voltage drops, short interruptions and fluctuations of the supply voltage according to IEC 61000-4-11	< 5 % $U_T$ (> 95 % drop of $U_T$ ) for 1/2 period 40 % $U_T$ (60 % drop of $U_T$ ) for 5 periods 70 % $U_T$ (30 % drop of $U_T$ ) for 25 periods < 5 % $U_T$ (> 95 % drop of $U_T$ ) for 5 s	N/A	The quality of the supply voltage should be equivalent to a typical commercial or hospital supply. If the user of the product 'BionicLink 60X3' demands continuous operation even in case of an interruption in the supply voltage, we recommend supplying the product 'BionicLink 60X3' with an uninterruptible power supply or battery.

Interference immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment directive
Magnetic field at mains frequency (50/60 Hz) according to IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields at the mains frequency should be equivalent to the typical levels for commercial or hospital environments.
Note: $U_T$ is the rated voltage before applying the test levels.			

**Table 4 – Electromagnetic interference immunity for non-life-sustaining devices and systems**

**Directives and manufacturer's declaration – electromagnetic interference immunity**

The product 'BionicLink 60X3' is designed for operation in an electromagnetic environment as specified below. The customer or user of the product 'BionicLink 60X3' must ensure that the device is operated in such an environment.

Interference immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment directive
Conducted HF interference according to IEC 61000-4-6	3 $V_{\text{eff}}$ 150 kHz to 80 MHz	Not applicable	Portable and mobile wireless equipment should not be used at
Radiated HF interference according to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	12 V/m 80 MHz to 1 GHz 3 V/m 1 GHz to 2.5 GHz	a lesser distance from the product 'BionicLink 60X3' including the lines than the recommended safety distance calculated using the equation applicable for the transmission frequency. <b>Recommended safety distance:</b> $d = 0.3 \sqrt{P}$ for 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ for 800 MHz to 2.5 GHz where P is the nominal output of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended safety distance in metres (m). The field strength of stationary radio transmitters at all frequencies according to an

Interference immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment directive
			on-site investigation <sup>a</sup> should be less than the compliance level. Interference is possible in the vicinity of devices that bear the following symbol. 

**Note 1:** The higher frequency range applies at 80 MHz and 800 MHz.

**Note 2:** These directives may not be applicable in all cases. The propagation of electromagnetic factors is influenced by absorption and reflection by buildings, objects and people.

<sup>a</sup> The field strength of stationary transmitters such as base stations of radio telephones and mobile land radio equipment, amateur radio stations, AM and FM radio and television stations cannot be precisely determined theoretically in advance. A site survey is recommended to establish the electromagnetic conditions as a result of stationary HF transmitters. If the field strength determined at the site of the product 'BionicLink 60X3' exceeds the compliance level specified above, the product 'BionicLink 60X3' has to be monitored in regards to normal operation at each place of use. If unusual performance characteristics are noted, taking additional steps may be necessary, for example changing the orientation or moving the product 'BionicLink 60X3'

**Table 6 – Recommended safe distances between portable and mobile HF telecommunication equipment and the product 'BionicLink 60X3' – for non-life-sustaining devices and systems**

**Directives and manufacturer's declaration – Recommended safe distances between portable and mobile HF telecommunication equipment and the product 'BionicLink 60X3'**

The product 'BionicLink 60X3' is designed for operation in an electromagnetic environment where HF interference is monitored. The customer or user of the product 'BionicLink 60X3' can help prevent electromagnetic interference by complying with the minimum distances between wearable and mobile HF communication equipment (transmitters) and the product 'BionicLink 60X3', as below according to the maximum output of the communication equipment.

Nominal output of transmitter [W]	Safety distance [m] by transmission frequency		
	150 kHz to 80 MHz	80 MHz to 800 MHz $d = 0.3 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	N/A	0.03	0.23
0.1	N/A	0.09	0.73
1	N/A	0.3	2.3
10	N/A	0.95	7.3
100	N/A	3	23

For transmitters with no maximum nominal output specified in the table above, the recommended safety distance  $d$  in metres (m) can be determined using the equation in the respective column, where  $P$  stands for the maximum nominal output of the transmitter in watts (W) according to the transmitter manufacturer's information.

**Note 1:** The higher frequency range applies at 80 MHz and 800 MHz.

Nominal output of transmitter [W]	Safety distance [m] by transmission frequency
<b>Note 2:</b> These directives may not be applicable in all cases. The propagation of electromagnetic factors is influenced by absorption and reflection by buildings, objects and people.	

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### **Caution: Exposure to Radio Frequency Radiation.**

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Responsible party:  
Otto Bock Health Care, LP  
3820 West Great Lakes Drive  
Salt Lake City, Utah 84120-7205 USA  
Phone + 1-801-956-2400  
Fax + 1-801-956-2401

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The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website

<http://www.hc-sc.gc.ca/rpb>.

Responsible party:

Otto Bock Healthcare Canada Ltd.  
5470 Harvester Road  
L7L 5N5 Burlington, Ontario  
Canada  
Phone + 1-800-665-3327

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