














12S6=L/R







DE Gebrauchsanweisung	7
EN Instructions for Use	12
FR Instructions d'utilisation	16
IT Istruzioni per l'uso	21
ES Instrucciones de uso	26
PT Manual de utilização	31
NL Gebruiksaanwijzing	35
SV Bruksanvisning	40
DA Brugsanvisning	44
NO Bruksanvisning	49
PL Instrukcja użytkowania	53
CS Návod k použití	58
TR Kullanma talimatı	62
EL Οδηγίες χρήσης	66
RU Руководство по применению	72

INFORMATION

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INFORMATION

Date of the last update: 2020-05-19

- Please read this document carefully before using the product and observe the safety notices.
- Instruct the user in the safe use of the product.
- Please contact the manufacturer if you have questions about the product or in case of problems.
- Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- Please keep this document for your records.

Explanation of Symbols

 **WARNING** Warnings regarding possible risks of severe accident or injury.

 **CAUTION** Warnings regarding possible risks of accident or injury.

 **NOTICE** Warnings regarding possible technical damage.

 **INFORMATION** Additional information on the fitting / use.

The scope of delivery is shown on the front page.

1 Single Components**1.1 Single Components**

Available single components:

13D1 O-Ring Set for 12S8

13D2=L Shoulder Bracket Set for 12S6

13D2=R Shoulder Bracket Set for 12S6

2 Description**2.1 Intended Use**

The Ottobock MovoShoulder Swing 12S6 is to be used **solely** for the exoprosthetic fitting of upper limbs in conjunction with the Ottobock myo arm prostheses or Ottobock passive arm prostheses.

2.2 Field of Application

The Ottobock MovoShoulder Swing 12S6 can be used for prosthetic fitting after unilateral or bilateral shoulder amputations or shoulder girdle amputation.

2.3 Safety Instructions

Please forward the following safety instructions to your patients:

 WARNING

Risk of accident when driving a motor vehicle. The ability of prosthesis users to drive a vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities.

All persons are required to observe their country's national and state driving laws when operating motor vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorized test centre.

For maximum safety and convenience, Ottobock recommend that, at the very least, a specialist evaluate the need for any adaptations to the car (such as by installing a steering fork, automatic shift). Risk-free driving should be ensured even when the prosthesis is not functioning.

CAUTION

Risk of injury as a result of excessive strain. The Ottobock prosthetic components described in these instructions for use have been developed for activities of daily living and must not be used for unusual activities such as extreme sports (free climbing, paragliding, etc.).

Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's safety!

Should prosthetic components be subjected to unusual stresses (such as a fall), immediately contact your prosthetist and have the components inspected for any damage. If necessary, the responsible prosthetist will pass the prosthesis on to the Ottobock Service Centre.

NOTICE

Risk of corrosion. Prosthetic components must not be exposed to corrosive elements such as fresh water, salt water, and acids.

Using this medical product in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

NOTICE

Damage caused by inappropriate environmental conditions. Prosthetic components must not be subjected to intense smoke, dust, vibrations, shocks, or high temperatures. Do not allow debris or liquids to get into the prosthetic components. Failure to follow these instructions can lead to malfunctioning and damage to the prosthesis.

NOTICE

Damage caused by improper cleaning agents. Do not use any aggressive cleaning agents, since they could cause damage to bearings, seals and plastic parts.

2.4 Function (Fig. 1–4)

The free swing behavior of the Ottobock MoveShoulder Swing promotes a natural gait pattern. It also reduces the stress on the skin in the contact area of the prosthetic socket as well as the spine and back muscles.

2.4.1 Locked Position

To lock the shoulder joint, move the prosthesis to approximately 15° of abduction using the sound hand (see Fig. 1). Move the whole prosthesis forward by about 30° and allow it to lock.

The prosthesis can also be positioned without the assistance of the sound hand by proceeding as follows: The patient bows his/her torso laterally (Fig. 2) so as to tilt the prosthesis away from

his body and then leans ahead to swing the prosthesis forward (Fig. 3). Then the patient moves back into an upright position while the prosthesis remains locked in the forward position (Fig. 4). In this position, the patient can perform simple gripping functions in combination with myo-electric components (e.g. System Electric Hand, System Electric Greifer).

The O-rings serve to dampen the tilt behaviour during abduction and prevent the prosthesis from unintentionally tilting away from the body. It is permissible to use a single O-ring resulting in reduced abduction dampening.

INFORMATION

Regularly check the O-rings for any wear and tear and have them replaced by a prosthetist if necessary.

2.4.2 Free Swing Position

To reach the free swing position, perform the reverse sequence of movements from the locked position. The patient leans forward, bows his/her torso laterally, and moves to the upright position. Now the shoulder joint is released for the free swing function.

The swinging behaviour of the prosthesis significantly influences the gait pattern of the patient. It is individually adapted to the patient by adjusting the friction of the shoulder joint.

The free swing friction can be adjusted by turning the cap screw in or out (Fig. 5, arrow). Turning the cap screw in will increase free swing dampening.

INFORMATION

The cap screw (Fig. 5, arrow) must not be removed under any circumstances. Apply Loctite® 241 (Article number 636K13) to the cap screw and turn the screw in until the rotation axis of the shoulder joint exhibits noticeable friction. Make sure that the cap screw never protrudes out of the countersink.

The amount of friction can be adjusted at any time through turning the cap screw. Apply Loctite® 241 (article number 636K13) to the cap screw after each adjustment.

3 Technical Data

Article number	12S6=*
Weight [g]	330
Anteversión, locked [°]	30
Max. anteversión in free swing [°]	28
Max. retroversión in free swing [°]	12
Max. abduction [°]	20

4 Handling

4.1 Assembling the Shoulder Joint (Fig. 6 to Fig. 13)

The shoulder joint is available in versions for left- and right-sided amputees. The shoulder joint is attached to the socket with the enclosed shoulder bracket.

INFORMATION

For patients with subtotal shoulder disarticulation, we recommend mounting the shoulder joint distally and slightly below shoulder height.

NOTICE

Damage of O-ring due to improper assembly. Make sure that there is no contact between the O-ring and socket when mounting the shoulder joint. O-ring friction against the socket may lead to inadmissible wear of the O-ring.

- 1) Mount the shoulder bracket to the shoulder joint with the enclosed oval head screws. Position the mounted shoulder joint on the patient (Fig. 6).
Left body half see Fig. 7, right body half see Fig. 8.
- 2) Adjust the shoulder bracket arm on the socket by bending using the Bending Iron 711S5 or 711S4. See Fig. 9 for an example.
- 3) Mark mounting holes for shoulder joint on the socket. If the \varnothing 1 mm holes are required, they must be bored up to \varnothing 4.2mm.
- 4) Dismantle the shoulder bracket from the shoulder joint.
- 5) Attach the shoulder bracket to the socket (Fig. 10).
- 6) Mark the positions of the four mounting holes on the socket (Fig. 11). Remove the shoulder bracket.
- 7) Bore up the marked locations on the socket with a \varnothing 8 mm drill. The bores are used to access to the oval head screws. This allows the shoulder joint to be dismantled from the shoulder bracket at any time.
- 8) Attach the shoulder bracket to the socket. Attach the shoulder joint to the shoulder bracket with four oval head screws M4 \times 12 and tighten them to a max. torque of 2 Nm. Secure the screws by applying Loctite® 241 (636K13)! (Fig. 12).

INFORMATION

While attaching the shoulder bracket, make sure it is in contact with the socket in the location of the four mounting holes.

- 9) Apply the Ottobock Special Lubricant 633G6 onto the contact surfaces of the O-ring (Fig. 13) and mount the O-ring.
- 10) We recommend covering the shoulder joint with a removable cosmetic to prevent the special lubricant from soiling the patient's clothing.

HINWEIS

Damage of electrode cable due to improper assembly. Damage of electrode cable due to improper assembly. Make sure that the electrode cable is not kinked or pinched. The electrode cable must be installed so that it is guided over the shoulder joint.

4.2 Replacing the O-rings

When the O-rings of the shoulder joint show signs of wear, they must be replaced using the 13D1 O-ring set (Fig. 14) and according to the following instructions:

- 1) Cut through the porous O-rings (Fig. 15) and dispose of them (Fig. 16).

12S6=L/R

- 2) Use Allen wrenches suitable to the screws (Fig. 17).
- 3) If necessary, warm the glued screw connections using a hot air gun to facilitate removal of the previously used adhesive.
- 4) Loosen the four screws (Fig. 18), remove (Fig. 19) and dispose of them.
- 5) Remove the steel bracket by pulling it downwards (Fig. 20).
- 6) Put the new O-rings over the steel bracket (Fig. 21).
- 7) Place the steel bracket on the shoulder bracket (Fig. 22) and pull the O-rings onto the designated position.
- 8) Apply Loctite® 241 (636K13) to the upper half of the large, middle screw only (Fig. 23). Do not apply any adhesive to the lower half of the screw. Apply a thin coat of adhesive to the thread of the other three screws as well (Fig. 23).
- 9) Insert the screws (Fig. 24) and tighten them with an installation torque of 2 Nm (Fig. 25).

5 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

6 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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