
2C1,2C3,2C6

DE	Gebrauchsanweisung	4
EN	Instructions for use	10
FR	Instructions d'utilisation	16
IT	Istruzioni per l'uso	22
ES	Instrucciones de uso	28
PT	Manual de utilização	34
NL	Gebruiksaanwijzing	40
SV	Bruksanvisning	46
DA	Brugsanvisning	51
NO	Bruksanvisning	57
PL	Instrukcja użytkowania	62
HU	Használati utasítás	68
CS	Návod k použití	74
RO	Instructiuni de utilizare	80
HR	Upute za uporabu	86
TR	Kullanma talimatı	91
EL	Οδηγίες χρήσης	97
RU	Руководство по применению	103
JA	取扱説明書	109
ZH	使用说明书	115

①

DE	Einzelteile-Pack	Mindestmenge
EN	Single-Component Pack	Minimum order quantity
FR	Kit de pièces de rechange	Quantité minimum
IT	Imballo dei singoli componenti	Quantità minima
ES	Kit componentes	Cantidad mínima
PT	Conjunto de peças de reposo	Quantidade mínima
NL	Onderdelenpakket	Minimikvantitet
SV	Service-Set	Minimum aantal
DA	Komponentpakke	Mindste mængde
NO	Enkeltkomponentpakke	Minsteordremengde
PL	Zestaw naprawczy	Ilość minimalna
HU	Alkatrész-csomag	Minimum mennyiség
CS	Sada jednotlivých dílů	Minimální množství pro objednání
RO	Set componente individuale	Cantitate minimă
HR	Pakiranje pojedinačnih dijelova	Minimalna količina
TR	Münferit parça ambalajı	Minimum miktar
EL	Σετ μεμονωμένων εξαρτημάτων	Ελάχιστη ποσότητα
RU	Комплект деталей	Минимальное количество
JA	セットで発注いただく部品	複数入パックで発注いただく部品
ZH	维修组件	最低起订量

②

③

④

⑤

⑥

⑦

⑧

INFORMATION

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INFORMATION

Last update: 2017-04-07

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

Explanation of symbols

NOTICE Warnings regarding possible technical damage.

INFORMATION Additional information on the fitting/use.

INFORMATION

Please read the Instructions for Use carefully. Pay special attention to the safety instructions!

1 Single components

1.1 Single component ■ (Fig. 1; for reordering)

(1) **2C19=connection cap for 2C1/2C6 normal footshells** (can be ordered separately)

Ordering example: 2C19=L23-25/4

Sides	Left (L), right (R)			
Sizes [cm]	22	23–25	26–28	29–30
Colors	beige (4), light brown (15)			

(2) 2C10=connection cap for 2C3 normal footshell

Order example: 2C10=L23-25/4

Sides	left (L), right (R)			
Sizes [cm]	21-22	23-25	26-28	29-31
Colours	beige (4), light brown (15)			

(3) 2C20=connection cap for 2C1/2C3/2C6 slim footshells (can be ordered separately)

Ordering example: 2C20=L22/4

Sides	Left (L), right (R)						
Sizes [cm]	21	22	23	24	25	26	27
Colors	beige (4), light brown (15)						

2 Description

2.1 Intended use

The 2C1/2C3/2C6 footshells were developed as individual components for the following Ottobock prosthetic feet:

2C1=* for the 1M10 Adjust

2C3=* for the 1C30 Trias

2C6=* for the Triton family

They are to be used **exclusively** with the corresponding prosthetic foot. The footshells are designed to fit most types of shoes.

2.2 Service Life

The product is a wear part, which means it is subject to normal wear and tear.

2.3 Safety instructions

Please forward the following safety instructions to your patients:

NOTICE

Damage due to incorrect environmental conditions. The footshell material is not UV- or water-resistant.

Do not grind the surface of this prosthetic foot, as this can lead to premature wear.

Using this medical product under such conditions will render all claims against Otto Bock HealthCare null and void.

3 Technical data

2C1=*N footshell, normal

Sides	left (L), right (R)								
Sizes [cm]	22	23	24	25	26	27	28	29	30
Weight [~g]	150	165	165	185	195	230	240	260	280
Colors	beige (4), light brown (15)								
Heel height [mm]	10 +/- 5								

2C1=*S footshell, narrow

Sides	left (L), right (R)				
Sizes [cm]	22	23	24	25	26
Weight [~g]	120	130	135	150	165
Colors	beige (4), light brown (15)				
Heel height [mm]	20 +/- 5				

2C3= * footshell, normal

Sides	left (L), right (R)									
Sizes [cm]	21	22	23	24	25	26	27	28	29	30
Weight [~g]	130	145	155	175	185	205	210	240	245	260
Colours	beige (4), light brown (15)									
Heel height [mm]	10 +/- 5									

2C3= *S footshell, slim

Sides	left (L), right (R)					
Sizes [cm]	21	22	23	24	25	26
Weight [~g]	110	120	130	135	150	165
Colours	beige (4), light brown (15)					
Heel height [mm]	20 +/- 5					

2C6= *footshell, normal

Sides	left (L), right (R)					
Sizes [cm]	24	25	26	27	28	29
Weight [~g]	190	200	220	225	245	270
Colors	beige (4), light brown (15)					
Heel height [mm]	10 +/- 5					

2C6= *S footshell, narrow

Sides	left (L), right (R)					
Sizes [cm]	21	22	23	24	25	26
Weight [~g]	135	145	155	180	190	200
Colors	beige (4), light brown (15)					
Heel height [mm]	15 +/- 5					

4 Handling

NOTICE

2C1=*: Damage caused by improper footshell assembly. Please do not use any tool to insert/remove the foot into/from the footshell as this can cause damage to the elements of the foot.

NOTICE

Risk of damage as a result of improper use. The footshell must not be dusted with talcum powder inside as this would increase abrasion. Replace the footshell if it becomes damaged or worn.

4.1 Inserting the foot into the footshell

Before inserting the foot into the footshell, attach the tube adapter to the prosthetic foot first. This will facilitate the procedure.

2C1=*

Slide the prosthetic foot with your hands into the forefoot section of the footshell (Fig. 2). Then press the foot down in the footshell (Fig. 3).

Make sure that the foot's base spring is located in front of the toe wedge of the footshell (Fig. 4) and that the foot adapter has engaged below the locking wedge (Fig. 5).

2C3=*/2C6=*

Guide the prosthetic foot into the footshell by hand and position the toe of the foot in the toe area of the footshell.

Hold the foot in the toe area of the footshell. Ensure that the foot heel spring has engaged below the heel wedges of the footshell.

4.2 Removing the footshell

2C1=*

Push the locking of the footshell to the rear (Fig. 6) and pull up the heel section of the prosthetic foot (Fig. 7). Now you can remove/pull the prosthetic foot out of the footshell completely (Fig. 8).

2C3=*/2C6=*

Push back on the toe area of the footshell and lever the heel out of the shell. Now the prosthetic foot can be removed completely from the footshell.

4.3 Care and use instructions

Please inform your patients.

Clean the footshell with a mild soap solution.

Wearing the footshell always in the foot increases the lifespan of the footshell. Advise the patient not to walk barefoot.

5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

5.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH