ottobock.

1WR95 Aqua-Fuß

DE	Gebrauchsanweisung (Fachpersonal)	3
EN	Instructions for use (Qualified Personnel)	10
FR	Instructions d'utilisation (Personnel spécialisé).....	17
IT	Istruzioni per l'uso (Personale tecnico specializzato)	24
ES	Instrucciones de uso (Personal técnico especializado)	31
PT	Manual de utilização (Pessoal técnico)	38
NL	Gebruiksaanwijzing (Vakmensen).....	45
SV	Bruksanvisning (Fackpersonal).....	52
DA	Brugsanvisning (Faguddannet personale)	58
NO	Bruksanvisning (Fagpersonell)	65
PL	Instrukcja użytkowania (Personel fachowy).....	71
HU	Használati utasítás (szakszemélyzet).....	78
CS	Návod k použití (Odborný personál)	85
RO	Instrucțiuni de utilizare (Personal de specialitate).....	92
HR	Upute za uporabu (Stručno osoblje)	99
TR	Kullanma talimatı (Uzman personel).....	106
EL	Οδηγίες χρήσης (Τεχνικό προσωπικό).....	113
RU	Руководство по применению (Квалифицированный персонал)	121
JA	取扱説明書 (有資格担当者)	129
ZH	使用说明书 (专业人员).....	136

INFORMATION

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INFORMATION

Last update: 2016-02-10

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

1.1 Function

The functional characteristics of the 1WR95 Aqua Foot are achieved through a specially designed rollover element and differentiated foam technology.

Anti-slip/waterproof properties

The 1WR95 Aqua Foot has particularly good anti-slip properties and is distinguished through good traction. The modular foot version of the 1WR95 Aqua Foot is sealed against water at the factory. The screw access on the sole is sealed with a plug.

Cosmetic appearance

The 1WR95 Aqua Foot has a natural shape, defined toes, and a separate big toe.

2 Application

2.1 Indications for use

The 1WR95 Aqua Foot is intended exclusively for the exoprosthetic fitting of the lower extremities.

2.2 Area of application

The 1WR95 Aqua Foot is intended exclusively for use in a waterproof prosthesis used as a secondary prosthesis.

The 1WR95 Aqua Foot is water resistant.

1WR95=L/R-0-P/4

The 1WR95=L/R-0-P/4 Aqua Foot is designed for use in modular prostheses and includes an installed titanium modular adapter.

It may only be used in combination with Ottobock components for water-proof walking devices.

1WR95=L/R-0-W/4

CAUTION

Risk of falls due to incorrect use

The 1WR95=L/R-0-W/4 is intended exclusively for prosthesis alignment using an exoskeletal design. The foot may not be combined with a modular foot adapter.

2.3 Conditions of use

The 1WR95 is approved for a body weight up to 150 kg.

3 Safety

3.1 Explanation of the ANSI symbols

 CAUTION Warnings regarding possible risks of accident or injury.

 NOTICE Warnings regarding possible technical damage.

 INFORMATION Additional information on the fitting / use.

3.2 General safety instructions

CAUTION

Non-observance of safety instructions

Non-observance of the following safety instructions can lead to malfunction of the 1WR95 Aqua Foot, resulting in a risk of injury to the patient and in destruction of the product.

3.3 Fitting safety instructions

CAUTION

Incorrect bonding of the moulded ankle part and open access to the screw in the sole may result in material damage and loss of functionality of the product.

Risk of a loss of functionality if water gets into the 1WR95=L/R-0-W/4

- Install the included plug with 636W17 adhesive in order to close the screw access in the sole of the foot and prevent water from entering.

3.4 Safety instructions for patients/users

Please make your patients aware of the following safety instructions:

⚠ CAUTION

Risk of falls due to incorrect use

The waterproof walking device is not intended to replace the initial prosthesis and should only be used in wet areas. It may not be used for extreme activities such as extreme sports (water skiing, kite surfing, etc.).

⚠ CAUTION

Risk of falls due to incorrect use

Walking speed must be adjusted accordingly in wet areas and in areas where a higher risk of falls can be expected.

⚠ CAUTION

Risk of burns in cases of heat exposure

The Ottobock components for waterproof walking devices are not suitable for use in saunas or for activities during which the prosthesis is constantly exposed to high heat above 50 °C.

NOTICE

Risk of loss of function due to incorrect storage/transport

The waterproof walking device may only be stored or transported at temperatures between -10 °C and +50 °C.

NOTICE

Risk of loss of function due to incorrect care

The waterproof walking device must be cleaned with clean fresh water every time after contact with saltwater or water containing chlorine/soap or contact with abrasive substances (e.g. sand). To clean the prosthesis, it should be placed in a sink (tub) with clean freshwater and moved around to rinse the corners, edges, and the tube.

The prosthesis must also be dried immediately after every use.

Contact with abrasive substance (e.g. sand) can lead to premature wear.

4 Scope of delivery

The scope of delivery is shown on the cover page and includes either the modular version with installed plug or the exoskeletal version with included plug.

5 Preparation for use

INFORMATION

According to the Ottobock alignment recommendations, bench alignment of the waterproof walking device with the 1WR95 Aqua Foot on 0 mm heel height is usually only possible with the aid of the angled tube adapter. Therefore you should use the 2WR95=1 tube adapter which is angled at the top. If necessary, the 2WR95 straight tube adapter can be used in the thigh area for TF amputees.

5.1 Bench alignment

Align TF modular prostheses in compliance with the „Alignment Recommendations for TF Modular Lower Limb Prostheses according to MOBIS®“ (order no.: 646F219):

- 1) Move the mid-foot forward 30 mm relative to the alignment reference line.
- 2) Set the effective heel height of the foot and add 5 mm.
- 3) Set the outward rotation of the foot.

Follow the instructions in the respective knee joint's instructions for use for the remainder of the bench alignment.

5.2 Static alignment

Significant optimisation of the bench alignment is possible using the L.A.S.A.R. Posture.

5.2.1 TT prostheses

The following alignment is used to achieve the highest possible degree of physiological knee stability for TT prostheses.

Positioning the load line in the sagittal plane:

The foot must be positioned so that the load line runs 15 mm in front of the knee's compromise pivot point.

5.2.2 TF Modular prostheses

Follow the „Alignment Recommendations for TF Modular Lower Limb Prostheses according to MOBIS®“ (order no.: 646F219) and the instructions on bench alignment and static alignment contained in the instructions for use of the 3WR95 Aqua Knee.

5.3 Dynamic alignment optimisation

Dynamic optimisation must be performed during trial walking.

5.4 Adjustment and final assembly

Please observe the following additional steps during adjustment and final assembly of the 1WR95=L/R-0-W/4 Aqua Foot:

- 1) To bond and seal at the same time, coat the glue surfaces of the foot and the 2K34 molded ankle part twice with ORTHOCRYL sealing resin.
- 2) Screw the foot and the moulded ankle part together immediately after applying the second coat.
- 3) Remove any sealing resin which escapes from the glued joint immediately.
- 4) Use 636W17 adhesive to install the included plug in order to seal the screw access hole in the sole of the foot, making sure to push the plug all the way in.

6 Maintenance / repair

INFORMATION

The fundamental design of all Ottobock foot components allows for a possible usage duration of 2 to 3 years, depending on the patient's activity level.

We recommend always performing regular yearly safety inspections.

The safety inspections must be performed by a prosthetist.

6.1 Spare parts

Single components which can be ordered separately

Single components

None

Minimum order quantity

None

Single-component pack

None

7 Disposal

7.1 Disposal information

INFORMATION

These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

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Otto Bock HealthCare GmbH