4R69=5


**2R50
2R49**


4R69=5, 2R50, 2R49

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Information

Date of the last update: 2013-03-20

- Please read this document carefully.
- Follow the safety instructions.

Note: The scope of delivery is shown on the front page.

1 Single Components (fig. 1)

1.1 Single Component ■

–

1.2 Single-Component Pack ●

–

1.3 Minimum Order Quantity ▲

- (1) 506G3=M8×14 Set Screw
- (2) 501Z2=M6×25 Cap Screw

2 Description

2.1 Intended Use

The structural parts shown on the front page are to be **solely** used for the prosthetic fitting of the lower limb.

2.2 Field of Application

Field of application according to the **Ottobock MOBIS Mobility System**:

2R49
2R50

Approved for a **patient weight of up to 100 kg**
(220 lbs)

4R69=5

Approved for a **patient weight of up to 125 kg**
(275 lbs)

Attention!

Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water, and acids. Using this medical product under such conditions will render all claims against Otto Bock HealthCare null and void.

Please also inform your patients!

2.3 Function

The tube adapters are used in combination with the tube clamp adapters as connection elements for modular lower extremity prostheses. These adapter combinations allow for controlled angle and length adaptation in the sagittal and frontal plane as well as adjustment of inward and outward rotation.

Attention!

To achieve maximum stability for the transtibial prostheses, we recommend using a Ø 34 mm tube adapter (e. g. 2R57/2R76) together with a Ø 34 mm tube clamp adapter (e. g. 4R82/4R91). This applies in particular if the prosthesis is to withstand increased pressure due to higher activity levels, longer foot or transtibial lever arms or similar factors.

3 Technical Data

	2R50	2R49	4R69=5
Weight	155 g	240 g	75 g
System height	–	–	33 mm
System height min.	97 mm	97 mm	–
System height max.	232 mm	432 mm	–
Material	Aluminum	Aluminum	Aluminum
Diameter	30 mm	30 mm	30 mm
Patient weight max.	100 kg (220 lbs)	100 kg (220 lbs)	125 kg (275 lbs)
Mobility grade	–	–	–

4 Handling

4.1 Assembly

The tube adapter must be shortened according to the measurements of the amputee.

Attention!

The tube adapter must always be installed distally.

Attention!

Do not clamp the tube into a vice ! It may re-shape or damage the tube.

The section surface must be perpendicular to the tube longitudinal axis (Fig. 2). We recommend using 719R2 Tube Cutter or equivalent.

Carefully debur the tube cutting edge inside and outside. For this purpose use the 718R1 Interior and Exterior Pipe Reamer or equivalent (Fig. 3).

Loosen the cap screw (Fig. 4, item 2) of the tube clamp adapter by two turns and insert the tube 5 mm deep into the tube clamp adapter while slightly rotating the tube. Next, place the tube clamp adapter with the tube inserted vertically on the worktop and fully insert the tube adapter into the tube clamp adapter until the tube is in contact with the base rim inside the clamp.

Attention!

For reasons of structural strength, the residual pieces of cut tube adapters must not be used as components in definitive prostheses.

4.2 Adjustments and Final Assembly

Attention!

Prior to final assembly (i.e. prior to inserting the tube adapter into the tube clamp adapter), the inside of the tube clamp adapter and the outside of the tube adapter must be cleaned with 634A3 (Acetone)!

For final assembly, tighten the cap screw with 710D4 Torque Wrench.

Torque: **13 Nm** (Fig. 4, item 2).

Attention!

The clamping slot of the tube clamp adapter must always be positioned anteriorly.

Static corrections are possible at any time during alignment and trial walking as well as after completion of the prosthesis using the four 506G3 Set Screws. To exchange a module or during disassembly the previous position may be retained by loosening only the two most deeply screwed in, adjacent set screws.

Replace set screws that appear too long (to prevent foam cover damage) or too short (stability) with the appropriate screw length.

Set screws are available in the following sizes.

<i>Article number</i>	<i>Length</i>
-----------------------	---------------

506G3=M8x10	10 mm
--------------------	-------

506G3=M8x12-V	12 mm
----------------------	-------

506G3=M8x14	14 mm
--------------------	-------

506G3=M8x16	16 mm
--------------------	-------

Use 710D4 Torque Wrench to tighten the set screws.

For trial fitting: Tighten the 506G3 Set Screws clockwise to **10 Nm**.

To finish the prosthesis: Secure the 506G3 Set Screws with 636K13 Loctite and pre-tighten to **10 Nm**. Then complete assembly by tightening to **15 Nm** (Fig. 5, item 1).

4.3 Maintenance Instructions

Note:

As a basic principle, all Ottobock modular adapters are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

We recommend carrying out regular safety checks once a year.

5 Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with tested modular components in accordance with the Ottobock Mobility System. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

6 Warranty

Ottobock warrants this product/these products from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction. Because warranty conditions may differ by country, please contact your local Ottobock distribution company for more information.

7 CE conformity

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH