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## EN Manufacturing Instructions for Ottobock System Orthotic Components and Leg Braces –Aluminium–

**Surface processing for light metal**

NOTICE
<p><b>Damage due to improper processing.</b> Excessive heat should not be used when forming the aluminium side bars. To prevent marring the surface, provide bending of these bars with well rounded edges. Use our 711S5 Special Bending Iron. Do not bend the side bars beyond a radius of 10 mm. Nicks and grooves from forming should be removed completely by grinding, brushing or polishing to prevent subsequent breaking in that area.</p>

For surface coating by sintering, the 618T40 Coating Powder is recommended. Heating of the material should not be longer than 5 minutes at a temperature of max. 230 °C.

**Connecting the Side Bars to System Joints**

After adapting the system bars to the system joints, the side bars must be bonded into the connection plugs of the system joints for reasons

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of stability. To do so, degrease the surfaces to be glued with 634A1 thinner or solvent, thoroughly mix the contents of the two cartridges of 636W28 special adhesive and hardener and apply the mixture into the connection plugs. Insert the side bars and secure with 501T7=7,5x9xM5 Phillips Head Screws. Drying time of the glue is 3 to 4 hours. Final strength is achieved after 16 hours.

**Liability**

The manufacturer's warranty applies only if the product has been used under the conditions and for the purposes described. The manufacturer recommends that the product be used and maintained according to the instructions for use.

**CE conformity**

This product meets the requirements of the 93/42/EWG guidelines for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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© Ottobock · 647G78-23-1701
Otto Bock HealthCare GmbH