

AxonEnergy Integral 757B501

DE Gebrauchsanweisung (Fachpersonal)	3
EN Instructions for use (qualified personnel)	10
FR Instructions d'utilisation (Personnel spécialisé)	18
IT Istruzioni per l'uso (Personale tecnico specializzato)	27
ES Instrucciones de uso (Personal técnico especializado)	35
PT Manual de utilização (Pessoal técnico)	43
NL Gebruiksaanwijzing (Vakmensen)	51
SV Bruksanvisning (Fackpersonal)	59
DA Brugsanvisning (Faguddannet personale)	67
NO Bruksanvisning (Fagpersonell)	75
FI Käyttöohje (Ammattihenkilöstö)	82
CS Návod k použití (Odborný personál)	90

INFORMATION

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INFORMATION

Last update: 2015-04-30

- ▶ Please read this document carefully before using the product.
- ▶ Follow the safety instructions to avoid injuries and damage to the product.
- ▶ Instruct the user in the proper and safe use of the product.

- ▶ Please contact the manufacturer if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Contact information can be found on the back page.
- ▶ Please keep this document in a safe place.

These instructions for use provide you with important information on the use, adaptation and handling of the product.

Only put the product into use in accordance with the information contained in the accompanying documents supplied.

2 Application

2.1 Intended Use

The 757B501 AxonEnergy Integral is intended **exclusively** for exoprosthetic fittings of the upper limbs and serves to supply energy to all active Ottobock Axon-Bus prosthesis components.

2.2 Use/Field of Application

The 757B501 AxonEnergy Integral may only be used in conjunction with components of the Ottobock Axon-Bus prosthetic system. **Only** the 757L500 AxonCharge Integral may be used for charging.

2.3 Conditions of Use

The 757B501 AxonEnergy Integral is to be used only for everyday activities.

2.4 User/Qualifications

Fitting a patient with the 757B501 AxonEnergy Integral may only be carried out by prosthetists who have been authorised by Ottobock after completion of a corresponding training course.

3 Safety

3.1 Explanation of Warning Symbols

WARNING

Warnings regarding possible risks of severe accident or injury.

CAUTION

Warnings regarding possible risks of accident or injury.

NOTICE

Warnings regarding possible technical damage.

3.2 General Safety Instructions

CAUTION

Non-observance of safety instructions

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.
 - ▶ Destruction of the battery.
- Follow the safety instructions in this accompanying document.

CAUTION

Penetration of liquid

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.

- Ensure that no liquid penetrates into the battery.
- Do not use the battery for bathing prostheses.

⚠ CAUTION

Disconnecting or establishing electrical connections when switched on

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.
- Power down the Axon-Bus prosthetic system by pressing the button in the charging receptacle before changing any Axon-Bus components (e.g. Axon-Bus gripping component).

⚠ CAUTION

Improper handling of the battery

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.
- Instruct the patient in the proper use of the battery.

3.3 Information to be Provided to the User

⚠ CAUTION

Penetration of dirt and moisture

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.
- Ensure that neither solid particles nor liquids can penetrate into the Axon-Bus prosthetic system.

⚠ CAUTION

Mechanical overload caused by mechanical influences or external loads

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.
- Do not subject the Axon-Bus components to mechanical vibrations or impacts.
- Check the Axon-Bus components for any visible damage before each use.

⚠ CAUTION

Wearing in extreme ambient temperatures

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.
- Avoid wearing the prosthesis in areas with temperatures outside of the permissible range. See the section "Technical Data".

⚠ CAUTION

Electromagnetic interference

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.
- Avoid areas where there are sources of strong electromagnetic interference (e.g. EAS systems in department stores).

⚠ CAUTION

Unauthorised manipulation of Axon-Bus components

- ▶ Injury due to faulty control or malfunction of the Axon-Bus prosthetic system.
- Manipulations to the Axon-Bus prosthetic system other than the tasks described in these instructions for use are not permitted.
- Do not damage the battery or separate the connection cables between the battery packs.

→ The Axon-Bus prosthetic system and any damaged Axon-Bus components may only be opened and/or repaired by certified Ottobock Myo-Service technicians.

⚠ CAUTION

Dirtying of the electrical contacts

- ▶ Risk of injuries due to faulty control or malfunction of the Axon-Bus prosthetic system caused by an insufficient supply of power to the Axon-Bus components.
- Ensure that the contacts of the charging receptacle are always clean and free of grease.
- Clean the electrical contacts of the battery charger and charging receptacle regularly using cotton swabs and a mild soap solution.
- Avoid damaging the contact surface coating with pointed or sharp objects.

NOTICE

Improper care of the housing

- ▶ Damage to the housing through the use of acetone, benzine or similar solvents.
- Only use a damp cloth and mild soap to clean the housing.

4 Delivery/Storage

4.1 Scope of Delivery

- 1 AxonEnergy Integral 757B501
- 2 Lamination Dummy – Battery
- 1 Lamination Dummy – Charging Receptacle
- 1 Drilling Template for Charging Receptacle
- 1 Component Holder Set
- 1 647G789 Instructions for Use

4.2 Storage

Observe the specifications in the chapter "Technical Data".

5 Product Description

5.1 Function

The 757B501 AxonEnergy Integral serves to provide power to the Axon-Bus prosthetic system.

6 Components

6.1 Rechargeable Battery

The rechargeable battery (Fig. 1, item 1) consists of Li-ion cells with integrated electronics. This protects the battery against short circuits, overvoltage, undervoltage and charging outside the allowable temperature range.

6.2 Charging Receptacle

The charging receptacle (Fig. 1, item 2) has the following functions:

- Contacts for charging the battery
- LED display to indicate the current battery charge level
- LED display for providing feedback on operating states

- Button to turn the Axon-Bus prosthetic system on and off, display the battery charge level and open the Axon-Bus gripping component in an emergency
- Beeper for providing feedback on operating states
- Button for activating the Bluetooth function.

6.3 Axon-Bus Cable

The Axon-Bus cable with the three-pin receptacle (Fig. 1, item 3) is used to exchange data and connects the respective Axon-Bus component to the battery.

6.4 Connection for 13E500 AxonMaster

Receptacle (Fig. 1, Item 4) for attaching the 13E500 AxonMaster. When using the battery in an elbow joint, e.g., 12K501=* AxonArm Ergo, the forearm cable leading to the flexion aid (AFB) is plugged into this receptacle.

7 Preparation for Use

7.1 Installation

The components are wired ready for installation and cannot be separated.

In combination with the component holder, the AxonEnergy Integral 757B501 can only be mounted in the forearm socket. The charging receptacle is installed in the outer socket through bonding with 617H46 Bonding Agent. Depending on how much space is available, the battery is secured in the outer socket or between the outer and inner socket. When little space is available, the battery is laminated into the outer socket.

NOTICE

Gluing the charging receptacle with an unsuitable adhesive

- ▶ Damage to the outer socket during removal of the charging receptacle due to the use of other adhesives (e.g. Cyamet quick-drying adhesive).
- 617H46 Bonding Agent must be used to glue the charging receptacle because the battery and charging receptacle must be removed from the outer socket in case of a defect.

7.2 Charging the Battery

The 757B501 AxonEnergy Integral must be charged prior to initial use. The 757B501 AxonEnergy Integral may **only** be charged with the 757L500 AxonCharge Integral. The charging process can only be performed with the 13E500 AxonMaster connected.

- 1) Connect the charging plug of the 757L500 AxonCharge Integral to the charging receptacle.
 - The beeper sounds briefly 2x.
 - The Axon-Bus prosthetic system is deactivated and the charging process starts automatically.
- 2) To activate the Axon-Bus prosthetic system, disconnect the charging plug and switch on the Axon-Bus prosthetic system.

INFORMATION: The prosthesis cannot be used during the charging process.

8 Handling

8.1 Switching On the Axon-Bus Prosthetic System

- ▶ Press the charging receptacle button and hold for one second (Fig. 2, arrow).
- The beeper emits two short signals.

8.2 Switching Off the Axon-Bus Prosthetic System

- ▶ Press the charging receptacle button and hold for more than one second (Fig. 2, arrow).
- The beeper emits one long signal.

8.3 Determining the Battery Charge Level

The battery charge level can be queried at any time.

1. With the Axon-Bus prosthetic system switched on, press the charging receptacle button (Fig. 2, arrow) and hold for less than one second.
2. The LED is illuminated. The colour of the LED indicates the current battery charge level.

Battery empty	Orange LED is illuminated
Battery 50% charged	Yellow LED is illuminated
Battery fully charged	Green LED is illuminated


8.4 Safety Shutoff

The purpose of the safety shutoff is to protect the battery; it is triggered in case of:

- Excessively high or low temperature
- Overvoltage or undervoltage
- Short circuit.

After a short circuit the battery charger must be connected and then disconnected in order for the electronics to be activated.

8.5 Open the Axon-Bus Gripping Component in an Emergency

This safety function allows the Axon-Bus gripping component to be opened regardless of the control signals present.

- 1) With the Axon-Bus prosthetic system switched on, press the charging receptacle button and hold for approximately three seconds until the Axon-Bus gripping component begins to open. A pulsating beep sounds as it opens.
- 2) Releasing the button immediately stops the process of opening the Axon-Bus gripping component and turns the entire Axon-Bus prosthetic system off.

8.6 Bluetooth Function

Pressing the button in the charging receptacle and holding for more than four seconds with the Axon-Bus prosthetic system turned off activates the Bluetooth function of the connected Axon-Master 13E500. The blue LED lights up. Once a connection to the PC has been successfully established, the LED flashes blue.

8.7 Beep Function

Action	Beep 1x long	Beep 2x short	Beep 3x short	Beep 1x short, 1x long	Beep pulsating
Switch on the Axon-Bus prosthetic system Press and hold the button for one to four seconds.					
Switch off the Axon-Bus prosthetic system Press the button and hold for more than one second.					
Charge the Axon-Bus prosthetic system Connect charging plug, charging process starts.					

Action	Beep 1x long	Beep 2x short	Beep 3x short	Beep 1x short, 1x long	Beep pulsating
Stop charging Disconnect charging plug. Charging ends and the Axon-Bus prosthetic system is switched off.					
Empty battery indicator					
Switch on the Bluetooth function (if the Axon-Bus prosthetic system was previously switched off) Press the button and hold for more than four seconds.					
Open the Axon-Bus gripping component in an emergency Press the button and hold for approximately three seconds, until the Axon-Bus gripping component opens.					

9 Fault Display and Fault Elimination

9.1 LED flashes red and beeper sounds for 3 seconds

Probable cause: Error in the Axon-Bus prosthetic system.

- Remedy:**
1. Switch the Axon-Bus prosthetic system off and back on.
 2. Check the plug connections in the socket.
 3. Connect the charger. If the charger displays the service symbol (yellow spanner) or the system error symbol (red spanner) send the Axon-Bus prosthetic system to Ottobock Myo-Service.

9.2 The LED flashes three times after pressing the button in the colour of the current battery capacity

Probable cause: The connection from the AxonEnergy Integral to the Axon-Master has been interrupted.

- Remedy:**
1. Check the connection (e.g. plug connections, flex cable, EasyPlug etc.).
 2. If no fault can be found, send the Axon-Bus prosthetic system to Ottobock Myo-Service.

9.3 LED alternately flashes yellow/red and beeper sounds for 3 seconds

Probable cause: The Axon-Bus prosthetic system is overheated.

- Remedy:**
1. Allow the Axon-Bus prosthetic system to cool for a few minutes.
 2. Do not use the Axon-Bus prosthetic system at high ambient temperatures.

9.4 The movements of the Axon-Bus gripping component are becoming increasingly slower

The patient is alerted to the diminishing battery charge level as the Axon-Bus gripping component increasingly slows down and the gripping force decreases. Finally, when there is very little battery capacity remaining, the Axon-Bus prosthetic system switches off automatically to protect the battery against harmful deep discharge.

10 Disposal

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

11 Legal Information

11.1 Liability

Otto Bock Healthcare Products GmbH, hereafter referred to as manufacturer, assumes liability only if the user complies with the processing, operating and maintenance instructions as well as the service intervals. The manufacturer explicitly states that this device may only be used in combination with components that were authorised by the manufacturer (see instructions for use and catalogues). The manufacturer does not assume liability for damage caused by component combinations which it did not authorise. The device may only be opened and repaired by authorised Ottobock technicians.

11.2 CE Conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

11.3 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

12 Appendices

12.1 Symbols Used

In some jurisdictions it is not permissible to dispose of these products with unsorted household waste. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the instructions of your national authority pertaining to return and collection.

Declaration of conformity according to the applicable European directives

Legal manufacturer

12.2 Technical Data

Capacity	approx. 1150 mAh
Output voltage	approx. 11.1 V DC
Charging voltage	approx. 12.3 V DC
Charging time	max. 3.5 h
Operating temperature	0 °C to +60 °C (+32 °F to +140 °F)
Storage temperature	-20 °C to +40 °C (-4 °F to +104 °F)
Charging temperature	0 °C to +40 °C (+32 °F to +104 °F)
Transport temperature	-20 °C to +40 °C (-4 °F to +104 °F)
Relative humidity	max. 80 % relative humidity, non-condensing
Weight	90 g
Dimensions (LxWxH)	55 x 35 x 23 mm

12.3 Axon

The term "Axon" stands for **A**daptive **eX**change **O**f **N**europlacement data. The Axon-Bus is an Ottobock innovation for the field of exoprosthetics: a data transmission system, derived from safety-related bus systems in the aviation and automobile industries. For the user this means enhanced safety and reliability because of a considerably reduced sensitivity to electromagnetic interference in comparison with conventional systems.

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