






## 4R204, 4R206

<b>DE</b> Gebrauchsanweisung .....	5
<b>EN</b> Instructions for use .....	12
<b>FR</b> Instructions d'utilisation .....	19
<b>IT</b> Istruzioni per l'uso .....	26
<b>ES</b> Instrucciones de uso .....	34
<b>PT</b> Manual de utilização .....	42
<b>NL</b> Gebruiksaanwijzing .....	49
<b>SV</b> Bruksanvisning .....	56
<b>DA</b> Brugsanvisning .....	64
<b>NO</b> Bruksanvisning .....	71
<b>FI</b> Käyttöohje .....	78
<b>HR</b> Upute za uporabu .....	85
<b>JA</b> 取扱説明書 .....	92

## INFORMATION

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### INFORMATION

Date of last update: 2020-09-11

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

### 1.1 Construction and Function

The 4R204 and 4R206 TF sport foot adapters are mounted on the 1E90 Sprinter prosthetic foot.

- **Selecting the prosthetic foot:** TF test sport foot adapter **4R206**
- **Definitive sport prosthesis:** TF definitive sport foot adapter **4R204**

To improve the readability of the texts, the TF test sport foot adapter is referred to as the test adapter below. The name definitive adapter is used for the TF definitive sport foot adapter.

The prosthetic foot and prosthetic sport knee joint are joined by the TF sport foot adapter with the help of a suitable socket adapter. The TF sport foot adapter permits anterior or posterior placement of the prosthetic foot.

The **4R206 test adapter** is used for the selection of a suitable prosthetic foot. The optimum position and height of the prosthetic foot in the prosthesis are determined with the help of the clamp. The clamp can be reopened as often as desired.

When the definitive sport prosthesis is fabricated, the prosthetic foot is adapted and the test adapter is replaced by the **4R204 definitive adapter**.

## 1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

### Limited combination options for Ottobock components

Designation	Reference Number
Socket Adapter	4R51, 4R54, 4R55, 4R73=A, 4R73=D, 4R77
Prosthetic Feet	1E90

## 2 Intended use

### 2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

### 2.2 Area of application

Approved for a body weight of **up to 100 kg**.

The product must be used in TF prostheses only.

### 2.3 Environmental conditions

Storage and transport
Temperature range $-20\text{ }^{\circ}\text{C}$ to $+60\text{ }^{\circ}\text{C}$ ( $-4\text{ }^{\circ}\text{F}$ to $+140\text{ }^{\circ}\text{F}$ ), relative humidity 20 % to 90 %, no mechanical vibrations or impacts
Allowable environmental conditions
<b>Temperature range:</b> $-10\text{ }^{\circ}\text{C}$ to $+45\text{ }^{\circ}\text{C}$
<b>Moisture:</b> relative humidity: 20% to 90%, non-condensing
<b>Chemicals/liquids:</b> fresh water as dripping water, occasional contact with salty air (e.g. near the ocean)
<b>Solids:</b> dust
Prohibited environmental conditions
<b>Chemicals/moisture:</b> salt water, perspiration, urine, acids, soapsuds, chlorine water
<b>Solids:</b> dust in high concentrations (e.g. construction site), sand, highly hygroscopic particles (e.g. talcum)

### 2.4 Service life

As part of a sport prosthesis, an exact service life cannot be determined for the product since strain varies widely depending on the sport and intensity of use.

## 3 Safety

### 3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

### 3.2 General safety instructions

 **CAUTION**

**Unallowable combination of prosthetic components**

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

 **CAUTION**

**Excessive strain on the product**

Risk of injury due to breakage of load-bearing components

- ▶ Use the product according to the specified area of application (see page 13).

 **CAUTION**

**Exceeding the service life**

Risk of injury due to change in or loss of functionality and damage to the product

- ▶ Ensure that the approved service life is not exceeded.

 **CAUTION**

**Improper assembly of the screw connections**

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.
- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

**NOTICE****Use under unallowable environmental conditions**

Damage to product due to unallowable environmental conditions

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop, etc.).

**Signs of changes in or loss of functionality during use**

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

**4 Scope of delivery**

4R206 TF test sport foot adapter					
Fig.	Item No.		Quantity	Designation	Reference Number
–	–	■	1	Instructions for Use	647G839
1	1		1	TF test sport foot adapter	–
1	2	■	1	Spacer Plate 3 mm	4G791
–	–	▲	4	Set Screw	4X298
–	–	▲	4	Countersunk Head Screw	501S41=M6x20

4R204 TF definitive sport foot adapter					
Fig.	Item No.		Quantity	Designation	Reference Number
–	–	■	1	Instructions for Use	647G839
2	1		1	TF definitive sport foot adapter	–
2	2	▲	4	Set Screw	506G3=M8x10
2	3	▲	1	Hexagon Head Screw	501K7=M10x30
2	4	▲	1	Flat Disc with Bevel	507U1=10
2	5	■	1	Spacer Plate 3 mm	4G791
2	6	■	1	Angle Plate	4X297
–	–	▲	4	Countersunk Head Screw	501S41=M6x20

## 5 Preparing the product for use

### CAUTION

#### **Incorrect alignment or assembly**

Risk of injury due to damaged prosthetic components

- ▶ Observe the alignment and assembly instructions.

### INFORMATION

- ▶ Not all of the materials may be available in your country. In this case, please contact the local branch of the manufacturer to obtain information on alternative materials.

### 5.1 Installing the Socket Adapter on the Test Adapter or Definitive Adapter

#### INFORMATION

#### **Positioning the socket adapter**

The socket adapter can be positioned on the test adapter or the definitive adapter in three different positions (offset **12 mm**).

#### **Required tools and materials:**

710D4 torque wrench, socket adapter, 4G791 spacer plate, 501S41=M6x20 countersunk head screws, degreasing cleaner (e.g. 636A58), 636K13 Loctite® 241

- 1) For definitive mounting, clean screw connection area using a degreasing cleaner.
- 2) Set the socket adapter with the intermediate pressure plate and the spacer plate onto the sport foot adapter in the desired position (see fig. 3).
- 3) Secure the screw connections with Loctite.
- 4) Mount the socket adapter on the sport foot adapter with the countersunk head screws (**torque: 12 Nm**).

### 5.2 Installing the Prosthetic Foot in the Test Adapter

#### > **Required tools and materials:**

710D4 Torque Wrench

- 1) Slide the prosthetic foot into the adapter from the distal direction and move it to the desired position (see fig. 4).
- 2) Tighten the 4 set screws to secure the prosthetic foot (see fig. 5) (torque: **15 Nm**).

### 5.3 Shifting the Prosthetic Foot

> **Required tools and materials:**

710D4 torque wrench

> No load on the prosthesis.

- 1) Loosen the 4 set screws (see fig. 4).
- 2) Shift the prosthetic foot in the test adapter.
- 3) Secure the prosthetic foot in the test adapter by tightening the 4 set screws (see fig. 5) (torque: **15 Nm**).

### 5.4 Marking and Processing the Prosthetic Foot

After determining the definitive prosthetic foot height, the position of the prosthetic foot is marked on the test adapter and the prosthetic foot is processed accordingly.

> **Required tools and materials:**

Saw with metal or diamond blade, carbide or diamond drill bit Ø 10 mm

- 1) Mark the cutting edge of the prosthetic foot with a pen (see fig. 6).
- 2) Mark the position of the bore hole with a pen (see fig. 7).
- 3) Remove the test adapter from the prosthetic foot.
- 4) Clamp the prosthetic foot in a vice with jaw pads.
- 5) Saw through the prosthetic foot at the marking (see fig. 8) and deburr the cut edge.
- 6) Drill the hole at the marked position and deburr it.

### 5.5 Installing the Prosthetic Foot in the Definitive Adapter

> **Required tools and materials:**

Degreasing cleaner (e.g. 636A58), 636K13 Loctite® 241, 710D4 torque wrench

- > The prosthetic foot is shortened to the definitive length and has a bore hole.
- > The socket adapter is mounted on the definitive adapter (see page 16)
- 1) For definitive mounting, clean screw connection area using a degreasing cleaner.
  - 2) Turn out the set screws of the definitive adapter so that they do not protrude into the insertion space.
  - 3) Set the angle plate onto the prosthetic foot (see fig. 9) and slide the prosthetic foot into the clamping space up to the stop (see fig. 10).
  - 4) Face the bevel on the washer toward the screw head and slide the washer onto the hexagon head screw.
  - 5) Secure the screw connections with Loctite.

- 6) Slide the hexagon head screw through the prosthetic foot into the definitive adapter and tighten it (see fig. 11) (torque: **30 Nm**).
- 7) Tighten the 2 posterior set screws (see fig. 12) (torque: **4 Nm**).
- 8) Tighten the 2 proximal set screws (torque: **4 Nm**).

## 6 Maintenance

### CAUTION

#### **Failure to follow the maintenance instructions**

Risk of injuries due to changes in or loss of functionality and damage to the product

- ▶ Observe the following maintenance instructions.
- ▶ Arrange regular maintenance intervals with the patient depending on the level of use.
- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.

## 7 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

## 8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

### 8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

## 9 Technical data

Reference Number	4R204	4R206
Weight [g]	440	580
Overall height [mm]	68	
System height [mm]	2	
Material	Aluminium	
Max. body weight [kg]	100	

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