



4R153*

DE Gebrauchsanweisung	4
EN Instructions for use	11
FR Instructions d'utilisation	17
IT Istruzioni per l'uso	24
ES Instrucciones de uso	31
PT Manual de utilização	38
NL Gebruiksaanwijzing	46
SV Bruksanvisning	53
DA Brugsanvisning	59
NO Bruksanvisning	66
FI Käyttöohje	72
TR Kullanma talimatı	79
RU Руководство по применению	85
JA 取扱説明書	93

INFORMATION

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INFORMATION

Date of last update: 2020-05-06

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

1.1 Construction and Function

The 4R153 adapter plate or 4R153=1 adapter plate with offset are installed between the prosthetic socket and the distal components of the modular prosthesis. They are mounted together with the socket adapter. The adapter plate features two connectors for generating a vacuum in the prosthetic socket and one connector for connecting to a vacuum pump. A PE sinter filter protects the vacuum pump from contamination. Additionally, the adapter plate is equipped with a check valve, which retains the vacuum generated in the prosthetic socket.

The adapter plate with offset can be used when a greater lateral clearance to the prosthetic knee joint or prosthetic foot is required.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

Approved for a body weight of **up to 150 kg**.

2.3 Environmental conditions

Storage and transport

Temperature range $-20\text{ }^{\circ}\text{C}$ to $+60\text{ }^{\circ}\text{C}$ ($-4\text{ }^{\circ}\text{F}$ to $+140\text{ }^{\circ}\text{F}$), relative humidity 20 % to 90 %, no mechanical vibrations or impacts

Allowable environmental conditions

Temperature range: $-10\text{ }^{\circ}\text{C}$ to $+45\text{ }^{\circ}\text{C}$

Moisture: relative humidity: 20% to 90%, non-condensing

Unacceptable environmental conditions

Chemicals/liquids: fresh water, salt water, perspiration, urine, acids, soapsuds, chlorine water

Solids: dust, sand, highly hygroscopic particles (e. g. talcum)

2.4 Service life

This prosthetic component has been tested by the manufacturer for 3 million load cycles according to ISO 10328. Depending on the patient's activity level, this corresponds to a service life of 3 to 5 years.

3 Safety

3.1 Explanation of warning symbols

 **CAUTION** Warning regarding possible risks of accident or injury.

 **NOTICE** Warning regarding possible technical damage.

3.2 General safety instructions

CAUTION

Excessive strain on the product

Risk of injury due to breakage of load-bearing components

- ▶ Use the prosthetic components according to the MOBIS classification.

CAUTION

Unallowable combination of prosthetic components

Risk of injury due to breakage or deformation of the product

- ▶ Only combine the product with prosthetic components that are approved for that purpose.
- ▶ Based on the instructions for use of the prosthetic components, verify that they may be combined with each other.

⚠ CAUTION

Use under unallowable environmental conditions

Risk of injury due to damage to the product

- ▶ Do not expose the product to unallowable environmental conditions.
- ▶ If the product has been exposed to unallowable environmental conditions, check it for damage.
- ▶ If damage is apparent or in case of doubt, do not continue using the product.
- ▶ Take suitable measures if required (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop etc.).

⚠ CAUTION

Exceeding the service life and reuse on another patient

Risk of injury due to loss of functionality as well as damage to the product

- ▶ Ensure that the approved service life is not exceeded.
- ▶ Only use the product on a single patient.

⚠ CAUTION

Mechanical damage to the product

Risk of injury due to change in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ If the product is damaged, check it for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, do not continue using the product (see "Signs of changes in or loss of functionality during use" in this section).
- ▶ Take any necessary measures (e.g. repair, replacement, inspection by the manufacturer's customer service, etc.).

Signs of changes in or loss of functionality during use

Among other factors, changes in functionality can be indicated by an altered gait pattern, a change in the positioning of the prosthetic components relative to each other and by the development of noises.

4 Scope of delivery

4R153* Adapter plate (see fig. 1)			
Item	Quantity	Designation	Reference number
1	1	Adapter plate (with or without offset)	–
2	1	Torx offset screwdriver	709S51=T8

4R153* Adapter plate (see fig. 1)			
Item	Quantity	Designation	Reference number
Single-component pack, consisting of:			4X305
3	1	PE solid rod (PE sinter filter)	–
4	1	TPU insert for 5R2 (thin)	–
		TPU insert for TEC plate (thick)	
5	1	Alternative connection for valve plate	–
6	1	O-ring	–
7	1	Locking screw	–
8	4	M6x38 mm screw	–
9	4	M6x35 mm screw	–

5 Preparing the product for use

CAUTION

Incorrect alignment or assembly

Risk of injury due to damaged prosthetic components

- Observe the alignment and assembly instructions.

5.1 Installing the Adapter Plate

CAUTION

Using a wrong prosthetic socket

Contusions and pressure sores due to wrong socket fit

- Use the vacuum pump only in combination with a prosthetic socket made with Harmony technology.

5.1.1 Preparing the Adapter Plate

A vacuum can be generated in the socket either via the proximal or the lateral connector of the adapter plate.

INFORMATION

Using the proximal connector of the adapter plate is recommended only when using a 5R2 lamination disc or the 6A94=3 TEC plate (with TPU insert respectively). The lateral connector for the intake hose can be recognised by the internal thread for the hose connector.

Proximal connection

INFORMATION

To ensure an optimum sealed connection to the socket, the following approach is recommended:

- ▶ Drill the distal socket bore with a **6 mm** bit.
- ▶ Apply PU adhesive inside the bore and insert the included hose piece (see fig. 2).
- ▶ Remove excess adhesive and allow the adhesive to set.
- ▶ Trim the hose piece flush on both ends.

> Required tools and materials:

710D4 torque wrench, 501Z61=M12x1.5 locking screw, 4X282 TPU insert (thin) for 5R2 or 4X279 TPU insert (thick) for TEC plate, 4X304 PE solid rod (PE sinter filter)

- 1) Screw the locking screw into the lateral connector of the adapter plate and tighten (mounting torque: **3 Nm**) (see fig. 3).
- 2) Press the PE sinter filter into the TPU insert to the stop (see fig. 4).
- 3) Slide the TPU insert with the PE sinter filter into the connector on the proximal side of the adapter plate (see fig. 5).

Lateral connection

> Required tools and materials:

710D4 Torque Wrench, 2 mm Allen Wrench, 501Z61=M12x1.5 Locking Screw, 4X304 PE Solid Rod (PE Sinter Filter), 627F13=7X1 O-Ring, 4X296 Hose Connector

- 1) Screw the locking screw into the proximal connector of the adapter plate and tighten (mounting torque: **3 Nm**) (see fig. 6).
- 2) Press the PE sinter filter into the lateral connector of the adapter plate up to the stop (see fig. 6).
- 3) Insert the O-Ring into the groove of the hose connector.
- 4) Use the Allen wrench to screw the hose connector into the lateral connector and tighten (installation torque: **1 Nm**) (see fig. 6).

5.1.2 Installing the Adapter Plate

⚠ CAUTION

Improper assembly of the screw connections

Risk of injury due to breakage or loosening of the screw connections

- ▶ Clean the threads before every installation.
- ▶ Apply the specified torque values.

- ▶ Follow the instructions regarding the length of the screws and about how to secure the screws.

INFORMATION

Note the position of the Harmony E2 vacuum pump in the prosthesis. The vacuum pump must not bump against the prosthesis or project beyond the flexed prosthetic knee joint. Use the 4R153=1 adapter plate with offset if you require greater lateral clearance (see fig. 10, see fig. 11).

Selection table for countersunk head screws

Chosen socket adapter	Screws to be used
4R22, 4R23, 4R54, 4R55, 4R73=A, 4R73=D, 4R74, 4R95	501S128=M6X35
4R37, 4R51, 4R77	501S128=M6X38

> **Required tools and materials:**

710D4 Torque Wrench, 636K13 Loctite® 241

- 1) Position the adapter plate on the lamination disc with the desired alignment (see fig. 7).
- 2) Set the chosen socket adapter onto the adapter plate (see fig. 8).
- 3) **INFORMATION: Secure the screws with Loctite®.**

Insert the screws through the socket adapter and the adapter plate into the lamination disc and, in diagonal sequence, tighten them using a torque wrench (mounting torque: **12 Nm**) (see fig. 9).

5.1.3 Optional: Installing the Intake Hose

An intake hose only needs to be installed if the lateral connector of the adapter plate will be used.

- 1) Connect the intake hose to the lateral connector of the adapter plate.
- 2) Run the intake hose to the intended connector of the prosthetic socket and connect it.

6 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Conduct annual safety inspections.

7 Disposal

This product may not be disposed of with regular domestic waste in all jurisdictions. Disposal that is not in accordance with the regulations of the country where the product is used may have a detrimental impact on health and

the environment. Please observe the information provided by the responsible authorities in the country of use regarding return, collection and disposal procedures.

8 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

8.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

8.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

9 Technical data

Reference number	4R153	4R153=1
Weight [g]	125	160
System height [mm]	22	22
Material	Aluminium	
Max. body weight [kg]	150	

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