ottobock.

50S16-7 Malleo Sprint light

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INFORMATION

Last update: 2013-07-22

- ▶ Please read this document carefully.
- ▶ Follow the safety instructions.

These instructions for use contain important information on the fitting and application of the 50S16-7 Malleo Sprint Light ankle orthosis.

2 Application

2.1 Intended Use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

2.2 Indications

- Functional treatment of ligament instability
- Post-operative treatment of ligament ruptures
- Ankle sprains
- Functional, conservative therapy of lateral ligament ruptures

Indications must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

Not known.

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb; pronounced varicose veins, especially with return flow impairment, lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the support will be applied; sensory and circulatory disorders in the foot area, e.g. in case of diabetic neuropathy.

2.4 Effects

The orthosis stabilises the ankle joint using straps. Eversion and inversion of the ankle joint are restricted without negatively affecting plantar and dorsal flexion.

3 Safety

3.1 Explanation of Warning Symbols

CAUTION

Warning regarding possible risks of accident or injury.

NOTICE

Warnings regarding possible technical damage.

3.2 General Safety Instructions

CAUTION

Using the product without following the instructions for use

Deteriorating health condition and damage to the product because of failure to follow the safety information

- ▶ Observe the safety information in these instructions for use.
- ▶ Make the patient aware of all safety instructions marked with "**Inform the patient**".

CAUTION

Reuse on other patients and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ Only use the product for a single patient.
- ▶ Clean the product regularly.
- ▶ **Inform the patient.**

⚠ CAUTION

Contact with heat, embers or fire

Injuries (e.g. burns) caused by melting of the material

- ▶ Keep the product away from open flame, embers and other sources of heat.
- ▶ **Inform the patient.**

NOTICE

Contact with salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to salves, lotions or other products that contain oils or acids.
- ▶ **Inform the patient.**

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Size Selection

- ▶ Select the orthosis size based on the circumference of the ankle (see size chart).

4.2 Fitting and Application

CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the orthosis is applied properly and fits correctly.
- ▶ **Inform the patient.**

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Instruct the patient to inspect the product for proper functionality, wear and damage prior to each use.
- ▶ Advise the patient that the product must not be used if the product or part of the product exhibits signs of wear (e.g. cracks, deformation, poor fit) or damage.

NOTICE

Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Only use this product carefully and for its intended purpose.
- ▶ Do not make any improper changes to the product.
- ▶ **Inform the patient.**

The orthosis is delivered ready for trial fitting. No adaptation is required.

- > The strap on the orthosis is open and the fasteners are loosened.
- 1) Pull the orthosis over your foot (see Fig. 1) and tie it like a normal shoe (see Fig. 2).
 - 2) Pass the outer end of the strap over the instep and pull it under the foot from the inner (medial) side.
 - 3) Fasten the outer end of the strap over the ankle on the outer (lateral) side (see Fig. 3).
 - 4) Pass the inner end of the strap over the instep and pull it under the foot from the outer (lateral) side.
 - 5) Fasten the inner end of the strap over the ankle on the inner (medial) side (see Fig. 4).
 - 6) If the physician has allowed weight-bearing on the foot:
Loose the ends of the strap and re-tighten (see Fig. 5).

4.3 Cleaning

- 1) Close all hook-and-loop closures.
- 2) Hand wash the orthosis in warm water at 40°C with mild laundry detergent.
- 3) Rinse well.
- 4) Allow to air dry. Do not expose to direct heat sources (e.g. sunlight, stove or radiator).

5 Disposal

Dispose of the product in accordance with national regulations.

6 Legal Information

6.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

6.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

Otto Bock HealthCare GmbH
Max-Näder-Straße 15 · 37115 Duderstadt · Germany
T +49 (0) 5527 848-0 · F +49 (0) 5527 72330
healthcare@ottobock.de · www.ottobock.com

Ottobock has a certified Quality Management System in accordance with ISO 13485.