


## 4R69, 2R50, 2R49

<b>DE</b>	Rohradapter Ø 30 mm, Schraubadapter Ø 30 mm .....	4
<b>EN</b>	Tube adapter Ø 30 mm, tube clamp adapter Ø 30 mm .....	6
<b>FR</b>	Adaptateur tubulaire Ø 30 mm, adaptateur à vis Ø 30 mm .....	9
<b>IT</b>	Tubo modulare Ø 30 mm, giunto modulare Ø 30 mm .....	11
<b>ES</b>	Adaptador de tubo Ø 30 mm, adaptador con rosca Ø 30 mm .....	14
<b>PT</b>	Adaptador com tubo Ø 30 mm, adaptador roscado Ø 30 mm .....	16
<b>NL</b>	Buisadapter Ø 30 mm, klemadapter Ø 30 mm .....	19
<b>SV</b>	Röradapter Ø 30 mm, Skruvadapter Ø 30 mm .....	21
<b>DA</b>	Røradapter Ø 30 mm, Skrueadapter Ø 30 mm .....	24
<b>NO</b>	Røradapter Ø 30 mm, Rørklemmeadapter Ø 30 mm .....	26
<b>PL</b>	Aluminiowe adaptery rurowe .....	28
<b>HU</b>	Csőadapter Ø 30 mm, csavaradapter Ø 30 mm .....	31
<b>CS</b>	Trubkový adaptér Ø 30 mm, šroubovací adaptér Ø 30 mm .....	33
<b>RO</b>	Adaptor tubular Ø 30 mm, Dispozitiv de adaptare cu șuruburi Ø 30 mm .....	35
<b>HR</b>	Prilagodnik cijevi Ø 30 mm, prilagodnik vijka Ø 30 mm .....	38
<b>TR</b>	Boru adaptörü Ø 30 mm, vidalama adaptörü Ø 30 mm .....	41
<b>EL</b>	Προσαρμογέας σωλήνα Ø 30 mm, βιδωτός προσαρμογέας Ø 30 mm .....	43
<b>RU</b>	Несущие модули Ø 30 мм, втулки Ø 30 мм .....	46
<b>JA</b>	チューブアダプター直径30mm, チューブクランプアダプター直径30mm .....	48
<b>ZH</b>	管连接器Ø30 毫米, 带管连接器Ø 30 毫米 .....	51

■	●	▲
Ⓓ Einzelteile	Einzelteile-Pack	Mindestmenge
ⒼⒷ Single components	Single-Component Pack	Minimum order quantity
Ⓕ Pièces à l'unité	Kit de pièces de rechange	Quantité minimum
Ⓘ Singoli componenti	Imballo dei singoli componenti	Quantità minima
Ⓔ Componentes	Kit componentes	Cantidad mínima
Ⓟ Componentes	Conjunto de peças de reposo	Quantidade mínima
ⓃⓁ Onderdelen	Onderdelenpakket	Minimikvantitet
Ⓢ Separata delar	Service-Set	Minimum aantal
ⒹⓀ Komponenter	Komponentpakke	Minste mængde
Ⓝ Enkle komponenter	Enkeltkomponentpakke	Minste ordremængde
ⓅⓁ Podzespoły	Zestaw naprawczy	Ilość minimalna
Ⓗ Alkatrészek	Alkatrész-csomag	Minimum mennyiség
ⒸⓏ Jednotlivé díly	Sada jednotlivých dílů	Minimální množství pro objednání
Ⓡⓐ Componente individuale	Set componente individuale	Cantitate minimă
ⒽⓇ Pojedinačni dijelovi	Pakiranje pojedinačnih dijelova	Minimalna količina
ⓉⓇ Münferit parçalar	Münferit parça ambalajı	Minimum miktar
ⒸⓇ Μεμονωμένα εξαρτήματα	Σετ μεμονωμένων εξαρτημάτων	Ελάχιστη ποσότητα
ⓇⓁⓈ Отдельные детали	Комплект деталей	Минимальное количество
Ⓝ 単一成分パック	単一成分パック	最低発注量
ⓃⓈ 单个部件	维修组件	最低起订量

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English

Last update: 2015-09-30

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Instruct the user in the proper and safe use of the product.
- Please keep this document in a safe place.

**Note:** The scope of delivery is shown on the front page.

## 1 Single components (fig. 1)

### 1.1 Single component

–

### 1.2 Single-component pack

–

### 1.3 Minimum order quantity

- |     |             |           |
|-----|-------------|-----------|
| (1) | 506G3=M8×14 | Set Screw |
| (2) | 501Z2=M6×25 | Cap Screw |

## 2 Description

### 2.1 Intended use

The structural parts shown on the front page are to be **solely** used for the prosthetic fitting of the lower limb.

### 2.2 Field of application

Field of application according to the **Ottobock MOBIS Mobility System**:

2R49, 2R50, 4R69

Approved for a body weight of up to 125 kg

**Attention!**

Please avoid exposing prosthetic components to corrosive elements such as fresh water, salt water, and acids. Using this medical product under such conditions will render all claims against Otto Bock HealthCare null and void.

**Please also inform your patients!**

### 2.3 Function

The tube adapters are used in combination with the tube clamp adapters as connection elements for modular lower extremity prostheses. These adapter combinations allow for controlled angle and length adaptation in the sagittal and frontal plane as well as adjustment of inward and outward rotation.

**Attention!**

To achieve maximum stability for the transtibial prostheses, we recommend using a  $\varnothing$  34 mm tube adapter (e. g. 2R57/2R76) together with a  $\varnothing$  34 mm tube clamp adapter (e. g. 4R82/4R91). This applies in particular if the prosthesis is to withstand increased pressure due to higher activity levels, longer foot or transtibial lever arms or similar factors.

## 3 Technical data

	2R50	2R49	4R69
Weight	155 g	240 g	75 g
System height	–	–	33 mm
System height min.	97 mm	97 mm	–
System height max.	232 mm	432 mm	–
Material	Aluminum	Aluminum	Aluminum
Diameter	30 mm	30 mm	30 mm
Patient weight max.	125 kg (275 lbs)	125 kg (275 lbs)	125 kg (275 lbs)

## 4 Handling

### 4.1 Assembly

The tube adapter must be shortened according to the measurements of the amputee.

**Attention!**

The tube adapter must always be installed distally.

**Attention!**

Do not clamp the tube into a vice ! It may re-shape or damage the tube.

**The section surface must be perpendicular to the tube longitudinal axis (Fig. 2).** We recommend using 719R2 Tube Cutter or equivalent.

**Carefully debur the tube cutting edge inside and outside.** For this purpose use the 718R1 Interior and Exterior Pipe Reamer or equivalent (Fig. 3).

Loosen the cap screw (Fig. 4, item 2) of the tube clamp adapter by two turns and insert the tube 5 mm deep into the tube clamp adapter while slightly rotating the tube. Next, place the tube clamp adapter with the tube inserted vertically on the worktop and fully insert the tube adapter into the tube clamp adapter until the tube is in contact with the base rim inside the clamp.

**Attention!**

For reasons of structural strength, the residual pieces of cut tube adapters must not be used as components in definitive prostheses.

## 4.2 Adjustments and final assembly

### Attention!

Prior to final assembly (i.e. prior to inserting the tube adapter into the tube clamp adapter), the inside of the tube clamp adapter and the outside of the tube adapter must be cleaned with 634A3 (Acetone)!

For final assembly, tighten the cap screw with 710D1 Torque Wrench.  
Torque: **13 Nm** (Fig. 4, item 2).

### Attention!

The clamping slot of the tube clamp adapter must always be positioned anteriorly.

Static corrections are possible at any time during alignment and trial walking as well as after completion of the prosthesis using the four 506G3 Set Screws. To exchange a module or during disassembly the previous position may be retained by loosening only the two most deeply screwed in, adjacent set screws.

Replace set screws that appear too long (to prevent foam cover damage) or too short (stability) with the appropriate screw length.

Set screws are available in the following sizes.

Article number	Length	Article number	Length
506G3=M8x10	10 mm	506G3=M8x14	14 mm
506G3=M8x12-V	12 mm	506G3=M8x16	16 mm

Use 710D4 Torque Wrench to tighten the set screws.

**For trial fitting:** Tighten the 506G3 Set Screws clockwise to **10 Nm**.

**To finish the prosthesis:** Secure the 506G3 Set Screws with 636K13 Loctite and pre-tighten to **10 Nm**. Then complete assembly by tightening to **15 Nm** (Fig. 5, item 1).

## 4.3 Maintenance instructions

### Note:

As a basic principle, all Ottobock modular adapters are subjected to tests involving three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

We recommend carrying out regular safety checks once a year.

## 5 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

### 5.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregard of this document, particularly due to improper use or unauthorised modification of the product.

### 5.2 CE conformity

This product meets the requirements of the European Directive 93/42/EEC for medical devices. This product has been classified as a class I device according to the classification criteria outlined in Annex IX of the directive. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to Annex VII of the directive.

### 5.3 Warranty

The manufacturer warrants this device from the date of purchase. The warranty covers defects that can be proven to be a direct result of flaws in the material, production or construction and that are reported to the manufacturer within the warranty period.

Further information on the warranty terms and conditions can be obtained from the competent manufacturer distribution company.

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Otto Bock HealthCare GmbH