## 743Y50 Anatomical SIT-Cast

<b>DE</b> Gebrauchsanleitung (Benutzer / Techniker) .....	3
<b>EN</b> Instructions for Use (User / Technician) .....	19
<b>FR</b> Mode d'emploi (Utilisateur / Technicien) .....	35
<b>IT</b> Istruzioni d'uso (utilizzatore / tecnici) .....	51
<b>ES</b> Instrucciones de uso (usuario / técnico) .....	67

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DE

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# Table of contents

EN

<b>1</b>	<b>Basic Information.....</b>	<b>20</b>
1.1	Publication Date .....	20
1.2	Product Data .....	20
1.3	Foreword .....	20
<b>2</b>	<b>Product Description.....</b>	<b>21</b>
2.1	Intended Use .....	21
2.2	Scope of Delivery .....	21
2.3	Interchangeable Components.....	21
2.4	Additional Components Required.....	22
2.5	Control Elements .....	22
<b>3</b>	<b>Safety.....</b>	<b>23</b>
3.1	Explanation of Warning Symbols .....	23
3.2	General Safety Instructions .....	23
<b>4</b>	<b>Technical Data.....</b>	<b>24</b>
4.1	Dimensions .....	24
4.2	Mass .....	24
<b>5</b>	<b>Assembly.....</b>	<b>24</b>
5.1	Mounting on the SIT-Cast Casting Device.....	24
5.2	Frontal Support Replacement and Installation.....	24
<b>6</b>	<b>Use.....</b>	<b>25</b>
6.1	Plaster Casting Technique .....	25
6.1.1	General Information.....	25
6.1.2	Advantages of the Anatomical Shape.....	25
6.1.3	Adjustment Possibilities for Taking Casts of Anatomical Structures .....	25
6.1.4	Taking Patient Measurements .....	26
6.1.5	Transferring the Measurements to the Anatomical SIT-Cast.....	27
6.1.6	Preparing to Take a Plaster Cast.....	29
6.1.7	Taking the Plaster Cast.....	29
6.1.7.1	Applying the Plaster Longuettes and Bandages, Setting .....	29
6.1.7.2	Modelling the Negative .....	31
6.1.7.3	Checking the Plaster Negative.....	32
6.1.7.4	Trial Fitting of the Plaster Socket.....	34
<b>7</b>	<b>Legal Information.....</b>	<b>34</b>
7.1	Liability .....	34

# 1 Basic Information

## 1.1 Publication Date

### INFORMATION

These instructions have been approved since the following date: 2012-06-12

## 1.2 Product Data

These instructions apply to the following products:

Designation	Article number	Comments
Anatomical SIT-Cast	743Y50	-

## 1.3 Foreword

These instructions contain important and useful information for all persons who use the product.

### ⚠ WARNING

#### **Risk to the user and the product due to failure to observe the instructions**

Risk of injury and / or death (user) – damage (product)

- ▶ Only use the product according to the information in these instructions.
- ▶ Always keep these instructions readily accessible in the vicinity of the product.

### Layout of the instructions

The instructions are divided into the following main sections:

- Basic Information
  - This section contains information about the instructions.
- Product Description
  - This section contains basic information about the product.
- Safety
  - This section contains basic safety information to protect persons and the product. The subsequent sections contain specific safety information on the topics in the respective section.
- Technical Data
  - The most important technical data are presented in this section in table form.
- Main sections on all important life cycle phases of the product (e.g. assembly, operation / use / handling etc.)
  - The main sections provide detailed information on the use of the product, depending on the product and the type of instructions.

## 2 Product Description

### 2.1 Intended Use

The Anatomical SIT-Cast in combination with the casting apparatus and the pneumatic SIT-Cast casting device or the SIT-Cast casting device is to be used exclusively for the plaster casting technique on patients with transfemoral amputations. It serves to create a plaster negative of a thigh socket with anatomical casting to achieve an exact fit around the ischium (tuber os ischii) under load.

The Anatomical SIT-Cast is a practical alternative to the manual plaster casting technique and can also be used on site, e.g. in the clinic, at the patient's home, etc.

### 2.2 Scope of Delivery

The Anatomical SIT-Cast is supplied as a base device with 2 adjustable frontal supports (left / right), with adjustable casting forms for the perineum area, for installation on the 743G10 Pneumatic Sit-Cast Casting Device or the 743G1 SIT-Cast Casting Device.


### 2.3 Interchangeable Components

#### Frontal supports


The 743Y50 Anatomical SIT-Cast is supplied with frontal supports for the right and left. Select the corresponding frontal support according to the patient's amputated side.

## 2.4 Additional Components Required

	Item no.	Designation	Article number
	1	Anatomical SIT-Cast	743Y50
	2	Casting Apparatus	743A11
	3	Pneumatic SIT-Cast Casting Device	743G10
	-	SIT-Cast Casting Device (same as 743G10 without pneumatics)	743G1

The required 743A11 Casting Apparatus and the SIT-Cast Casting Device, which is available in the versions 743G10 and 743G1, are not included in the scope of delivery. These must be ordered separately.

## 2.5 Control Elements

The adjustment wheels (1 + 5) are used to adjust the position and angle of the tuberosity supports (4). Consult the illustration above for the effect that turning the adjustment wheels has on the tuberosity supports.

The frontal support (3) can be adjusted after loosening the knurled screw (2). Afterwards re-tighten the knurled screw.

When adjusting the position of the tuberosity supports (4), first use the adjustment wheel (5) to set the approximate position. Then use the other adjustment wheel (1) to fine-tune the angle. Finally, adjust the position of the frontal support (3).

## 3 Safety

### 3.1 Explanation of Warning Symbols

#### **DANGER**

**This warning identifies an immediately impending hazard.**

If the hazard is not avoided, death or serious injury (irreversible) will result.

#### **WARNING**

**This warning identifies a possibly hazardous situation.**

If the hazard is not avoided, death or serious injury (irreversible) may result.

#### **CAUTION**

**This warning identifies a possibly hazardous situation.**

If the hazard is not avoided, slight or minor injuries (reversible) may result.

#### **NOTICE**

**This warning identifies a possibly detrimental situation.**

If the hazard is not avoided, the product or something in the vicinity may be damaged.

### 3.2 General Safety Instructions

#### **WARNING**

**Risk of falling due to the patient's lack of ability to stand**

Injuries caused by falling

- ▶ Do not use the product if the patient does not have a sufficient ability to stand.
- ▶ Use a standing aid (e.g. standing support, forearm crutches) during the entire process.
- ▶ Having two prosthetists in attendance while taking the plaster cast is recommended.

#### **WARNING**

**Risk of falling due to unsecured screw connections**

Injuries caused by falling

- ▶ Check the tightness of all screw connections.
- ▶ Tighten all screw connections that are loose.
- ▶ Secure the locking screw after adjusting the height of the 743A11 Casting Apparatus to the tuberosity-floor measurement of the patient.
- ▶ Be sure to tighten the clamping screw of the support after completing the adjustment.

#### **CAUTION**

**Risk of health impairment due to improper hygiene**

Skin irritation, formation of eczema or infections

- ▶ Clean and disinfect the product after every use.

## 4 Technical Data

### 4.1 Dimensions

<b>Length</b>	444 mm / 17 1/2"
<b>Width</b>	140 mm / 5 1/2"
<b>Height</b>	181 mm / 7 1/8"

### 4.2 Mass

<b>Total weight</b>	2.5 kg / 5.5 lbs
---------------------	------------------


## 5 Assembly

### 5.1 Mounting on the SIT-Cast Casting Device

The 743Y50 Anatomical SIT-Cast is mounted on the 743G10 Pneumatic SIT-Cast Casting Device or the 743G1 SIT-Cast Casting Device. Firmly tighten the knurled screws.

### 5.2 Frontal Support Replacement and Installation

Assignment of the frontal supports to the amputated side:

1. Frontal support with amputation on the right
2. Frontal support with amputation on the left

To change the frontal support, remove the knurled screw and replace the frontal support. Firmly tighten the knurled screws after reassembly.

## 6 Use

### 6.1 Plaster Casting Technique

#### 6.1.1 General Information

The custom-made residual limb socket is of the utmost importance for the quality and wearer comfort of a prosthesis. For transfemoral sockets in particular, the lack of residual limb load-bearing capacity due to anatomical reasons requires the application of special design criteria.

The connecting element between the body of the amputee and the prosthesis has to fulfil the following functions:

- Correct containment of the residual limb volume
- Secure adhesion of the prosthesis on the residual limb
- Favourable load and force transmission
- Secure guidance of the prosthesis

Arterial, venous and lymphatic circulation must not be impaired; non-physiological pressure also has to be avoided.

The various parts of the residual limb socket have to conform to anatomical features; the circumference is reduced based on functional criteria in order to compress the residual limb.

Using the Anatomical SIT-Cast Casting Apparatus allows a systematic approach for taking plaster casts of the transfemoral residual limb to be implemented without disregarding individual criteria. In combination with the 743A11 Otto Bock Casting Apparatus and the 743G10 Pneumatic SIT-Cast Casting Device or the 743G1 SIT-Cast Casting Device, plaster casts are taken under load.

#### 6.1.2 Advantages of the Anatomical Shape

The socket shape has the following advantages:

- Sitting is more comfortable thanks to the gluteus cutout, since the patient is not sitting on part of the socket (no interference by the plastic material). This offers freedom of movement, an improved and more functional gait pattern and, very important for the patient, an improved cosmetic appearance.
- The lower brim significantly enhances freedom of movement. Full hip flexion is possible and, depending on the length of the residual limb and hip flexibility, the patient can achieve virtually normal external rotation up to 90° (sitting cross-legged).
- The appearance of the hip and buttocks while walking is more physiological. There is virtually no apparent difference compared to the sound side.

#### 6.1.3 Adjustment Possibilities for Taking Casts of Anatomical Structures

The precisely adjustable casting elements serve to precisely capture anatomical structures already while taking the plaster cast. This makes it easier to produce a custom plaster negative for the patient.

In order to follow the individual anatomy, the elements of the ramus embracing are continuously adjustable in two dimensions.

1. The ramus angle is precisely captured by the elements of the ramus embracing. These should run parallel to the ramus branch and are positioned more anterior compared to classic CAT-CAM sockets. At this position, the angle is approximately 30° for the male pelvis and 40° to 45° for the female pelvis. These values are guidelines only and should be verified on a case by case basis. During advance adjustment and while taking the plaster cast, ensure that the elements of the ramus embracing are flush with the ischial tuberosity on the posterior side.
2. The distance of the elements of the ramus embracing in the medio-lateral direction can be varied, and depends on the a-p position of the elements as well as the ramus angle.

The frontal support has a recess in the medial direction to relieve pressure on the area over the adductor tendons.

It is also continuously adjustable in the anterior-posterior direction in order to set the precise a-p measurement (between the adductor tendons and the posterior part of the ischial tuberosity). During advance adjustment ensure that the elements of the ramus embracing are flush with the ischial tuberosity on the posterior side.

### 6.1.4 Taking Patient Measurements

To protect against wet plaster and for soft tissue compression, apply 2 layers of body protection stockinette. Measure the width of the pelvis below the trochanter from the dorsal direction.

Possible positions for marks are the proximal socket brim, the course of the femur and the trochanter major as well as the circumference measurement at intervals of 5 cm.

Measure the tuberosity-sub-trochanter distance.

The lateral calliper arm faces in the walking direction while the medial calliper arm bends in the joint.

Ensure that the anterior part of the arm is positioned flat on the tuberosity.

Measure the distance from the medial ischium branch to below the trochanter with the callipers from the dorsal direction ("bony" medio-lateral residual limb dimensions).

Read the angle.

Measure the diagonal distance.

The medial calliper arm has to be flat on the tuberosity / ramus.

Measure the a-p distance.

1. Localise the adductor tendon.
2. The anterior calliper arm is applied to the adductor, the posterior arm to the tip of the tuberosity.

**Attention:** Do not exert compression force!

The diagonal a-p measurement is used for additional verification.


Measure the residual limb length from the ischial tuberosity to the end of the residual limb.

Take the measurement along the residual limb axis.

If the residual limb is in a highly abducted position, tilt the callipers accordingly. The residual limb musculature must be relaxed when doing so, in order to keep the callipers from slipping off.


### 6.1.5 Transferring the Measurements to the Anatomical SIT-Cast

Use the front and rear adjustment wheels to set the tuberosity width. First adjust the posterior and then the anterior width. This results in the angle.

Fine-tune the angle of the tuberosity supports.

Adjust the a-p width.

Conduct a sitting test as described below.  
Palpate the tip of the tuberosity. It should be at the level of the posterior support tip.  
Verify support contact.  
Readjust the width and angle if required.

The patient must be sitting comfortably and has to confirm that the supports are snug against the tuberosity / ramus.

Readjust the a-p socket width at the front if required.

## 6.1.6 Preparing to Take a Plaster Cast

Prepare the plaster longuettes according to the SIT-Cast plaster casting technique:

1x 20 m, 4 layers: lateral section

2x 15 cm, 2 layers: residual limb end

1x 10 cm, 3 layers: medial socket brim

4 – 5 elastic plaster bandages, 15 cm

Select a 99B4 Rubber Cuff of suitable size and pull it over the brim.

The rubber cuff compresses the residual limb and protects the device against wet plaster.

Apply plaster insulation cream to the rubber cuff.

## 6.1.7 Taking the Plaster Cast

### 6.1.7.1 Applying the Plaster Longuettes and Bandages, Setting

#### Summary:

The transfemoral residual limb is covered with Cellona longuettes and elastic plaster bandages. Then a stockinette is applied for pulling into the SIT-Cast casting device.

Possible positions for marks are the proximal socket brim, the course of the femur and the trochanter major as well as the circumference measurement at intervals of 5 cm.

Apply and mould the 4-layer plaster longuette (20 cm) from the residual limb end to the iliac crest on the lateral side.

Cover the residual limb end with another plaster longuette from medial to lateral and anterior to posterior.

To reinforce the proximal brim, apply a 3-layer plaster longuette (10 cm) in the pubic bone area. It runs from medial along the groin to the iliac crest.

Wrap three or four elastic plaster bandages around the residual limb from the proximal end. Start dorsally at the level of the iliac crest and work quickly. Be sure to keep the tension consistent.

Cut 12 cm wide perlon stockinette to double the length of the residual limb.

Pull the perlon stockinette over the residual limb and secure it with another elastic plaster bandage.

Have the patient get into the device.

Briefly put the maximum weight on the residual limb side and tighten the perlon stockinette.

Apply the lateral support at approximately 0.5 to 1.2 bar.

Slightly model the gluteus pocket.  
Allow the plaster negative to harden under load.

Cut the upper body protection stockinette and have the patient exit the device.

#### **INFORMATION**

The same plaster casting technique can also be used for fittings with a liner.

#### **6.1.7.2 Modelling the Negative**


Fill the negative with plaster on the outside and reinforce it with a plaster bandage.

Fill the frontal and dorsal bulges in the plaster negative with mixed plaster.  
Carry out additional modelling above the trochanter major and the dorsal-lateral areas.  
Reduce the circumference by 3 % – 5 %.

**INFORMATION**

In case of a pronounced flexion position of the residual limb, the lower corner of the frontal area can develop a large bulge. Use a knife to cut away this corner so that a smooth surface results.

Deepen the formed channel in the adductor area so that the patient does not perceive pressure in this area.  
Smooth the surface and finely sand it with wet sanding paper until the model is ready for the trial fitting.

**6.1.7.3 Checking the Plaster Negative**

Check the plaster negative.

Trim the socket brim.


### Checking the measurements

Check the m-l measurement.

Check the a-p measurement.

Check the diagonal measurement and circumference measurements.

#### 6.1.7.4 Trial Fitting of the Plaster Socket

For the trial fitting, drill a pull-in hole in the middle of the socket bottom.

Set the plaster socket on the gimbal ring of the plaster casting apparatus and adjust the height.

Pull the residual limb into the plaster socket with a bandage.

Ask the patient to distribute his or her weight evenly and verify that the pelvis is straight.

Once a satisfactory fit is achieved, control bores can be drilled into the plaster negative in order to check the residual limb volume.

Check the following points in the course of final verification:

- Embedding of the medial ischium sections
- Course of the socket brim
- Fit via control bores
- Bony tension in the medio-lateral direction

## 7 Legal Information

### 7.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

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Otto Bock HealthCare GmbH