F23 Maverick Vertical Shock

DE Gebrauchsanweisung	7
EN Instructions for use	16
FR Instructions d'utilisation	26
IT Istruzioni per l'uso	36
ES Instrucciones de uso	46
PT Manual de utilização	56
NL Gebruiksaanwijzing	66
SV Bruksanvisning	76
DA Brugsanvisning	85
NO Bruksanvisning	94
FI Käyttöohje	103
RU Руководство по применению	112

ISO 22675 - P n - m kg*)


*) Body mass limit not to be exceeded!


Stiffness	Weight (kg)	Label text
1	52	ISO 22675-P3-52-kg
2	59	ISO 22675-P3-59-kg
3	68	ISO 22675-P3-68-kg
4	77	ISO 22675-P4-77-kg
5	88	ISO 22675-P4-88-kg
6	100	ISO 22675-P5-100-kg
7	116	ISO 22675-P5-116-kg
8	130	ISO 22675-P6-130-kg
9	147	ISO 22675-P7-147-kg

1

2

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INFORMATION

Date of last update: 2022-03-21

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

1.1 Construction and Function

This document applies to the following products: F23 Maverick Vertical Shock

The prosthetic foot has fibreglass spring elements. The heel stiffness can be increased with heel wedges.

The prosthetic foot has a functional ring. The functional ring absorbs vertical impact loads and enables slight torsional movements.

1.2 Combination possibilities

This prosthetic component is compatible with Ottobock's system of modular connectors. Functionality with components of other manufacturers that have compatible modular connectors has not been tested.

Limited combination options for Ottobock components

The prosthetic foot generates high moments in the ankle area. Use structural components with higher weight limits:

Body weight (high activity) [kg]	up to 59	up to 77		up to 100		up to 130
Foot size [cm]	up to 26	up to 28	from 29	up to 28	from 29	up to 31
Structural component weight limit [kg]	≥75	≥100	≥125	≥125	≥150	≥150
Body weight (normal activity) [kg]				up to 100		up to 130
Foot size [cm]				from 29		up to 31
Structural component weight limit [kg]				≥125		≥150

2 Intended use

2.1 Indications for use

The product is intended exclusively for lower limb exoprosthetic fittings.

2.2 Area of application

Our components perform optimally when paired with appropriate components based upon weight and mobility grades identifiable by our MOBIS classification information and which have appropriate modular connectors.

The product is recommended for mobility grade 3 (unrestricted outdoor walker) and mobility grade 4 (unrestricted outdoor walker with particularly high demands).

The table that follows shows the suitable spring stiffness of the prosthetic foot, matching the patient's body weight and activity.

Spring stiffness relative to body weight and activity level		
Body weight [kg]	Normal	High
44 to 52	1	2
53 to 59	2	3
60 to 68	3	4
69 to 77	4	5
78 to 88	5	6
89 to 100	6	7
101 to 116	7	8
117 to 130	8	9
131 to 147	9	–

2.3 Environmental conditions

Storage and transport
Temperature range $-20\text{ }^{\circ}\text{C}$ to $+60\text{ }^{\circ}\text{C}$ ($-4\text{ }^{\circ}\text{F}$ to $+140\text{ }^{\circ}\text{F}$), relative humidity 20 % to 90 %, no mechanical vibrations or impacts
Allowable environmental conditions
Temperature range: $-10\text{ }^{\circ}\text{C}$ to $+45\text{ }^{\circ}\text{C}$ ($14\text{ }^{\circ}\text{F}$ to $113\text{ }^{\circ}\text{F}$)
Moisture: relative humidity: 20% to 90%, non-condensing
Chemicals/liquids: fresh water as dripping water, occasional contact with salty air (e.g. near the ocean)
Solids: dust
Prohibited environmental conditions
Chemicals/moisture: salt water, perspiration, urine, acids, soapsuds, chlorine water
Solids: dust in high concentrations (e.g. construction site), sand, highly hygroscopic particles (e.g. talcum)

2.4 Lifetime

Prosthetic foot

Depending on the patient's activity level, the maximum lifetime of the product is 3 years.

Footshell, protective sock, functional ring and O-rings

The product is a wear part, which means it is subject to normal wear and tear.

3 Safety

3.1 General safety instructions

 **CAUTION!**

Risk of injury and risk of product damage

- ▶ Note the combination possibilities/combination exclusions in the instructions for use of the products.
- ▶ Comply with the product's field of application and do not expose it to excessive strain (see page 17).
- ▶ To avoid the risk of injury and product damage, do not use the product beyond the tested lifetime.
- ▶ To avoid the risk of injury and product damage, only use the product for a single patient.
- ▶ To prevent mechanical damage, use caution when working with the product.
- ▶ If you suspect the product is damaged, check it for proper function and readiness for use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

NOTICE!

Risk of product damage and limited functionality

- ▶ Check the product for damage and readiness for use prior to each use.
- ▶ Do not use the product if its functionality is restricted. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).
- ▶ Do not expose the product to prohibited environmental conditions.
- ▶ Check the product for damage if it has been exposed to prohibited environmental conditions.
- ▶ Do not use the product if it is damaged or in a questionable condition. Take suitable measures (e.g. cleaning, repair, replacement, inspection by the manufacturer or a specialist workshop).

Signs of changes in or loss of functionality during use

Reduced spring effect (e.g. decreased forefoot resistance or changed rollover behaviour) or delamination of the spring are indications of loss of functionality. Unusual noises can indicate a loss of functionality.

4 Scope of delivery

Quantity	Designation	Reference number
1	Instructions for use	-
1	Prosthetic foot	-
1	Protective sock (black)	SL=Spectra-Sock-7
1	Heel wedge set	2F11=1

Additional accessories/spare parts (not included in the scope of delivery)

Designation	Reference number
Footshell with sandal toe (without cap)	FTC-2A-1*
Tool to change the footshell	ACC-00-10300-00
Protective sock (white)	SL=Spectra-Sock

5 Preparing the product for use

CAUTION

Incorrect alignment, assembly or adjustment

Injury due to incorrectly installed or adjusted as well as damaged prosthetic components

- ▶ Observe the alignment, assembly and adjustment instructions.

NOTICE

Altering the prosthetic foot or footshell

Premature wear resulting from damage to the product

- ▶ Do not alter the prosthetic foot or footshell.

5.1 Alignment

INFORMATION

There is a plastic protector on the product's pyramid. It protects the connecting section from scratches during the alignment and fitting of the prosthesis.

- ▶ Remove the protector before the patient leaves the workshop/fitting area.

5.1.1 Pre-compressing the functional ring

In order to achieve its full flexibility, the functional ring has to be compressed before the first use.

- 1) Using the pre-compression kit, compress the functional ring for **15 minutes** (see fig. 4, see fig. 5).
- 2) Disassemble the pre-compression kit.

5.1.2 Applying/removing the footshell

INFORMATION

- ▶ Pull the protective sock over the prosthetic foot to prevent noises in the footshell.
 - ▶ Always use the prosthetic foot with the footshell.
- ▶ Apply or remove the footshell as described in the footshell instructions for use.

5.1.3 Bench Alignment

TT bench alignment

Bench alignment process	
Required materials: 662M4 goniometer, 743S12 heel height measuring device, 743A80 50:50 gauge, alignment tool (e.g. 743L200 L.A.S.A.R. Assembly or 743A200 PROS.A. Assembly)	
Assemble and align the prosthetic components in the alignment apparatus according to the following specifications:	
①	Dorsiflexion/plantar flexion: Adjust the prosthetic foot to the heel height of the shoe. (The prosthetic foot is designed for shoes with a heel height of 10 mm.)
②	Adduction/abduction: Adjust the angle of the prosthetic socket in the frontal plane.
③	Flexion/extension: Adjust the angle of the prosthetic socket in the sagittal plane.
④	Linear shifting of the prosthetic socket: Shift the prosthetic socket so that the plumb line runs along the front edge of the adapter on the prosthetic foot (see fig. 2).

TF bench alignment

- ▶ Observe the information in the prosthetic knee joint instructions for use.

5.1.4 Static Alignment

- Ottobock recommends checking the alignment of the prosthesis using the L.A.S.A.R. Posture and adapting it as needed.

5.1.5 Dynamic Trial Fitting

- Adapt the alignment of the prosthesis in the frontal plane and the sagittal plane (e.g. by making angle or slide adjustments) to ensure an optimum gait pattern.
- **TT fittings:** Make sure that physiological knee movement in the sagittal and frontal plane is achieved when the leg begins to bear weight after the heel strike. Avoid medial movement of the knee joint. If the knee joint moves in the medial direction in the first half of the stance phase, move the prosthetic foot in the medial direction. If the medial movement occurs in the second half of the stance phase, reduce the exterior rotation of the prosthetic foot.
- Remove the plastic protector from the pyramid after completing the dynamic fitting and the walking exercises.

Heel too soft	
Symptoms	Possible solutions
<ul style="list-style-type: none">• Foot makes full ground contact too quickly• Forefoot feels too stiff• Knee goes into hyperextension	<ul style="list-style-type: none">• Shift prosthetic socket forward relative to the foot• Use a heel wedge

Heel too hard	
Symptoms	Possible solutions
<ul style="list-style-type: none">• Fast knee flexion, low stability• Transition from heel strike to toe-off too fast• Energy return feels low	<ul style="list-style-type: none">• Shift prosthetic socket back relative to the foot• Reduce the heel stiffness (shift or remove heel wedge)

Prosthetic foot too stiff	
Symptoms	Possible solutions
<ul style="list-style-type: none">• Minimal rollover movement of the prosthetic foot at low walking speed (long full ground contact)	<ul style="list-style-type: none">• Select prosthetic foot in lower stiffness

Prosthetic foot too soft	
Symptoms	Possible solutions
<ul style="list-style-type: none">• Clicking noise at the start of ground contact.• Very pronounced deformation of the forefoot at high activity levels	<ul style="list-style-type: none">• Select prosthetic foot in higher stiffness

5.1.5.1 Optimising the heel characteristics

The behaviour of the prosthetic foot at heel strike and during heel contact in the mid-stance phase can be adapted with the placement of a heel wedge.

The heel wedge is attached with adhesive tape for testing. For definitive assembly, it is glued to the prosthetic foot.

Fitting

- 1) Use the adhesive tape on the underside of the heel wedge.
- 2) Place the heel wedge on the base spring in the recommended position (see fig. 3).
- 3) Set the desired stiffness by shifting the heel wedge (anterior = harder, posterior = softer).

Definitive assembly

- 1) Remove the adhesive tape from the heel wedge. This can be done using a cloth with acetone.
- 2) Roughen the contact surface on the prosthetic foot slightly using sandpaper. Remove the dust after sanding.
- 3) Glue the heel wedge to the prosthetic foot with cyanoacrylate adhesive.

5.1.5.2 Checking the stiffness of the functional ring

The stiffness of the functional ring influences the level of comfort while walking. The stiffness is specified on the functional ring (the higher the number, the greater the stiffness).

Functional ring while walking	Cause	Correction/solution
Soft damping for the user. The functional ring is not compressed to the stop.	The functional ring has the appropriate stiffness	–
User hits the stop during damping. The functional ring is fully compressed.	Functional ring too soft	Install a harder functional ring (see page 24)
No damping for the user. The functional ring is not compressed.	Functional ring too hard	Install a softer functional ring (see page 24)
	Possible: Heel wedge too soft	Install a harder heel wedge (see page 22)

6 Cleaning

- > **Allowable cleaning agent:** pH neutral soap (e.g. 453H10 Derma Clean)
- 1) **NOTICE! To avoid product damage, only use the allowable cleaning agents.**
Clean the product with clear fresh water and a pH neutral soap.

- 2) **If present:** Remove dirt from water drainage contours using a toothpick and rinse.
- 3) Rinse the soap away with clear fresh water. In doing so, rinse the footshell until all dirt has been removed.
- 4) Dry the product with a soft cloth.
- 5) Allow to air dry in order to remove residual moisture.

7 Maintenance

- ▶ A visual inspection and functional test of the prosthetic components should be performed after the first 30 days of use.
- ▶ Inspect the entire prosthesis for wear during normal consultations.
- ▶ Inspect the product for visible damage every six months. Make additional inspection appointments as needed (e.g. for highly active or heavy users).
- ▶ Check the condition of the functional ring as part of each inspection. If the user frequently compresses the functional ring to the stop while walking, or if the functional ring exhibits visible damage, replace the functional ring (see page 24).

7.1 Replacing the functional ring

To replace the functional ring, the adapter with pyramid is removed from the base body. The O-rings on the adapter and the mounting screw are also replaced together with the functional ring.

- > **Required materials:** Pre-compression kit, torque wrench (e.g. 710D20), degreasing cleaner (e.g. 634A58 isopropyl alcohol), Loctite 243 or Loctite 246
- 1) Remove the footshell and the prosthetic foot from the prosthesis.
 - 2) Insert an Allen key (**6 mm**) into the hole in the pyramid and unscrew the mounting screw in the distal end of the base body (see fig. 6).
 - 3) Pull the adapter with pyramid up and out of the base body (see fig. 7).
 - 4) Remove the functional ring.
 - 5) **If the O-rings are worn:** Remove the O-rings from the adapter (see fig. 8).
 - 6) **If the O-ring is worn:** Remove the O-ring from the mounting screw (see fig. 8).
 - 7) Clean the mounting screw and adapter including thread with degreasing cleaner (see fig. 8).
 - 8) Slide new O-rings onto the adapter and apply silicone grease. Make sure that no silicone grease gets into the threads (see fig. 8).
 - 9) Slide a new O-ring onto the mounting screw (see fig. 8).

- 10) Apply Loctite to the thread of the mounting screw.
- 11) Slide the new functional ring onto the adapter of the prosthetic foot. The structure of the functional ring only permits installation with the correct alignment.
- 12) Insert the adapter of the prosthetic foot into the base body and attach it with the mounting screw (tightening torque: **15 Nm**). To do so, insert the torque wrench through the hole in the pyramid.
- 13) Pre-compress the new functional ring (see page 21).
- 14) Reinstall the prosthetic foot on the prosthesis and pull on the footshell.

8 Disposal

In some jurisdictions it is not permissible to dispose of the product with unsorted household waste. Improper disposal can be harmful to health and the environment. Observe the information provided by the responsible authorities in your country regarding return, collection and disposal procedures.

9 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

9.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

9.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

10 Technical data

F23 Maverick VS								
Sizes [cm]	23	24	25	26	27	28	29	30
Heel height [mm]	10							
System height [mm]	181		189			202		
Build height [mm]	209		217			230		
Average weight with footshell [g]	920		1072			1166		
Max. body weight [kg]	147							

F23 Maverick VS								
Sizes [cm]	23	24	25	26	27	28	29	30
Mobility grade	3, 4							

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