

### Flat seam

To ensure the proximal end is barely felt and the distal seal is extremely precise.

### Pre-shaped kneecap area

Reduces pressure on the kneecap when the knee bends.

### Smooth, elastic fabric

Makes it easier to put on, facilitates mobility and ensures a long service life.

### Medical white oil and material quality

The white oil is released gradually. This helps nourish the skin and retains the wall thickness for a long time.

### 15° preflexion

Makes it easier to bend the knee, and allows fewer creases to form in the hollow of the knee.

### Conical shape

The sealing sleeve is tapered in two directions: from the distal end in the direction of the knee, and from the knee in the direction of the proximal end. This increases the grip on the prosthetic socket, and at the same time reduces pressure on the thigh.







## ProFlex Plus sealing sleeve

### Flex it!

The ProFlex is one of the most widely used sealing sleeves in the world and now we have made it better!

Introducing the new ProFlex Plus. Designed with a more flexible fabric and smoother proximal seam the ProFlex Plus offers a soft and supple but tough sealing sleeve. ProFlex Plus is an extension of the ProFlex line and is now available in Tan!

Technical data			
<b>Article number</b>	<b>453A30</b>	<b>453A40</b>	
<b>Thigh length from MPT to proximal brim</b>	Regular approx. 26 cm	Short approx. 18 cm	
<b>Size</b>	1	2	3
<b>Knee center circumference</b>	24–32 cm	30–40 cm	36–47 cm
<b>Circumference 20 cm proximal to knee center</b>	34–46 cm	40–54 cm	48–66 cm
<b>Colors</b>	tan (0)		
<b>Order example</b>	453A30=2-0		

<b>Liner</b>	 TPE 6Y92	 DVS Liner 6Y94	 Anatomic 3D 6Y512
<b>Sleeve</b>	 ProFlex Plus 453A30/453A40		
<b>Suspension</b>	 Dynamic Vacuum 4R220	 Harmony P3 4R147	