

ottobock.



bebionicTM standard gloves

Crafted using multiple layers of advanced silicone material, and constructed with an integrated reinforcing mesh, the bebionic standard gloves are soft and durable, but remain easy to remove and clean.

Available in 19 different lifelike color shades, the unique TrueFinishTM micropigmentation makes the glove the ultimate in authenticity. Additional detailing on the palms, knuckles, nails, and joints enhances the natural appearance of the hand. These gloves are also available in a jet black shade.

Quality for life

bebionic standard gloves

Part number	Side	Description
BBGSML/*	Left	bebionic small glove – left
BBGSMR/*	Right	bebionic small glove – right
BBGMDL/*	Left	bebionic medium glove – left
BBGMDR/*	Right	bebionic medium glove – right
BBGLGL/*	Left	bebionic large glove – left
BBGLGR/*	Right	bebionic large glove – right
BBGTSML/BK	Left	bebionic small glove, jet black – no definition – left
BBGTSMR/BK	Right	bebionic small glove, jet black – no definition – right
BBGTMDL/BK	Left	bebionic medium glove, jet black – no definition – left
BBGTMDR/BK	Right	bebionic medium glove, jet black – no definition – right
BBGTLGL/BK	Left	bebionic large glove, jet black – no definition – left
BBGTLGR/BK	Right	bebionic large glove, jet black – no definition – right



19 TrueFinish™ color shades

Swatches are available on request:
E24655 E0-E18 TrueFinish colors

Ordering examples
BBGMDL/E5
Orders a bebionic medium glove – left in TrueFinish color E5

Swatch colors are intended for illustrative purposes only.
*Black color only available for bebionic standard gloves.

To request swatch samples or to place an order for bebionic standard gloves, please contact your local customer representative.