

Foot Family Portfolio



Image	K-Level	Material Number	Foot Name	Foot Description	Max Pt Weight	Inches - Clearance (size 26)	L-Codes
	K3	1B1-2	Meridium®	Microprocessor Feet	275 lbs	6 1/4" (160 mm)	L5973
	K3 (K4)	1A1-1	Empower®	Microprocessor Feet	287 lbs	8 3/4" (222 mm)	L5999
	K3	1C66	Triton® smart ankle	Microprocessor Feet	220 lbs	6 1/4" (158 mm)	L5973
	K3-K4	1C60	Triton®	Moderate-High Activity	330 lbs	6 1/4" (158 mm)	L5980 + L5986
	K3-K4	1C61	Triton® Vertical Shock	Moderate-High Activity	330 lbs	8 1/16" (205 mm)	L5986 - L5987
	K3-K4	1C62	Triton Harmony®	Moderate-High Activity	330 lbs	8 1/16" (205 mm)	L5986 + L5987 + L5781
	K3-K4	1C63	Triton® Low Profile	Moderate-High Activity	330 lbs	2 7/8" (73 mm)	L5981 + L5986
	K3-K4	1C64	Triton® HD	Moderate-High Activity	330 lbs	5 7/8" (149 mm)	L5980 + L5986
	K3-K4	1E57	Lo Rider™	Moderate-High Activity	300 lbs	1 3/8" (36 mm) (with pyramid)	L5981
	K3-K4	1C20	ProSymes™	Specialty	275 lbs	1 3/4" (43 mm)	L5981
	K3-K4	1E80	Chopart™ 00mm HH	Specialty	300 lbs	3/4" (20 mm)	L5976
	K3	1C30	Trias®	Moderate-High Activity	275 lbs	4 3/8" (111 mm)	L5981 + L5986
	K3	1D35	Dynamic Motion™	Moderate-High Activity	220 lbs	3 3/8" (87 mm)	L5979
	K1-K2	1C11	Terion™ K2	Lower Activity	385 lbs	3" (77 mm)	L5972 + L5986
	K1-K2	1M10	Adjust®	Lower Activity	275 lbs	3" (74 mm) (without Footshell)	L5972 + L5986
	K1-K2	1D10/1D11	Dynamic Foot™ with Toes	Lower Activity	330 lbs	3 1/2" (90 mm) (without adapter)	L5972
	K1-K2	1A30	Greissinger™ Plus Foot HH 10mm	Lower Activity	165 lbs - 220 lbs	3 3/8" (87 mm)	L5972 + L5986
	K1	1G6	Cosmetic™ Light Foot with Toes	Lower Activity	165 lbs	3 3/8" (85 mm) (w/ 2R54/2R31)	5971 (Replacement) Included in base code
	K1	1H38	Single Axis Foot only Men 10mm Toes	Lower Activity	220 lbs	2 1/2" (65 mm) (w/ 2R33/2R10) 2 3/4" (69 mm) (w/ 2R51)	L5974
	K1-K2	1S101/102 /103	SACH+™ Feet	Lower Activity	275 lbs	3 3/8" (85 mm)	5971 (Replacement) Included in base code
	K1-K2	1S30	SACH™ Foot Children	Pediatric	77 lbs	2 1/8" (55 mm) (size 12)	5971 (Replacement) Included in base code
	K1-K2	1K10	Dynamic Foot™ for Children	Pediatric	99 lbs	2 5/8" with 2R40 (66 mm) (size 18)	L5972
	K1-K3	1WR95	Foot for water resistant prosthesis	Specialty	330 lbs	3 3/8" (85 mm)	L5972