

# Ottobock Software System Requirements

## Supported Operating Systems and Minimum System Requirements

This document contains the minimum system requirements and supported operating systems for:

- AxonSoft 2.2
- C-Soft 2.8
- C-Soft Plus 1.4
- K-Soft 1.2
- ElbowSoft 1.6
- MyolinoSoft 1.0
- M-Soft 1.2
- X-Soft 1.6

*Note: Ottobock Software is not compatible with Apple products (Mac, iPad, iPhone, etc.)*

### Supported Operating Systems (32-bit and 64-bit)

- Windows 7, SP1 or higher, all editions except Windows 7 Starter
- Windows 8
- Windows 10

### Minimum PC Hardware Requirements (32-bit platforms)

- PC with a Pentium IV processor
- At least 1 GB RAM
- At least 1.5 GB free hard drive space
- Screen resolution of 1024 x 768 or higher
- Display scale 100% or 96 DPI
- 1 available USB port

### Minimum PC Hardware Requirements (64-bit platforms)

- PC with an Intel Core2 (x64) or faster
- At least 2 GB RAM
- At least 1.5 GB free hard drive space
- Screen resolution of 1024 x 768 or higher
- Display scale 100% or 96 DPI
- 1 available USB port