

Aluminum adapters

New design, lightweight and higher weight limit.



Technical data:

Article number	4R72=32AL	4R72=45AL	4R72=60AL	4R72=75AL	4R72=90AL	2R49=AL	2R50=AL	4R69=AL	4R74=AL
Product description	Double adapter					Pylon Long	Pylon Short	Tube clamp adapter	Socket adapter with pyramid adapter
Material						Aluminum			
Diameter	-	-	-	-	-	30 mm	30 mm	30 mm	-
Build height	32 mm	45 mm	60 mm	75 mm	90 mm	454 mm	214 mm	47 mm	26 mm
Weight	80 g	95 g	105 g	110 g	120 g	255 g	155 g	75 g	70 g
Max. body weight	136 kg / 300 lbs								

Today's state of the art lower limb prostheses not only benefit from sophisticated functional units, but also from high-performance modular adapters. They are key elements when it comes to lightweight design and stability leading to a high degree of wearing comfort for the user.

Benefits at a glance

- Aluminum lightweight design with high degree of stability
- Bodyweight now up to 136 kg / 300 lbs
- Reliably tested and certified in accordance with the most recent ISO 10328:2006 standards