

Caleo 3D 6Y95

TPE Comfort for Active Users

The liner is the most intimate component of the prosthesis and can have considerable impact on the success of the fitting. Caleo liners made of copolymer are valued by many users thanks to their protective and nurturing qualities.

The Caleo range has now been expanded to include a further version, the Caleo 3D 6Y95 liner. This liner with new gel formula is particularly suited to more active users due to its high level of durability and resistance to wear.

All benefits at a glance

- The 3D design features a thicker wall (6 mm) in the front for better protection, and a thinner wall (3 mm) at the rear for more flexibility
- Medical white oil for nurturing dry and sensitive residual limbs
- New, improved gel formula combines durability and comfort, and can also be thermoformed
- Knee flexion zone for easy knee flexibility with simultaneous relief of the patella
- Monoelastic fabric prevents pistoning and eliminates the need for a distal matrix
- Softer flexible cap cover for locking liners or no cap for cushion liners



Technical data

Article number	6Y95=L6 (size)	6Y95=C6*
Amputation level	transtibial (TT)	transtibial (TT)
Distal connection	with	without
Textile cover	with	with
Wall thickness	6 mm (front), 3 mm (back)	6 mm (front), 3 mm (back)
Socket design	specific or total surface weight-bearing socket	specific or total surface weight-bearing socket
Matrix	integrated in textile	integrated in textile
Example for ordering	6Y95=L6-M	6Y95=C6-M

Sizes (circumference)

	4 cm from distal		30 cm from distal	
	min	max	min	max
S	150	260	S	200 330
M	180	310	M	230 410
MP	200	310	MP	330 560
L	230	350	L	360 560
LP	280	430	LP	400 630
XL	330	500	XL	430 690