

623T1 - Cotton Stockinette

Material number 623T 1

Page: 1 of 6

1. Product and company identification

Product identifier

Trade name: 623T1 - Cotton Stockinette

Relevant identified uses of the substance or mixture and uses advised against

General use: Article for orthopedic procedures. For commercial user only.

Details of the supplier of the safety data sheet

Company name: Otto Bock Health Care

Street/POB-No.: 3820 W. Great Lakes Drive

Postal Code, city: Salt Lake City, UT 84120
USAWWW: www.ottobockus.com

Telephone: +1 (801) 956-2400

Telefax: +1 (801) 956-2401

Dept. responsible for information:

Quality Department,

Telephone: +1 (801) 954-2304 (7 AM – 3 PM, Mountain Time),

Email: USRegulatory@ottobock.com

Additional information:

Corporate headquarters:

Ottobock SE & Co. KGaA

Max-Näder-Straße 15

Duderstadt

Germany

Emergency phone number**CHEMTREC, Telephone: +1 (800) 424-9300**

2. Hazards identification

Emergency overview

Appearance: Form: solid, fibers

Color: white

Odor: odorless

Classification: This material is classified as not hazardous.

Regulatory status

This material is not considered hazardous by the U.S. OSHA Hazard Communication Standard (29 CFR 1910.1200) and SIMDUT in Canada.

Hazards not otherwise classified

Processing, e.g. by cutting, sawing or grinding, can produce particles and dust. For risks which have to be observed thereby, see section 7: Handling, section 8: Exposure controls / personal protection and section 11: Toxicology.

see section 11: Toxicological information

3. Composition / Information on ingredients

Chemical characterization: Cotton-fabric (99%-100%), contains in traces coating agent (paraffin)

4. First aid measures

General information: For mechanical processing: dust formation.

In case of inhalation: In case of troubles after inhalation of dust: Provide fresh air. Seek medical attention.

Following skin contact: Dust: Remove residues with water. In case of skin irritation, consult a physician.

After eye contact: Dust: Immediately flush eyes with plenty of flowing water for 10 to 15 minutes holding eyelids apart.
In case of troubles or persistent symptoms, consult an ophthalmologist.

After swallowing: Not a probable route of exposure

Most important symptoms/effects, acute and delayed

Processing, e.g. by cutting, sawing or grinding, can produce particles and dust.
Dust: May cause irritations.

Information to physician

Treat symptomatically.

5. Fire fighting measures

Flash point/flash point range: No data available

Auto-ignition temperature: No data available

Suitable extinguishing media: Water spray jet, foam, dry chemical powder, carbon dioxide.

Extinguishing media which must not be used for safety reasons: Full water jet

Specific hazards arising from the chemical

Combustible. May form dangerous gases and vapours in case of fire.
In case of fire may be liberated: carbon monoxide and carbon dioxide.

Protective equipment and precautions for firefighters: Wear self-contained breathing apparatus. Wear appropriate protective equipment.

Additional information: Cool endangered containers with water spray and, if possible, remove from danger zone.
Do not breathe fumes.
Do not allow fire water to penetrate into surface or ground water.

6. Accidental release measures

Personal precautions: At processing: Do not breathe dust. Ensure adequate ventilation, especially in confined areas.

Environmental precautions: Do not allow to penetrate into soil, waterbodies or drains.

Methods for clean-up: Take up mechanically, placing in appropriate containers for disposal.

7. Handling and storage

Handling

Advices on safe handling: Handle in accordance with good industrial hygiene and safety practice.

For mechanical processing: Provide adequate ventilation. Avoid generation of dust. Do not breathe dust.

Precautions against fire and explosion:

Take precautionary measures against static discharges.

Storage

Requirements for storerooms and containers:

Keep container dry.

Hints on joint storage:

Keep away from food, drink and animal feedingstuffs.

8. Exposure controls / personal protection

Engineering controls

For mechanical processing: Provide adequate ventilation.

See also information in chapter 7, section storage.

Personal protection equipment (PPE)

Eye/face protection: Tightly sealed safety glasses according to OSHA Standard - 29 CFR: 1910.133 or ANSI Z87.1-2003

Skin protection: For mechanical processing: Wear suitable protective clothing.

For machine processing:

Protective gloves against mechanical risks. OSHA Standard - 29 CFR: 1910.138

Observe glove manufacturer's instructions concerning penetrability and breakthrough time.

Respiratory protection: For mechanical processing:

Use self-contained respiratory apparatus if respirable dust type occurs.

Dust mask or particulates filter.

OSHA Standard - 29 CFR: 1910.134 or ANSI Z88.2

General hygiene considerations:

Avoid generation of dust.

Do not breathe dust. Wash hands before breaks and after work.

9. Physical and chemical properties

Information on basic physical and chemical properties

Appearance: Form: solid, fibers

Color: white

Odor: odorless

Odor threshold: No data available

pH value: No data available

Melting point/freezing point: No data available

Initial boiling point and boiling range: No data available

Flash point/flash point range: No data available

Evaporation rate: No data available

623T1 - Cotton Stockinette

Material number 623T 1

Flammability:	No data available
Explosion limits:	No data available
Vapor pressure:	No data available
Vapor density:	No data available
Density:	at 68 °F: 1.52 - 1.54 g/cm ³
Water solubility:	at 68 °F: insoluble
Partition coefficient: n-octanol/water:	No data available
Auto-ignition temperature:	No data available
Thermal decomposition:	No data available
Additional information:	No data available

10. Stability and reactivity

Reactivity:	No data available.
Chemical stability:	Stable under recommended storage conditions.
Possibility of hazardous reactions	No dangerous reactions with proper and specified storage and handling.
Conditions to avoid:	Protect from moisture contamination.
Incompatible materials:	Strong acids, strong bases
Hazardous decomposition products:	In case of fire may be liberated: carbon monoxide and carbon dioxide.
Thermal decomposition:	No data available

11. Toxicological information

Toxicological tests

Toxicological effects:	Acute toxicity (oral): Lack of data.
	Acute toxicity (dermal): Lack of data.
	Acute toxicity (inhalative): Lack of data.
	Skin corrosion/irritation: Lack of data.
	Serious eye damage/irritation: Lack of data.
	Sensitisation to the respiratory tract: Lack of data.
	Skin sensitisation: Lack of data.
	Germ cell mutagenicity/Genotoxicity: Lack of data.
	Carcinogenicity: Lack of data.
	Reproductive toxicity: Lack of data.
	Effects on or via lactation: Lack of data.
	Specific target organ toxicity (single exposure): Lack of data.
	Specific target organ toxicity (repeated exposure): Lack of data.
	Aspiration hazard: Lack of data.

623T1 - Cotton Stockinette

Material number 623T 1

Page: 5 of 6

Symptoms

In case of inhalation: Dust: After intake of large amounts: mild irritant
After eye contact: Dust: mild irritant

12. Ecological information**Ecotoxicity**

Further details: No data available

Mobility in soil

No data available

Persistence and degradability

Further details: Product is biodegradable.

Additional ecological information

General information: Do not allow to penetrate into soil, waterbodies or drains.

13. Disposal considerations**Product**

Recommendation: Incinerate according to applicable local, state and federal regulations.

Contaminated packagingRecommendation: Dispose of waste according to applicable legislation.
Non-contaminated packages may be recycled.**14. Transport information****USA: Department of Transportation (DOT)**

Proper shipping name: Not restricted

Sea transport (IMDG)

Proper shipping name: Not restricted

Marine pollutant: no

Air transport (IATA)

Proper shipping name: Not restricted

Further information

No dangerous good in sense of these transport regulations.

15. Regulatory information**National regulations - U.S. Federal Regulations**

cotton: TSCA Inventory: not listed

623T1 - Cotton Stockinette

Material number 623T 1

Page: 6 of 6

National regulations - Great Britain

Hazchem-Code: -

16. Other information

Hazard rating systems:



NFPA Hazard Rating:

Health: 0 (Minimal)

Fire: 1 (Slight)

Reactivity: 0 (Minimal)

HMIS Version III Rating:

Health: 0 (Minimal)

Flammability: 1 (Slight)

Physical Hazard: 0 (Minimal)

Personal Protection: X = Consult your supervisor

HEALTH	0
FLAMMABILITY	1
PHYSICAL HAZARD	0
	X

Reason of change: Changes in section 1.3: Corporate headquarters

Date of first version: 1/9/2008

Department issuing data sheet

Contact person: see section 1: Dept. responsible for information

The information in this data sheet has been established to our best knowledge and was up-to-date at time of revision. It does not represent a guarantee for the properties of the product described in terms of the legal warranty regulations.