

# RS2 Runway

**Product features**

- Adjust heel height from zero to 2 inches; option of push button or 4 mm hex key adjustable design
- Push button design offers 10 height settings
- Runway HX offers infinite adjustments on the ankle slide
- Anatomical “gliding” ankle design maintains appropriate prosthetic alignment, regardless of heel height



**Technical data**

K-Level	3
Max. body weight	255 lbs / 116 kg
Sizes	22-28 cm
Heel height	0"-2" / 0-50 mm
Weight with footshell*	27.94 oz / 792 g
Build height*	3 15/16" / 100 mm

\* Reference size is 26 cm.

## RS2 Runway

### Order form

Customer		Shipping address (if different from customer address)	
Customer no.	<input type="text"/>	Customer no.	<input type="text"/>
Company	<input type="text"/>	Company	<input type="text"/>
Street	<input type="text"/>	Street	<input type="text"/>
Postal code/city	<input type="text"/>	Postal code/city	<input type="text"/>
O&P professional	<input type="text"/>	Customer ref.	<input type="text"/>

#### Stiffness chart

Body weight	Activity level		
	Low	Moderate	High
100–115 lbs / 44–52 kg	1	1	2
116–130 lbs / 53–59 kg	1	2	3
131–150 lbs / 60–68 kg	2	3	4
151–170 lbs / 69–77 kg	3	4	5
171–195 lbs / 78–88 kg	4	5	6
196–220 lbs / 89–100 kg	5	6	–
221–255 lbs / 101–116 kg	6	–	–

#### Foot module part number

Reference number	Product ID	Stiffness	Size	Toe Configuration
RS2	-	0	A	-
	<b>Product ID</b>	<b>Stiffness</b>	<b>Size (cm)</b>	<b>Toe Configuration</b>
	00 Runway	1–6	22, 23, ..., 28	RL Regular toe, left
	K1 Runway HX			RR Regular toe, right
				SL Sandal toe, left
				SR Sandal toe, right

#### Footshell part number

Reference number	Generation	Profile	Cap	Size	Color	Toe	Side
FTC	-	-	-	-	4	-	-
	<b>Generation</b>	<b>Profile</b>	<b>Cap</b>	<b>Size (cm)</b>	<b>Color</b>	<b>Toe</b>	<b>Side</b>
	2 Sandal toe	D Runway	0 Yes	22, 23, ..., 28	L Light	R Regular	L Left
	3 Standard (mid profile)	E Runway HX	1 No		M Medium	S Sandal	R Right
					D Dark		

#### Scope of delivery

- RS2-0\* Runway or RS2-K1\* Runway HX
- KIT-00-113\* Set of heel wedges
- SO-NPS-200\* Spectra protective sock black

#### Accessories & spare parts

- SO-NPS-100\* Spectra protective sock white
- ACC-00-10300-00 Footshell removal tool
- ACC-00-10400-00 Shim adjustment wrench

#### Available footshell variants

- FTC-3D-1\*: Footshell Runway standard (without cap)
- FTC-2D-1\*-S\*: Footshell Runway sandal toe (without cap)
- FTC-3E-1\*: Footshell Runway HX standard (without cap)
- FTC-2E-1\*-S\*: Footshell Runway HX sandal toe (without cap)