


VS1 Senator

Product features

- Lightweight, energy-return design facilitates smooth, easy stride and durability

**Technical data**

K-Level	3
Max. body weight	300 lbs / 136 kg
Sizes	22-31 cm
Heel height	3/8" / 10 mm
Weight with footshell*	23.10 oz / 655 g
Build height*	5 1/2" / 140 mm

* Reference size is 26 cm.

VS1 Senator

Order form

Customer		Shipping address (if different from customer address)	
Customer no.		Customer no.	
Company		Company	
Street		Street	
Postal code/city		Postal code/city	
O&P professional		Customer ref.	

Stiffness chart

Body weight	Activity level	
	Low	Moderate
100–120 lbs / 45–55 kg	1	2
121–145 lbs / 56–66 kg	2	3
146–175 lbs / 67–80 kg	3	4
176–210 lbs / 81–96 kg	4	5
211–250 lbs / 97–114 kg	5	6
251–300 lbs / 115–136 kg	6	–

Foot module part number

Reference number	Stiffness	Size
VS1 - 00 - 0	A	RU
	Stiffness	Size (cm)
	1–6	22, 23, ..., 31

Footshell part number

Reference number	Generation	Profile	Cap	Size	Color	Side
FTC	-	-	-	4	-	R
	Generation	Profile	Cap	Size (cm)	Color	Side
	2–3	M Mid	0 Yes	22, 23, ..., 31	L Light	L Left
		H High	1 No		M Medium	R Right
					D Dark	

Scope of delivery

- VS1* Senator
- KIT-00-113* Set of heel wedges
- SO-NPS-200* Spectra protective sock black

Accessories & spare parts

- SO-NPS-100* Spectra protective sock white
- ACC-00-10300-00 Footshell removal tool

Available footshell variants

- FTC-3M-0*: Footshell standard (with cap)
- FTC-3M-1*: Footshell standard (without cap)
- FTC-2H-1*: Footshell high profile (without cap)